

MacroModel 9.1

User Manual

Copyright © 2006 Schrödinger, LLC. All rights reserved. CombiGlide, Epik, Glide, Impact, Jaguar, Liaison, LigPrep, Maestro, Phase, Prime, QikProp, QikFit, QikSim, QSite, SiteMap, and Strike are trademarks of Schrödinger, LLC.

Schrödinger and MacroModel are registered trademarks of Schrödinger, LLC.

Python is a copyrighted work of the Python Software Foundation. All rights reserved.

The C and C++ libraries for parsing PDB records are a copyrighted work (1989) of the Regents of the University of California. All rights reserved.

See the [Copyright Notices](#) for full copyright details.

To the maximum extent permitted by applicable law, this publication is provided “as is” without warranty of any kind. This publication may contain trademarks of other companies.

Please note that any third party programs (“Third Party Programs”) or third party Web sites (“Linked Sites”) referred to in this document may be subject to third party license agreements and fees. Schrödinger, LLC and its affiliates have no responsibility or liability, directly or indirectly, for the Third Party Programs or for the Linked Sites or for any damage or loss alleged to be caused by or in connection with use of or reliance thereon. Any warranties that we make regarding our own products and services do not apply to the Third Party Programs or Linked Sites, or to the interaction between, or interoperability of, our products and services and the Third Party Programs. Referrals and links to Third Party Programs and Linked Sites do not constitute an endorsement of such Third Party Programs or Linked Sites.

Revision A, April 2006

Contents

Document Conventions	xiii
Chapter 1: MacroModel Overview	1
1.1 MacroModel User Manual	1
1.2 MacroModel	1
1.3 MacroModel and Maestro Interaction	2
1.4 Calculation Preparation and Submission	2
1.5 Command-Line Utilities	3
1.6 Citing MacroModel in Publications	4
Chapter 2: Introduction to Maestro.....	5
2.1 General Interface Behavior	5
2.2 Starting Maestro	5
2.3 The Maestro Main Window	6
2.3.1 The Menu Bar	8
2.3.2 The Toolbar	9
2.3.3 Mouse Functions in the Workspace	12
2.3.4 Shortcut Key Combinations	13
2.4 Maestro Projects	13
2.4.1 The Project Table Toolbar	15
2.4.2 The Project Table Menus	16
2.4.3 Selecting Entries	17
2.4.4 Including Entries in the Workspace	18
2.4.5 Mouse Functions in the Project Table	18
2.4.6 Project Table Shortcut Keys	19
2.5 Building a Structure	20
2.5.1 Placing and Connecting Fragments	20
2.5.2 Adjusting Properties	22
2.5.3 The Build Panel Toolbar	22

2.6 Selecting Atoms	23
2.6.1 Toolbar Buttons	23
2.6.2 Picking Tools	24
2.6.3 The Atom Selection Dialog Box	25
2.7 Scripting in Maestro	25
2.7.1 Python Scripts	25
2.7.2 Command Scripts	26
2.7.3 Macros	27
2.8 Specifying a Maestro Working Directory	27
2.9 Undoing an Operation	28
2.10 Running and Monitoring Jobs	28
2.11 Getting Help	30
2.12 Ending a Maestro Session	30
Chapter 3: Basic Molecular Modeling	31
3.1 MacroModel Force Field Implementation	31
3.2 Parameter Quality Considerations	33
3.3 MacroModel Solvation Treatment	35
3.4 Truncation of Electrostatic Interactions	36
3.4.1 Energetics of Charges	37
3.4.2 Molecular Mechanics Description of Charge	37
3.4.3 Bond Dipole Cutoffs (BDCO)	38
3.4.4 BDCO and Molecular Modeling	38
3.5 General Guidelines for Convergence	39
3.5.1 Minimization	39
3.5.2 Conformational Searching	39
3.5.3 Molecular and Stochastic Dynamics	40
3.6 Problem Size: Limits of Operation	41
3.7 Nonbonded Cutoffs	41

3.8	Electrostatic Treatments	42
3.9	File Conversion	42
3.10	Modified Force Field, Solvent, or Atom Type Files	43
3.11	Hydrogen Treatment	43
3.11.1	Hydrogens on Hetero Atoms.....	43
3.11.2	Hydrogen Treatment in the AMBER* Force Field.....	43
Chapter 4: Running MacroModel From the Command Line.....		45
4.1	Preparing for MacroModel Calculations	45
4.1.1	Environment Variables	45
4.1.2	The Command File	45
4.1.3	The MacroModel Command File Format	46
4.2	Submitting Jobs From the Command Line.....	47
4.2.1	Using an Executable Command File.....	49
4.3	Running Remote and Distributed Calculations	50
4.4	Interacting With a Running MacroModel Job.....	51
Chapter 5: General Energetic Settings.....		53
5.1	General MacroModel Calculation Setup.....	53
5.1.1	Specifying Job Input Source	53
5.2	The Potential Folder	54
5.3	Fixing and Freezing Atoms in Large Systems	55
5.3.1	Methods for Freezing and Fixing Atoms	56
5.3.2	Freezing and Fixing Individual Atoms	57
5.3.3	The Substructure Facility	57
5.4	The Constraints Folder.....	59
5.4.1	Constraining Atoms, Distances, Angles, and Torsions.....	60
5.4.2	Freezing Atoms	61

5.5 The Substructure Folder	61
5.5.1 Defining a Substructure	61
5.5.2 Creating a Shell of Atoms	62
5.5.3 Reading and Writing Substructures	63
5.6 Setting up the Job	64
5.7 Writing Job Files	65
5.8 Running and Monitoring Jobs	65
Chapter 6: Current Energy Calculations	67
6.1 The Current Energy Panel	67
6.2 The ECalc Folder	67
6.3 Command File Examples	68
6.4 The Force Field Viewer	70
6.5 Checking and Interpreting Results	71
Chapter 7: Minimizations	73
7.1 The Minimization Panel	73
7.2 Performing a Minimization Calculation	73
7.2.1 Minimization Methods	74
7.2.2 Convergence Parameters	75
7.3 Command File Examples	76
7.4 Checking and Interpreting Results	77
Chapter 8: Multiple Minimizations	79
8.1 The Multiple Minimization Panel	79
8.2 Setting Up Multiple Minimization Calculations	79
8.3 Partition Coefficient Estimation	82
8.4 Automatic Setup (AUTO)	83

8.5 Preparation of Multi-Ligand Structure Files With premin	84
8.6 Command File Examples	85
8.6.1 Energy Minimization of Multiple Non-Conformers.....	85
8.6.2 Multiple Conformer Minimization With Automatic Redundant Conformer Elimination	86
8.6.3 Partition Coefficient Estimation	87
 Chapter 9: Geometry Scans	89
9.1 Performing a Dihedral Scan in Maestro	89
9.2 Command File Examples	90
9.3 Plotting Scan Results in Maestro	92
9.3.1 1D and 2D Plot Panel General Features.....	92
9.3.2 The 1D Plot Panel.....	93
9.3.3 2D Plot Panel	94
9.4 Checking and Interpreting Results	95
 Chapter 10: Conformational Searches	97
10.1 Conformational Search Methods.....	97
10.2 Performing Conformational Searches	98
10.2.1 Conformational Search Method	99
10.2.2 Global Search Restrictions	101
10.2.3 Low-Mode Parameters.....	102
10.2.4 Automatic Setup of Conformational Search Variables	103
10.2.5 Setting Conformational Search Variables Manually	104
10.2.5.1 Ring Closures	105
10.2.5.2 Torsion Rotations	106
10.2.5.3 Molecule Trans/Rot	107
10.2.5.4 Comparison Atoms	109
10.2.5.5 Chiral Atoms	110
10.2.5.6 Distance Check.....	111
10.2.5.7 Torsion Check.....	112
10.2.5.8 Ligand Bonds.....	113

10.3 Command File Examples	113
10.3.1 Conformational Search Using MCMM	114
10.3.2 Multi-structure Conformational Search Using LMOD	116
10.3.3 Mixed MCMM/Low-Mode Search Using a Substructure File	117
10.4 Checking and Interpreting Results	118
Chapter 11: Ligand Torsional Searches	119
11.1 The Ligand Torsion Search Panel	120
11.1.1 The Potential Folder	120
11.1.2 The Mini Folder	121
11.1.3 The LTSearch Folder	121
11.2 Setting up ConfGen Calculations	124
11.3 Command File Example	124
Chapter 12: Dynamics Calculations	129
12.1 The Dynamics Panel	129
12.2 Performing a Dynamics Calculation	129
12.2.1 The Monitor Folder	129
12.2.1.1 Surface Areas	131
12.2.1.2 Distances	131
12.2.1.3 Angles	132
12.2.1.4 Torsions	132
12.2.1.5 H-Bonds	133
12.2.2 The Dynamics Folder	134
12.3 Command File Examples	136
12.3.1 Stochastic Dynamics	136
12.3.2 Simulated Annealing	137
12.4 Checking and Interpreting Results	138
Chapter 13: MC/SD Calculations	139
13.1 The MC/SD Panel	139

13.2	Setting Up MC/SD Calculations	139
13.3	Command File Examples	140
13.4	Checking and Interpreting Results	142
Chapter 14: Minta Calculations.....		143
14.1	The MINTA Panel	143
14.2	Command File Examples	145
14.3	Checking and Interpreting Results	146
Chapter 15: Protein Loop Construction		147
15.1	Performing LOOP Calculations	147
15.1.1	Input Restrictions	147
15.1.2	Structure Handling Considerations	148
15.2	Example .lsq Input File.....	148
15.3	Command File Examples	149
15.3.1	LOOP Job Using the Input Structure Sequence	149
15.3.2	LOOP Job Using the Sequence From a File.....	151
15.4	Example LOOP Job Output	151
15.4.1	Atom Renumbering Using Input Protein Structure.....	151
15.4.2	Conformational Search Using Input Protein Structure	151
15.4.3	LOOP Run Using an .lsq File as Input	152
Chapter 16: eMBRAcE		153
16.1	Minimizations With eMBRAcE.....	153
16.2	Conformational Searches With eMBRAcE.....	155
16.3	Specifying a Substructure for eMBRAcE.....	158
16.4	Command File Examples	160
16.4.1	eMBRAcE Minimization Calculations	160
16.4.1.1	Interaction Mode	160
16.4.1.2	Energy Difference Mode.....	163

16.4.2	Conformational Searches With eMBrAcE	165
16.4.2.1	MCMM Conformational Search	165
16.4.2.2	Low-Mode Conformational Search	168
16.4.2.3	MCMM Conformation Search With COPY/ALGN	169
16.4.3	Distributed eMBrAcE Calculations	171
Chapter 17: Free Energy Simulations		175
17.1	Free Energy Perturbation	175
17.2	Setting Up FEP Calculations	176
17.2.1	The Structure File	176
17.2.2	The Command File	177
17.2.3	Other Possible Command Files	179
17.2.4	The Output File	180
17.3	Criteria for a Successful Simulation	182
17.4	Other Types of Free Energy Calculations	184
17.4.1	Hydrogen Bonding Preference of a Glycyl Lactam in Organic Solution	184
17.4.2	Conformational Free Energies for a Diamide in Solution	185
17.5	Literature	186
Chapter 18: Molecular Clustering with XCluster		189
18.1	The XCluster Panel	189
18.1.1	Selecting the Structure Source	189
18.1.2	Selecting the Clustering Criterion	190
18.1.3	Defining the Comparison Atoms or Torsions	191
18.2	Command File Examples	191
18.3	XCluster Output	192
Chapter 19: Redundant Conformer Elimination		193
19.1	Eliminating Redundant Conformers Using Maestro	193
19.2	Command File Examples	195

Chapter 20: Additional Features	197
20.1 Geometry Calculations.....	197
20.2 Calculating Interaction Energies Using ASET	197
20.3 Rewinding the Output File for Additional Minimization	198
20.4 Visualizing Vibrational Modes	199
20.5 Alignment of Structures.....	200
20.6 autoref: Restrained Minimizations	201
20.7 serial_split: Split Serial Search Output Structure Files	201
Chapter 21: Getting Help	203
References	205
Copyright Notices	209
Index.....	211

Document Conventions

In addition to the use of italics for names of documents, the font conventions that are used in this document are summarized in the table below.

Table 1.1.

Font	Example	Use
Sans serif	Project Table	Names of GUI features, such as panels, menus, menu items, buttons, and labels
Monospace	\$SCHRODINGER/maestro	File names, directory names, commands, environment variables, and screen output
Italic	<i>filename</i>	Text that the user must replace with a value
Sans serif uppercase	CTRL+H	Keyboard keys

In descriptions of command syntax, the following UNIX conventions are used: braces { } enclose a choice of required items, square brackets [] enclose optional items, and the bar symbol | separates items in a list from which one item must be chosen. Lines of command syntax that wrap should be interpreted as a single command.

In this document, to *type* text means to type the required text in the specified location, and to *enter* text means to type the required text, then press the ENTER key.

References to literature sources are given in square brackets, like this: [10].

MacroModel Overview

1.1 MacroModel User Manual

This manual contains an introduction to the MacroModel molecular mechanics program. For detailed information about command line MacroModel and MacroModel operation codes, see the *MacroModel Reference Manual*. For tutorial exercises that demonstrate Maestro and MacroModel functionality, see the *MacroModel Quick Start Guide*. For detailed information on using Maestro, the graphical user interface (GUI) for MacroModel, see the Maestro online help or the *Maestro User Manual*.

[Chapter 1](#) provides a brief introduction to MacroModel and describes how it interacts with Maestro, while [Chapter 2](#) provides a brief overview of how to use Maestro. A good general background to force field based molecular modeling, particularly as it is implemented in MacroModel, is given in [Chapter 3](#). Running MacroModel from the command-line (without Maestro) is described in [Chapter 4](#). Chapters [5-17](#) provide more information on how to run particular classes of MacroModel calculations. [Chapter 20](#) describes additional features: geometry queries and interaction energy calculations. Finally, [Chapter 21](#) describes how to get additional help for performing MacroModel calculations.

1.2 MacroModel

MacroModel is a general purpose, force-field-based molecular modeling program with applicability to a wide range of chemical systems. MacroModel provides researchers with multiple advanced methods to aid the understanding of chemical structure, energetics, and dynamics. A large selection of force fields is available in MacroModel, including the latest technical advances introduced into OPLS_2005, a force field that Schrödinger is actively developing. Numerous minimization methods are available, enabling geometry optimizations for a broad selection of structural classes. A wide range of methods is available for conformational searching, allowing you to efficiently sample the potential energy surface for low-energy structures, for systems ranging from small molecules to entire proteins. Solvation effects can be accounted for using the efficient continuum solvation models employed by MacroModel. Additional advanced features include molecular dynamics simulations, free energy perturbation simulations, and pure and mixed methods for ensemble sampling. MacroModel 9.1 contains several new features and performance enhancements, which reflect a commitment to provide the latest advancements in computational science.

1.3 MacroModel and Maestro Interaction

Maestro is the graphical user interface for MacroModel. MacroModel runs energy calculations as independent tasks and consequently does not tie up Maestro during lengthy computations. Maestro monitors MacroModel tasks so that both numerical and structural information may be viewed while the tasks are running. Such monitoring is the default mode of operation for newly started tasks, although monitoring can be broken off and reestablished at a later time. Thus, you can initiate several MacroModel tasks, disconnect from them, carry out graphical modeling operations, and periodically reconnect to and examine the progress of each of the previously submitted MacroModel tasks.

1.4 Calculation Preparation and Submission

MacroModel calculations can be prepared and launched from the Maestro GUI or from the command line. An overview of each job submission type is provided in this section.

From Maestro

Separate panels for current energy, energy minimization, dihedral driving, conformational search, ligand torsion search, multiple minimization, dynamics, MC/SD, and MINTA are available in Maestro. To set up a calculation, display the relevant panel and adjust the settings as desired, then click Start to set up and launch the job. Alternatively, you can choose to write out the structure and command files that you will need to launch the job later from the command line by clicking Write.

From the Command Line

Some types of MacroModel calculations cannot be prepared or submitted from Maestro. For these jobs, you must manually create the necessary command file, which provides MacroModel with the instructions it needs to perform the related calculation. For this task, you might find it easiest to use a Maestro-generated command file as a template.

Once you set up a job in Maestro and click either Start or Write, Maestro writes out a molecular structure file and a command file. These files are required to run MacroModel calculations. The molecular structure file contains the structures to be used as input in the calculation. The command file contains the name of the input structure file, the name of the output structure file, and an ordered list of operation codes (*opcodes*) for the calculations.

The *MacroModel Reference Manual* contains detailed information about all of the MacroModel commands and operational codes. Examples of command files for various standard modeling operations are provided at the end of relevant chapters in this manual. Because all

jobs are now handled by the Schrödinger Job Control facility, MacroModel jobs can be monitored from Maestro even when jobs are launched from the command line.

1.5 Command-Line Utilities

MacroModel is distributed with the following command-line utilities:

<code>autoref</code>	Performs a restrained minimization of a protein-ligand structure using MacroModel. See Section 20.6 on page 201 .
<code>para_bmin</code>	Provides an easy-to-use method of distributing serial MacroModel searches across multiple processors. See Section 4.3 on page 50 .
<code>premin</code>	Robustly minimizes ligand structures in a multistructure Maestro file using MacroModel. Problematic structures are separated out for examination. See Section 8.5 on page 84 .
<code>serial_split</code>	Splits the output structure file from a serial search into separate structure files for each individual search performed. See Section 20.7 on page 201 .

The following command-line utilities, which are installed with the software, may be useful in conjunction with MacroModel calculations:

<code>applyhtreat</code>	Adds or removes hydrogen atoms, dummy atoms, and lone pairs to structures in a MacroModel or Maestro file and produces a new Maestro file.
<code>maemmod, mmodmae</code>	Convert between Maestro and MacroModel file formats.
<code>mmodmol, molmmod</code>	Convert between MacroModel and Sybyl Mol2 file formats.
<code>maesubset</code>	Selects a subset of the structures from a multistructure Maestro file.
<code>propfilter</code>	Selects structures from a multistructure Maestro file based on properties.
<code>proplist</code>	Lists selected properties present in structures within a multistructure Maestro file.
<code>sdconvert</code>	Converts between SD and Maestro or MacroModel file formats.
<code>sdssubset</code>	Selects a subset of the structures from a multistructure SD file.

For more information on these utilities, see [Appendix D](#) of the *Maestro User Manual*.

1.6 Citing MacroModel in Publications

The use of this product should be acknowledged in publications as:

MacroModel, version 9.1, Schrödinger, LLC, New York, NY, 2005.

Introduction to Maestro

Maestro is the graphical user interface for all of Schrödinger's products: CombiGlide™, Epik™, Glide™, Impact™, Jaguar™, Liaison™, LigPrep™, MacroModel®, Phase™, Prime™, QikProp™, QSite™, and Strike™. It contains tools for building, displaying, and manipulating chemical structures; for organizing, loading, and storing these structures and associated data; and for setting up, monitoring, and visualizing the results of calculations on these structures. This chapter provides a brief introduction to Maestro and some of its capabilities. For more information on any of the topics in this chapter, see the [Maestro User Manual](#).

2.1 General Interface Behavior

Most Maestro panels are amodal: more than one panel can be open at a time, and a panel need not be closed for an action to be carried out. Each Maestro panel has a Close button so you can hide the panel from view.

Maestro supports the mouse functions common to many graphical user interfaces. The left button is used for choosing menu items, clicking buttons, and selecting objects by clicking or dragging. This button is also used for resizing and moving panels. The right button displays a shortcut menu. Other common mouse functions are supported, such as using the mouse in combination with the SHIFT or CTRL keys to select a range of items and select or deselect a single item without affecting other items.

In addition, the mouse buttons are used for special functions described later in this chapter. These functions assume that you have a three-button mouse. If you have a two-button mouse, ensure that it is configured for three-button mouse simulation (the middle mouse button is simulated by pressing or holding down both buttons simultaneously).

2.2 Starting Maestro

Before starting Maestro, you must first set the SCHRODINGER environment variable to point to the installation directory. To set this variable, enter the following command at a shell prompt:

```
csh/tcsh:      setenv SCHRODINGER installation-directory
bash/ksh:      export SCHRODINGER=installation-directory
```

You might also need to set the `DISPLAY` environment variable, if it is not set automatically when you log in. To determine if you need to set this variable, enter the command:

```
echo $DISPLAY
```

If the response is a blank line, set the variable by entering the following command:

```
csh/tcsh:      setenv DISPLAY display-machine-name:0.0
```

```
bash/ksh:      export DISPLAY=display-machine-name:0.0
```

After you set the `SCHRODINGER` and `DISPLAY` environment variables, you can start Maestro using the command:

```
$SCHRODINGER/maestro options
```

If you add the `$SCHRODINGER` directory to your path, you only need to enter the command `maestro`. Options for this command are given in [Section 2.1](#) of the *Maestro User Manual*.

The directory from which you started Maestro is Maestro's current working directory, and all data files are written to and read from this directory unless otherwise specified (see [Section 2.8 on page 27](#)). You can change directories by entering the following command in the command input area (see [page 8](#)) of the main window:

```
cd directory-name
```

where *directory-name* is either a full path or a relative path.

2.3 The Maestro Main Window

The Maestro main window is shown in [Figure 2.1 on page 7](#). The main window components are listed below.

The following components are always visible:

- **Title bar**—displays the Maestro version, the project name (if there is one) and the current working directory.
- **Auto-Help**—automatically displays context-sensitive help.
- **Menu bar**—provides access to panels.
- **Workspace**—displays molecular structures and other 3D graphical objects.

The following components can be displayed or hidden by choosing the component from the Display menu. Your choice of which main window components are displayed is persistent between Maestro sessions.




Figure 2.1. The Maestro main window.

- **Toolbar**—contains buttons for many common tasks and provides tools for displaying and manipulating structures, as well as organizing the Workspace.
- **Status bar**—displays information about a particular atom, or about structures in the Workspace, depending on where the pointer pauses (see [Section 2.5](#) of the *Maestro User Manual* for details):
 - **Atom**—displays the chain, residue number, element, PDB atom name, formal charge, and title or entry name (this last field is set by choosing Preferences from the Maestro menu and selecting the Feedback folder).

- **Workspace**—displays the number of atoms, entries, residues, chains, and molecules in the Workspace.
- **Clipping planes window**—displays a small, top view of the Workspace and shows the clipping planes and viewing volume indicators.
- **Sequence viewer**—shows the sequences for proteins displayed in the Workspace. See [Section 2.6](#) of the *Maestro User Manual* for details.
- **Command input area**—provides a place to enter Maestro commands.

When a distinction between components in the main window and those in other panels is needed, the term *main* is applied to the main window components (e.g., main toolbar).

You can expand the Workspace to occupy the full screen, by pressing CTRL+=. All other components and panels are hidden. To return to the previous display, press CTRL+= again.

2.3.1 The Menu Bar

The menus on the main menu bar provide access to panels, allow you to execute commands, and control the appearance of the Workspace. The main menus are as follows:

- **Maestro**—save or print images in the Workspace, execute system commands, save or load a panel layout, set preferences, set up Maestro command aliases, and quit Maestro.
- **Project**—open and close projects, import and export structures, make a snapshot, and annotate a project. These actions can also be performed from the Project Table panel. For more information, see [Section 2.4 on page 13](#).
- **Edit**—undo actions, build and modify structures, define command scripts and macros, and find atoms in the Workspace.
- **Display**—control the display of the contents of the Workspace, arrange panels, and display or hide main window components.
- **Tools**—group atoms; measure, align, and superimpose structures; and view and visualize data.
- **Applications**—set up, submit, and monitor jobs for Schrödinger’s computational programs. Some products have a submenu from which you can choose the task to be performed.
- **Scripts**—manage and install Python scripts that come with the distribution and scripts that you create yourself. (See [Chapter 13](#) of the *Maestro User Manual* for details.)
- **Help**—open the Help panel, the PDF documentation index, or information panels; run a demonstration; and display or hide Balloon Help (tooltips).

2.3.2 The Toolbar

The main toolbar contains three kinds of buttons for performing common tasks:



Action—Perform a simple task, like clearing the Workspace.



Display—Open or close a panel or open a dialog box, such as the Project Table panel.



Menu—Display a *button menu*. These buttons have a triangle in the lower right corner.

There are four types of items on button menus, and all four types can be on the same menu (see Figure 2.2):

- **Action**—Perform an action immediately.
- **Display**—Open a panel or dialog box.
- **Object types for selection**—Choose Atoms, Bonds, Residues, Chains, Molecules, or Entries, then click on an atom in the Workspace to perform the action on all the atoms in that structural unit.

The object type is marked on the menu with a red diamond and the button is indented to indicate the action to be performed.

- **Other setting**—Set a state, choose an attribute, or choose a parameter and click on atoms in the Workspace to display or change that parameter.

The toolbar buttons are described below. Some descriptions refer to features not described in this chapter. See the *Maestro User Manual* for a fuller description of these features.




Figure 2.2. The Workspace selection *button menu* and the Adjust distances, angles or dihedrals *button menu*.

Workspace selection

- Choose an object type for selecting
- Open the Atom Selection dialog box




Undo/Redo

Undo or redo the last action. Performs the same function as the Undo item on the Edit menu, and changes to an arrow pointing in the opposite direction when an Undo has been performed, indicating that its next action is Redo.

Open a project

Open the Open Project dialog box.



Import structures

Open the Import panel.

Open/Close Project Table

Open the Project Table panel or close it if it is open.



Save as

Open the Save Project As dialog box, to save the project with a new name.

Create entry from Workspace

Open a dialog box in which you can create an entry in the current project using the contents of the Workspace.




Delete

- Choose an object type for deletion
- Delete hydrogens and waters
- Open the Atom Selection dialog box
- Delete other items associated with the structures in the Workspace
- Click to select atoms to delete
- Double-click to delete all atoms

Open/Close Build panel

Open the Build panel or close it if it is open.



Add hydrogens

- Choose an object type for applying a hydrogen treatment
- Open the Atom Selection dialog box
- Click to select atoms to treat
- Double-click to apply to all atoms

Local transformation

- Choose an object type for transforming
- Click to select atoms to transform
- Open the Advanced Transformations panel



Adjust distances, angles or dihedrals

- Choose a parameter for adjusting
- Delete adjustments

Fit to screen

Scale the displayed structure to fit into the Workspace and reset the center of rotation.



Clear Workspace

Clear all atoms from the Workspace.

Set fog display state

Choose a fog state. Automatic means fog is on when there are more than 40 atoms in the Workspace, otherwise it is off.



Enhance depth cues

Optimize fogging and other depth cues based on what is in the Workspace.

Rotate around X axis by 90 degrees

Rotate the Workspace contents around the X axis by 90 degrees.



Rotate around Y axis by 90 degrees

Rotate the Workspace contents around the Y axis by 90 degrees.

Tile entries

Arrange entries in a rectangular grid in the Workspace.

**Save view**

Save the current view of the Workspace: orientation, location, and zoom.

**Display only selected atoms**

- Choose an object type for displaying
- Click to select atoms to display
- Double-click to display all atoms

**Also display**

- Choose a predefined atom category
- Open the Atom Selection dialog box

**Display residues within N angstroms of currently displayed atoms**

- Choose a radius
- Open a dialog box to set a value

**Draw bonds in wire**

- Choose an object type for drawing bonds in wire representation
- Open the Atom Selection dialog box
- Click to select atoms for representation
- Double-click to apply to all atoms

**Draw atoms in Ball & Stick**

- Choose an object type for drawing bonds in Ball & Stick representation
- Open the Atom Selection dialog box
- Click to select atoms for representation
- Double-click to apply to all atoms

**Color all atoms by scheme**


Choose a predefined color scheme.

**Label atoms**


- Choose a predefined label type
- Delete labels

**Reset Workspace**


Reset the rotation, translation, and zoom of the Workspace to the default state.

**Restore view**


Restore the last saved view of the Workspace: orientation, location, and zoom.

**Display only**


- Choose a predefined atom category
- Open the Atom Selection dialog box

**Undisplay**


- Choose a predefined atom category
- Open the Atom Selection dialog box

**Show, hide, or color ribbons**


- Choose to show or hide ribbons
- Choose a color scheme for coloring ribbons

**Draw atoms in CPK**


- Choose an object type for drawing bonds in CPK representation
- Open the Atom Selection dialog box
- Click to select atoms for representation
- Double-click to apply to all atoms

**Draw bonds in tube**


- Choose an object type for drawing bonds in tube representation
- Open the Atom Selection dialog box
- Click to select atoms for representation
- Double-click to apply to all atoms

**Color residue by constant color**

- Choose a color for applying to residues
- Click to select residues to color
- Double-click to color all atoms

**Label picked atoms**

- Choose an object type for labeling atoms
- Open the Atom Selection dialog box
- Open the Atom Labels panel at the Composition folder
- Delete labels
- Click to select atoms to label
- Double-click to label all atoms



Display H-bonds

- Choose bond type:
intra—displays H-bonds within the selected molecule
- inter—displays H-bonds between the selected molecule and all other atoms.
- Delete H-bonds
- Click to select molecule



Measure distances, angles or dihedrals

- Choose a parameter for displaying measurements
- Delete measurements
- Click to select atoms for measurement

2.3.3 Mouse Functions in the Workspace

The left mouse button is used for selecting objects. You can either click on a single atom or bond, or you can drag to select multiple objects. The right mouse button opens shortcut menus, which are described in [Section 2.7](#) of the *Maestro User Manual*.

The middle and right mouse buttons can be used on their own and in combination with the SHIFT and CTRL keys to perform common operations, such as rotating, translating, centering, adjusting, and zooming.

Table 2.1. Mapping of Workspace operations to mouse actions.

Mouse Button	Keyboard	Motion	Action
Left		click, drag	Select
Left	SHIFT	click, drag	Toggle the selection
Middle		drag	Rotate about X and Y axes Adjust bond, angle, or dihedral
Middle	SHIFT	drag vertically	Rotate about X axis
Middle	SHIFT	drag horizontally	Rotate about Y axis
Middle	CTRL	drag horizontally	Rotate about Z axis
Middle	SHIFT + CTRL	drag horizontally	Zoom
Right		click	Spot-center on selection
Right		click and hold	Display shortcut menu
Right		drag	Translate in the X-Y plane
Right	SHIFT	drag vertically	Translate along the X axis
Right	SHIFT	drag horizontally	Translate along the Y axis
Right	CTRL	drag horizontally	Translate along the Z axis
Middle & Right		drag horizontally	Zoom

2.3.4 Shortcut Key Combinations

Some frequently used operations have been assigned shortcut key combinations. The shortcuts available in the main window are described in [Table 2.2](#).

Table 2.2. Shortcut keys in the Maestro main window.

Keys	Action	Equivalent Menu Choices
CTRL+B	Open Build panel	Edit > Build
CTRL+C	Create entry	Project > Create Entry From Workspace
CTRL+E	Open Command Script Editor panel	Edit > Command Script Editor
CTRL+F	Open Find Atoms panel	Edit > Find
CTRL+H	Open Help panel	Help > Help
CTRL+I	Open Import panel	Project > Import Structures
CTRL+M	Open Measurements panel	Tools > Measurements
CTRL+N	Create new project	Project > New
CTRL+O	Open project	Project > Open
CTRL+P	Print	Maestro > Print
CTRL+Q	Quit	Maestro > Quit
CTRL+S	Open Sets panel	Tools > Sets
CTRL+T	Open Project Table panel	Project > Show Table
CTRL+W	Close project	Project > Close
CTRL+Z	Undo/Redo last command	Edit > Undo/Redo
CTRL+=	Enter and exit full screen mode (Workspace occupies full screen)	None

2.4 Maestro Projects

All the work you do in Maestro is done within a *project*. A project consists of a set of *entries*, each of which contains one or more chemical structures and their associated data. In any Maestro session, there can be only one Maestro project open. If you do not specify a project when you start Maestro, a *scratch* project is created. You can work in a scratch project without saving it, but you must save it in order to use it in future sessions. When you save or close a project, all the view transformations (rotation, translation, and zoom) are saved with it. When you close a project, a new scratch project is automatically created.

Likewise, if there is no entry displayed in the Workspace, Maestro creates a *scratch* entry. Structures that you build in the Workspace constitute a scratch entry until you save the structures as project entries. The scratch entry is not saved with the project unless you explicitly add it to the project. However, you can use a scratch entry as input for some calculations.

To add a scratch entry to a project, do one of the following:

- Click the Create entry from Workspace button:



- Choose Create Entry from Workspace from the Project menu.
- Press CTRL+C.

In the dialog box, enter a name and a title for the entry. The entry name is used internally to identify the entry and can be modified by Maestro. The title can be set or changed by the user, but is not otherwise modified by Maestro.

Once an entry has been incorporated into the project, its structures and their data are represented by a row in the Project Table. Each row contains the row number, an icon indicating whether the entry is displayed in the Workspace (the In column), the entry title, a button to open the Surfaces panel if the entry has surfaces, the entry name, and any entry properties. The row number is not a property of the entry.

Entries can be collected into groups, and the members of the group can be displayed or hidden. Most additions of multiple entries to the Project Table are done as entry groups.

You can use entries as input for all of the computational programs—Glide, Impact, Jaguar, Liaison, LigPrep, MacroModel, Phase, Prime, QikProp, QSite, and Strike. You can select entries as input for the ePlayer, which displays the selected structures in sequence. You can also duplicate, combine, rename, and sort entries; create properties; import structures as entries; and export structures and properties from entries in various formats.

To open the Project Table panel, do one of the following:

- Click the Open/Close Project Table button on the toolbar



- Choose Show Table from the Project menu
- Press CTRL+T.

The Project Table panel contains a menu bar, a toolbar, and the table itself.




Figure 2.3. The Project Table panel.

2.4.1 The Project Table Toolbar

The Project Table toolbar contains two groups of buttons and a status display. The first set of buttons opens various panels that allow you to perform functions on the entries in the Project Table. The second set of buttons controls the ePlayer, which “plays through” the selected structures: each structure is displayed in the Workspace in sequence, at a given time interval. See [Section 2.3.2 on page 9](#) for a description of the types of toolbar buttons. The buttons are described below.



Find

Open the Find panel for locating alphanumeric text in any column of the Project Table, except for the row number.



Sort

Open the Sort panel for sorting entries by up to three properties.



Plot

Open the Plot panel for plotting entry properties.



Import Structure

Open the Import panel for importing structures into the project.



Export Structure

Open the Export panel for exporting structures to a file.



Columns

Choose an option for adjusting the column widths.



Select only

Open the Entry Selection dialog box for selecting entries based on criteria for entry properties.



Go to start

Display the first selected structure.



Previous

Display the previous structure in the list of selected structures.



Play backward

Display the selected structures in sequence, moving toward the first.



Stop

Stop the ePlayer.



Play forward

Display the selected structures in sequence, moving toward the last.



Next

Display the next structure in the list of selected structures.



Go to end

Display the last selected structure.



Loop

Choose an option for repeating the display of the structures. **Single Direction** displays structures in a single direction, then repeats. **Oscillate** reverses direction each time the beginning or end of the list is reached.

The status display, to the right of the toolbar buttons, shows the number of selected entries. When you pause the cursor over the status display, the Balloon Help shows the total number of entries, the number shown in the table, the number selected, and the number included in the Workspace.

2.4.2 The Project Table Menus

- **Table**—find text, sort entries, plot properties, import and export structures, and configure the Project Table.
- **Select**—select all entries, none, invert your selection, or select classes of entries using the Entry Selection dialog box and the Filter panel.

- Entry—include or exclude entries from the Workspace, display or hide entries in the Project Table, and perform various operations on the selected entries.
- Property—display and manipulate entry properties in the Project Table.
- ePlayer—view entries in succession, stop, reverse, and set the ePlayer options.

2.4.3 Selecting Entries

Many operations in Maestro are performed on the entries selected in the Project Table. The Project Table functions much like any other table: select rows by clicking, shift-clicking, and control-clicking. However, because clicking in an editable cell of a selected row enters edit mode, you should click in the Row column to select entries. See [Section 2.4.5 on page 18](#) for more information on mouse actions in the Project Table. There are shortcuts for selecting classes of entries on the Select menu.

In addition to selecting entries manually, you can select entries that meet a combination of conditions on their properties. Such combinations of conditions are called *filters*. Filters are Entry Selection Language (ESL) expressions and are evaluated at the time they are applied. For example, if you want to set up a Glide job that uses ligands with a low molecular weight (say, less than 300) and that has certain QikProp properties, you can set up a filter and use it to select entries for the job. If you save the filter, you can use it again on a different set of ligands that meet the same selection criteria.

To create a filter:

1. Do one of the following:
 - Choose Only, Add, or Deselect from the Select menu.
 - Click the Entry selection button on the toolbar.



2. In the Properties folder, select a property from the property list, then select a condition.
3. Combine this selection with the current filter by clicking Add, Subtract, or Intersect. These buttons perform the Boolean operations OR, AND NOT, and AND on the corresponding ESL expressions.
4. To save the filter for future use click Create Filter, enter a name, and click OK.
5. Click OK to apply the filter immediately.

2.4.4 Including Entries in the Workspace

In addition to selecting entries, you can also use the Project Table to control which entries are displayed in the Workspace. An entry that is displayed in the Workspace is *included* in the Workspace; likewise, an entry that is not displayed is *excluded*. Included entries are marked by an X in the diamond in the In column; excluded entries are marked by an empty diamond. Entry inclusion is completely independent of entry selection.

To include or exclude entries, click, shift-click, or control-click in the In column of the entries, or select entries and choose Include or Exclude from the Entry menu. Inclusion with the mouse works just like selection: when you include an entry by clicking, all other entries are excluded.

It is sometimes useful to keep one entry in the Workspace and include others one by one: for example, a receptor and a set of ligands. You can fix the receptor in the Workspace by selecting it in the Project Table and choosing Fix from the Entry menu or by pressing CTRL+F. A padlock icon replaces the diamond in the In column to denote a *fixed* entry. To remove a fixed entry from the Workspace, you must exclude it explicitly (CTRL+X). It is not affected by the inclusion or exclusion of other entries. Fixing an entry affects only its inclusion; you can still rotate, translate, or modify the structure.

2.4.5 Mouse Functions in the Project Table

The Project Table supports the standard use of shift-click and control-click to select objects. This behavior applies to the selection of entries and the inclusion of entries in the Workspace. You can also drag to resize rows and columns and to move rows.

You can drag a set of non-contiguous entries to reposition them in the Project Table. When you release the mouse button, the entries are placed after the first unselected entry that precedes the entry on which the cursor is resting. For example, if you select entries 2, 4, and 6, and release the mouse button on entry 3, these three entries are placed after entry 1, because entry 1 is the first unselected entry that precedes entry 3. To move entries to the top of the table, drag them above the top of the table; to move entries to the end of the table, drag them below the end of the table.

A summary of mouse functions in the Project Table is provided in [Table 2.3](#).

Table 2.3. Mouse operations in the Project Table.

Task	Mouse Operation
Change a Boolean property value	Click repeatedly in a cell to cycle through the possible values (On, Off, Clear)
Display the Entry menu for an entry	Right-click anywhere in the entry. If the entry is not selected, it becomes the selected entry. If the entry is selected, the action is applied to all selected entries.
Display a version of the Property menu for a property	Right-click in the column header
Edit the text or the value in a table cell	Click in the cell and edit the text or value
Include an entry in the Workspace, exclude all others	Click the In column of the entry
Move selected entries	Drag the entries
Paste text into a table cell	Middle-click
Resize rows or columns	Drag the boundary with the middle mouse button
Select an entry, deselect all others	For an unselected entry, click anywhere in the row except the In column; for a selected entry, click the row number.
Select or include multiple entries	Click the first entry then shift-click the last entry
Toggle the selection or inclusion state	Control-click the entry or the In column

2.4.6 Project Table Shortcut Keys

Some frequently used project operations have been assigned shortcut key combinations. The shortcuts, their functions, and their menu equivalents are listed in [Table 2.4](#).

Table 2.4. Shortcut keys in the Project Table.

Keys	Action	Equivalent Menu Choices
CTRL+A	Select all entries	Select > All
CTRL+F	Fix entry in Workspace	Entry > Fix
CTRL+I	Open Import panel	Table > Import Structures
CTRL+N	Include only selected entries	Entry > Include Only
CTRL+U	Deselect all entries	Select > None
CTRL+X	Exclude selected entries	Entry > Exclude
CTRL+Z	Undo/Redo last command	Edit > Undo/Redo in main window

2.5 Building a Structure

After you start Maestro, the first task is usually to create or import a structure. You can open existing Maestro projects or import structures from other sources to obtain a structure, or you can build your own. To open the Build panel, do one of the following:

- Click the Open/Close Build panel button in the toolbar:



- Choose Build from the Edit menu.
- Press CTRL+B.

The Build panel allows you to create structures by drawing or placing atoms or fragments in the Workspace and connecting them into a larger structure, to adjust atom positions and bond orders, and to change atom properties. This panel contains a toolbar and three folders.

2.5.1 Placing and Connecting Fragments

The Build panel provides several tools for creating structures in the Workspace. You can place and connect fragments, or you can draw a structure freehand.

To place a fragment in the Workspace:

1. Select Place.
2. Choose a fragment library from the Fragments menu.
3. Click a fragment.
4. Click in the Workspace where you want the fragment to be placed.

To connect fragments in the Workspace, do one of the following:

- Place another fragment and connect them using the Connect & Fuse panel, which you open from the Edit menu on the main menu bar or with the Display Connect & Fuse panel on the Build toolbar.



- Replace one or more atoms in the existing fragment with another fragment by selecting a fragment and clicking in the Workspace on the main atom to be replaced.
- Grow another fragment by selecting Grow in the Build panel and clicking the fragment you want to add in the Fragments folder.

2.5.2 Adjusting Properties

In the Atom Properties folder, you can change the properties of the atoms in the Workspace. For each item on the Property option menu—Element, Atom Type (MacroModel), Partial Charge, PDB Atom Name, Grow Name, and Atom Name—there is a set of tools you can use to change the atom properties. For example, the Element tools consist of a periodic table from which you can choose an element and select an atom to change it to an atom of the selected element.

Similarly, the Residue Properties folder provides tools for changing the properties of residues: the Residue Number, the Residue Name, and the Chain Name.

To adjust bond lengths, bond angles, dihedral angles, and chiralities during or after building a structure, use the Adjust distances, angles or dihedrals button on the main toolbar:



You can also open the Adjust panel from this button menu, from the Display Adjust panel button on the Build panel toolbar (which has the same appearance as the above button) or from the Edit menu in the main window.

2.5.3 The Build Panel Toolbar

The toolbar of the Build panel provides quick access to tools for drawing and modifying structures and labeling atoms. See [Section 2.3.2 on page 9](#) for a description of the types of toolbar buttons. The toolbar buttons and their use are described below.



Free-hand drawing

Choose an element for drawing structures freehand in the Workspace (default C). Each click in the Workspace places an atom and connects it to the previous atom.



Delete

Choose an object for deleting. Same as the [Delete](#) button on the main toolbar, see [page 10](#).



Set element

Choose an element for changing atoms in the Workspace (default C). Click an atom to change it to the selected element.



Increment bond order

Select a bond to increase its bond order by one, to a maximum of 3.



Decrement bond order

Select a bond to decrease its bond order by one, to a minimum of 0.

**Increment formal charge**

Select an atom to increase its formal charge by one.

**Decrement formal charge**

Select an atom to decrease its formal charge by one.

**Move**

Choose a direction for moving atoms, then click the atom to be moved. Moves in the XY plane are made by clicking the new location. Moves in the Z direction are made in 0.5 Å increments.

**Label**

Apply heteroatom labels as you build a structure. The label consists of the element name and formal charge, and is applied to atoms other than C and H.

**Display Connect & Fuse panel**

Open the Connect & Fuse panel so you can connect structures (create bonds between structures) or fuse structures (replace atoms of one structure with those of another).

**Display Adjust panel**

Open the Adjust panel so you can change bond lengths, bond angles, dihedral angles, or atom chiralities.

**Add hydrogens**

Choose an atom type for applying the current hydrogen treatment. Same as the [Add hydrogens](#) button on the main toolbar, see [page 10](#).

**Geometry Symmetrizer**

Open the Geometry Symmetrizer panel for symmetrizing the geometry of the structure in the Workspace.

**Geometry Cleanup**

Clean up the geometry of the structure in the Workspace.

2.6 Selecting Atoms

Maestro has a powerful set of tools for selecting atoms in a structure: toolbar buttons, picking tools in panels, and the Atom Selection dialog box. These tools allow you to select atoms in two ways:

- Select atoms first and apply an action to them
- Choose an action first and then select atoms for that action

2.6.1 Toolbar Buttons

The small triangle in the lower right corner of a toolbar button indicates that the button contains a menu. Many of these buttons allow you to choose an object type for selecting: choose Atoms, Bonds, Residues, Chains, Molecules, or Entries, then click on an atom in the Workspace to perform the action on all the atoms in that structural unit.

For example, to select atoms with the Workspace selection toolbar button:

1. Choose Residues from the Workspace selection button menu:



The button changes to:



2. Click on an atom in a residue in the Workspace to select all the atoms in that residue.

2.6.2 Picking Tools

The picking tools are embedded in each panel in which you need to select atoms to apply an operation. The picking tools in a panel can include one or more of the following:

- Pick option menu—Allows you to choose an object type. Depending on the operation to be performed, you can choose Atoms, Bonds, Residues, Chains, Molecules, or Entries, then click on an atom in the Workspace to perform the action on all the atoms in that structural unit.

The Pick option menu varies from panel to panel, because not all object types are appropriate for a given operation. For example, some panels have only Atoms and Bonds in the Pick option menu.

- All button—Performs the action on all atoms in the Workspace.
- Selection button—Performs the action on any atoms already selected in the Workspace.
- Previous button—Performs the action on the most recent atom selection defined in the Atom Selection dialog box.
- Select button—Opens the Atom Selection dialog box.
- ASL text box—Allows you to type in an ASL expression for selecting atoms.

ASL stands for Atom Specification Language, and is described in detail in the [Maestro Command Reference Manual](#).

- Clear button—Clears the current selection



- Show markers option—Marks the selected atoms in the Workspace.

For example, to label atoms with the Label Atoms panel:

1. Choose Atom Labels from the Display menu.
2. In the Composition folder, select Element and Atom Number.
3. In the picking tools section at the top of the panel, you could do one of the following:
 - Click Selection to apply labels to the atoms already selected in the Workspace (from the previous example).
 - Choose Residues from the Pick option menu and click on an atom in a different residue to label all the atoms in that residue.

2.6.3 The Atom Selection Dialog Box

If you wish to select atoms based on more complex criteria, you can use the Atom Selection dialog box. To open this dialog box, choose Select from a button menu or click the Select button in a panel. See [Section 5.3](#) of the *Maestro User Manual* for detailed instructions on how to use the Atom Selection dialog box.

2.7 Scripting in Maestro

Although you can perform nearly all Maestro-supported operations through menus and panels, you can also perform operations using Maestro commands, or compilations of these commands, called *scripts*. Scripts can be used to automate lengthy procedures or repetitive tasks and can be created in several ways. These are summarized below.

2.7.1 Python Scripts

Python is a full-featured scripting language that has been embedded in Maestro to extend its scripting facilities. The Python capabilities within Maestro include access to Maestro functionality for dealing with chemical structures, projects, and Maestro files.

The two main Python commands used in Maestro are:

- `pythonrun`—executes a Python module. (You can also use the alias `pyrun`.) The syntax is:

```
pythonrun module.function
```

- `pythonimport`—rereads a Python file so that the next time you use the `pythonrun` command, it uses the updated version of the module. (You can also use the alias `pyimp`.)

From the Maestro Scripts menu you can install, manage, and run Python scripts. For more information on the Scripts menu, see [Section 13.1](#) of the *Maestro User Manual*.

For more information on using Python with Maestro, see *Maestro Scripting with Python*.

2.7.2 Command Scripts

All Maestro commands are logged and displayed in the Command Script Editor panel. This means you can create a command script by performing the operations with the GUI controls, copying the logged commands from the Command History list into the Script text area of the panel, then saving the list of copied commands as a script.

To run an existing command script:


1. Open the Command Script Editor panel from the Edit menu in the main window.
2. Click Open Local and navigate to the directory containing the desired script.
3. Select a script in the Files list and click Open.

The script is loaded into the Script window of the Command Script Editor panel.

4. Click Run Script.

Command scripts cannot be used for Prime operations.

The *Command History* window displays a log of all commands issued internally within Maestro when you interact with a panel, menu, or structure



Opens the *Show/Hide Command* panel, used to determine which commands are logged in the *Command History* list

Figure 2.5. The Command Script Editor panel.

2.7.3 Macros

There are two kinds of macros you can create: named macros and macros assigned to function keys F1 through F12.

To create and run a named macro:

1. Open the Macros panel from the Edit menu in the main window.
2. Click New, enter a name for the macro, and click OK.
3. In the Definition text box, type the commands for the macro.
4. Click Update to update the macro definition.
5. To run the macro, enter the following in the command input area in the main window:

```
macrorun macro-name
```

If the command input area is not visible, choose Command Input Area from the Display menu.

To create and run a function key macro:

1. Open the Function Key Macros panel from the Edit menu in the main window.
2. From the Macro Key option, select a function key (F1 through F12) to which to assign the macro.
3. In the text box, type the commands for the macro.
4. Click Run to test the macro or click Save to save it.
5. To run the macro from the main window, press the assigned function key.

For more information on macros, see [Section 13.5](#) of the *Maestro User Manual*.

2.8 Specifying a Maestro Working Directory

When you use Maestro to launch MacroModel jobs, Maestro writes job output to the directory specified in the Directory folder of the Preferences panel. By default, this directory (the file I/O directory) is the directory from which you started Maestro.

To change the Maestro working directory:

1. Open the Preferences panel from the Maestro menu.
2. Click the Directory tab.
3. Select the directory you want to use for reading and writing files.




Figure 2.6. T

You can also set other preferences in the Preferences panel. See [Section 12.2](#) of the *Maestro User Manual* for details.

2.9 Undoing an Operation

To undo a single operation, click the Undo button in the toolbar, choose Undo from the Edit menu, or press CTRL+Z. The word Undo in the menu is followed by text that describes the operation to undo. Not all operations can be undone: for example, global rotations and translations are not undoable operations. For such operations you can use the Save view and Restore view buttons in the toolbar, which save and restore a molecular orientation.

2.10 Running and Monitoring Jobs

Maestro has panels for each product for preparing and submitting jobs. To use these panels, choose the appropriate product and task from the Applications menu and its submenus. Set the appropriate options in the panel, then click Start to open the Start dialog box and set options for running the job. For a complete description of the Start dialog box associated with your computational program, see your product's User Manual. When you have finished setting the options, click Start to launch the job and open the Monitor panel.

The Monitor panel is the control panel for monitoring the progress of jobs and for pausing, resuming, or killing jobs. All jobs that belong to your user ID can be displayed in the Monitor panel, whether or not they were started from Maestro. Subjobs are indented under their parent in the job list. The text pane shows various output information from the monitored job, such as the contents of the log file. The Monitor panel opens automatically when you start a job. If it is




Figure 2.7. The Monitor panel.

not open, you can open it by choosing Monitor from the Applications menu in the Maestro main window.

While jobs are running, the Detach, Pause, Resume, Stop, Kill, and Update buttons are active. When there are no jobs currently running, only the Monitor and Delete buttons are active. These buttons act on the selected job. By default, only jobs started from the current project are shown. To show other jobs, deselect Show jobs from current project only.

When a monitored job ends, the results are incorporated into the project according to the settings used to launch the job. If a job that is not currently being monitored ends, you can select it in the Monitor panel and click Monitor to incorporate the results. Monitored jobs are incorporated only if they are part of the current project. You can monitor jobs that are not part of the current project, but their results are not incorporated. To add their results to a project, you must open the project and import the results.

Further information on job control, including configuring your site, monitoring jobs, running jobs, and job incorporation, can be found in the [Job Control Guide](#) and the [Installation Guide](#).

2.11 Getting Help

Maestro comes with automatic, context-sensitive help (Auto-Help), Balloon Help (tooltips), an online help facility, and a user manual. To get help, follow the steps below:

- Check the Auto-Help text box at the bottom of the main window. If help is available for the task you are performing, it is automatically displayed there. It describes what actions are needed to perform the task.
- If your question concerns a GUI element, such as a button or option, there may be Balloon Help for the item. Pause the cursor over the element. If the Balloon Help does not appear, check that Show Balloon Help is selected in the Help menu of the main window. If there is Balloon Help for the element, it appears within a few seconds.
- If you do not find the help you need using either of the steps above, click the Help button in the lower right corner of the appropriate panel. The Help panel is displayed with a relevant help topic.
- For help with a concept or action not associated with a panel, open the Help panel from the Help menu or press CTRL+H.

If you do not find the information you need in the Maestro help system, check the following sources:

- The *Maestro User Manual*
- The Frequently Asked Questions page, found at <http://www.schrodinger.com/Support/faq.html>

You can also contact Schrödinger by e-mail or phone for help:

- E-mail: help@schrodinger.com
- Phone: (503) 299-1150

2.12 Ending a Maestro Session

To end a Maestro session, choose Quit from the Maestro menu. To save a log file with a record of all operations performed in the current session, click Quit, save log file in the Quit panel. This information can be useful to Schrödinger support staff when responding to any problem you report.

Basic Molecular Modeling

The potential energy model used for MacroModel energy calculations is the classical one known as molecular mechanics. It is an empirical model which is parameterized to reproduce known data from experiment or quantum mechanical calculations. The equation and parameter sets which allow calculation of energy from a molecular geometry are known as force fields.

Generally, MacroModel uses equations and parameters from published standard force fields. However, the MacroModel implementations differ in various ways from the authentic force fields. We distinguish the MacroModel implementation from the original force fields by adding a “*” to the end of the force field name. MacroModel force field parameters and equations selectors are found in force field files having suffixes .fld (e.g., mm2.fld). Differences between the MacroModel force fields and the standard fields are summarized below.

3.1 MacroModel Force Field Implementation

All MacroModel force field files contain the authentic parameter set published by the original authors of the force field. In addition to these parameters are other parameters from other sources (e.g., the literature or work at Columbia). Parameters in the force field files are labeled as to their origin (O = original from the force field authors, M = modified from the original values, and A = added from some other source where no original parameter exists). They are also labeled by quality (1 = high quality final value, 2 = tentative value based on more than one experimental or quantum calculation, 3 = crude low quality parameter). Sources of A and M parameters are given at the ends of the lines in the force field files and recent additions to the force fields are documented in the *MacroModel Technical Manual*.

The MacroModel implementations of standard force fields differ from the authentic force fields in the following ways:

MM2*—All force field equations are identical with those of authentic MM2 from Allinger [1], with the exception of the following:

- The electrostatic equation (MM2* uses partial charges and Coulomb’s law, whereas MM2 uses bond dipoles and the Jeans equation).
- The out-of-plane bending equation (MM2* uses an improper torsion while MM2 uses a pyramidalization distance—the difference being insignificant except for substantially distorted sp^2 systems).

- Handling of conjugation (MM2* uses specific V1-V3 torsional terms for various conjugated systems, whereas MM2 uses an SCF π calculation in uncommon systems).

MM3*—All force field equations are identical to those of authentic MM3 from Allinger [2], except for those differences listed above for MM2.

AMBER*—All force field equations are identical to those of authentic AMBER from Kollman. The MacroModel default for hydrogen bonding uses Kollman's recent 6,12-Lennard Jones treatment [3] and an improved peptide backbone parameter set [4].

OPLS*—All force field equations are identical to those of OPLS/AMBER from Jorgensen [5].

OPLS_2001—Also referred to as OPLSAA, this force field, developed by Professor W. Jorgensen of Yale University, is probably the best available for condensed-phase simulations of peptides. All force-field equations are identical to those of authentic OPLSAA [6]. MacroModel's implementation has been validated by comparison to BOSS OPLSAA calculations for a wide variety of organic systems. Comparisons to ab initio calculations and experiment show that OPLS_2001 reproduces conformational energies well for systems for which it has been specifically parameterized. Especially good results can be expected for proteins. With the exception of improved charge, van der Waals and torsion parameters for sulfur in thiols and thiol ethers [7], all parameters are native OPLS_2001. The new thio parameters, which use appreciably smaller charges on sulfur and which have been validated in liquid-phase simulations on thiols and thiol ethers, significantly improve the conformational energetics of CYS and MET residues in proteins.

OPLS_2005—OPLS_2005 is an enhanced version of the OPLSAA all atom force field developed by Schrödinger to provide a larger coverage of organic functionality. In particular all torsional parameters have been refit to reproduce the conformational energetics derived at a higher level of quantum theory and additional charges have been fit to support additional organic functionality. The parameters for proteins have been updated to the ones published more recently (Kaminski, G. A.; Friesner, R. A.; Tirado-Rives, J.; Jorgensen, W. J. J. *Phys. Chem. B* 2001, 105, 6474).

AMBER94—All force field equations and parameters are the same as in Cornell et al. [8], with the following small exceptions:

- In MacroModel partial charges are specified by bond dipoles rather than as charge values. The partial charges may differ slightly between the two implementations; these differences are typically in the fifth significant figure.
- The atoms defining improper torsions are not specified by the AMBER protocol in situations of high local symmetry. This may sometimes give rise to small differences in molecular energies or geometries between the two programs.

- The paper gives the two nitrogen types different van der Waals parameters, but the AMBER 4.1 program uses the same parameters for both. We follow the program's convention.

MMFF—Our implementation is identical to that described by Halgren [9–15]. We supply both MMFF94 and MMFF94s; the latter enforces planarity about delocalized sp^2 nitrogens.

OPLS_2001 and OPLS_2005 have much in common. In this document we refer to them collectively as OPLSAA.

Whenever a current energy calculation (ECalc) is carried out, a listing file (*jobname.mmo*) can be produced which contains all parameters used in the calculation along with the origin and quality of each parameter. Note that any torsion parameter where V_1 - V_3 are all set to zero will not be included in the output. To include these in the listing it is necessary to include the DEBG 56 command in the command file: see the *MacroModel Reference Manual* for details. This automatic parameter referencing feature provides important information on the quality and reliability of the calculation.

Different force fields use different defaults for their electrostatic treatment (constant or distance-dependent dielectric) and their cutoff distances (van der Waals and electrostatic). It is possible to set such options exactly as in the authentic fields using the Electrostatic Treatment and Cutoff option menus in the Potential folder within Maestro's MacroModel panels.

3.2 Parameter Quality Considerations

Molecular mechanics force fields are empirical, and the accuracy of the results rest entirely on the ability of the parameters and functional forms used to mimic the real potential energy of molecules. If the parameters are deficient, is it impossible to obtain agreement with experiment or to make useful predictions.

When MacroModel performs an energy calculation, the program checks the quality of each parameter in use. Use of low quality (quality = 3) parameters, especially torsional ones, may result in inaccurate conformational energy differences and geometries. Low quality stretches often indicate crude partial charges since charge information often originates from bond dipoles. Using such low quality parameter values can cause charges and solvation energies to be inaccurate. Consequently, whenever MacroModel initiates an energy calculation, a warning and the numbers of low quality stretch, bend, and torsional parameters in use are listed in the MacroModel *jobname.log* file and in the Maestro Monitor panel viewing window. An example is shown below:

```
WARNING - Conformational Energies May Not Be Accurate
WARNING - Solvation Energies/Charges May Not Be Accurate
Low quality force field parameters in use:
```

Number of low quality stretches, bends & torsions = 1 1 8

The above message indicates that eight low quality torsions, one low quality bend, and one low quality stretch are in use in the calculation. Consequently, conformational energy differences and solvation energies may be unreliable. By looking in the job's *jobname.log* file, you can see which line(s) in the force field file is the source of the low quality parameters. Listings of the specific torsions having low quality, type 3 parameters can be found in the job's *.mmio* file after running the job with a "Complete" energy listing.

If the file notes that low quality parameters are being used, first attempt to find improved parameters by trying the system of interest with each MacroModel force field. MacroModel supports a number of different force fields, and each has its own strengths. If you are unable to find a force field with acceptable parameters for your particular system, you should seriously consider developing your own parameters, or try to obtain parameters from other practitioners in the field.

Typically, low quality torsional parameters are generalized ones (e.g., any bond between two sp^3 carbons has a 3-fold rotational barrier of ~ 3 kcal/mol). Thus for a specific torsion involving, for example $Cl-C(sp^3)-C(sp^3)-O$, a new torsional parameter line would be added to the force field file with specific V1-V3 terms which reproduce reliable experimental data or high quality ab initio calculations (e.g., 6-31G*) on, for example, β -chloroethyl methyl ether. To make the calculation more reliable, new parameters that fit experimental or high quality ab initio data need to be determined and added to the force field file. However, if you are comparing conformations that do not involve significant changes in the torsions having low quality parameters, then errors may not be large.

Obtaining correct torsional profiles around rotatable bonds is one of the most commonly encountered problems in reparameterizing a potential energy surface. For example, if a new charge set is introduced, the non-bonded interactions on either side of a rotatable bond may alter the torsional profile. It is important to test that the relative energies of the minima and maxima are in agreement with the results of high-level molecular-orbital calculations on model compounds. In order to make this sort of comparison, the molecular mechanics calculations should be performed in the gas phase; that is, with the solvation model turned off and with constant-dielectric electrostatics. Once the potential-energy surface is correctly described in the gas phase, the solvation model should correctly describe the potential-energy surface in solution.

To determine the force-field quality for the problem at hand, check the summary of low quality parameters in use by using Maestro's Force Field Viewer panel, which is opened from the Analysis menu. Additional information on developing parameters can be found in the *MacroModel Technical Manual*.

3.3 MacroModel Solvation Treatment

While many molecular modeling studies are carried out without including the effect of solvent, the omission is largely one of expedience. Most experimental studies are carried out in solvent, and the solvent medium can have a major effect on molecular structures and energies. For some molecular types such as a small organic molecule with only one polar functional group, solvation does not appear to be a major determinant of conformational energies. However, for molecules having several polar functional groups, the effect of solvent can be dramatic since the electrostatically least stable structures are often the most heavily solvated (stabilized) in a polar solvent.

Using an explicit solvent model is one approach; however, it has its own disadvantages. In particular, “explicit solvation” calculations run much more slowly because there are so many particles when hundreds of explicit solvent molecules are included. Furthermore, convergence is a problem in that longer simulations or different solvent starting configurations often give different final energies. Consequently, simple energy minimization is not useful in an explicit solvent.

MacroModel uses an alternative solution model which treats the solvent as a fully equilibrated analytical continuum starting near the van der Waals surface of the solute. The model is termed the GB/SA model and is described in the literature [16]. MacroModel is provided with parameter files for water (`water.slv`), octanol (`octanol.slv`), and chloroform (`chcl3.slv`). Using the GB/SA model slows calculations by a factor of approximately three relative to the gas phase. However, because of the increased accuracy of modeling in solvent, it is suggested that the GB/SA continuum solvation model be used in all calculations for molecules for aqueous solutions. Solvation is controlled by the Solvent: option menu, located in the Potential folder of the MacroModel panels within Maestro, or by the `SOLV` command in MacroModel command files.

The calculation of Born radii in the GB part of the GB/SA model is performed by doing a volume integral. Stretch, bend, and non-bonded pairs (including 1-4 interactions) contribute to this integral. For large systems, such as proteins in which non-bonded cutoffs typically are used and not all non-bonded pairs are included on the non-bonded pair list in order to expedite calculation, Born radii are subject to systematic error. This is because the volume integration is performed by using pairs on the stretch, bend, and non-bonded pair lists, and the latter excludes longer range interactions when non-bonded cutoffs are in use. To correct for this error, MacroModel calculates the contribution from such longer range non-bonded pairs every time a non-bonded pair list update is done. This contribution is taken to be a fixed value and is used in energetic and derivative calculations. Although not an exact correction, inclusion of this contribution in this manner does help significantly. The correction is enabled by default;

see the description for debug flags 830 and 832 under the DEBG opcode description in the *MacroModel Reference Manual* for information on disabling the correction.

In calculations using continuum solvation, MacroModel uses an approximate solvent accessible surface area function for derivatives and then computes final energies with a more exact area function at the end of the calculation. Consequently, the intermediate energies which are listed during energy minimization iterations will differ from the final energies.

See the *MacroModel Technical Manual* for a discussion on the performance of the GB/SA model with different sources of partial atomic charge.

3.4 Truncation of Electrostatic Interactions

For large systems, interactions between non-bonded pairs separated by more than a given distance need to be ignored in order to make simulations tractable. Such treatment for van der Waals interactions is less necessary since the interaction dies off as $1/r^6$ where r is the inter-atomic distance. Electrostatic interactions, on the other hand, die off as $1/r$, and pose a greater challenge because of their long range nature.

MacroModel versions prior to 8.1 used residue-based cutoffs for systems that had residue information in their coordinate files and atom-based cutoffs for systems with only one residue or without such information. MacroModel versions from 8.1 on include a new method for the truncation of electrostatic interactions. It is termed “Bond Dipole Cutoffs” (BDCO) and leverages the physics of charge-charge, charge-dipole, and dipole-dipole interactions to give very accurate absolute electrostatic and generalized Born energies. An example is the benzamidine/trypsin complex (PDB code 3ptb), which contains 3245 atoms in an all-atom representation, for which results are given in Table 3.1. In addition to a smaller error in energy by two orders of magnitude than that seen with residue-based cutoffs, the BDCO calculation uses fewer non-bonded pairs. Because of this, the BDCO calculation runs slightly faster than the residue-based cutoff calculation. Similar results are seen for generalized Born solvation polarization energies.

Table 3.1. Use of different cutoffs in a calculation on the benzamidine/trypsin complex.

Method	Number of Non-bonded Pairs in Energy Calculation	Total Electrostatic Energy (kJ/mol)	Error (kJ/mol)
BDCO	2731027	-28697.73	2.08
res ^a	3628479	-29179.45	479.64
all ^b	5254184	-28699.81	0.0

a. Residue-based cutoffs

b. All non-bonded pairs included in calculation

The following subsections provide an overview of the physical basis and implementation of this novel method for truncating electrostatic interactions.

3.4.1 Energetics of Charges

Coulomb's law states that the potential energy of two point charges separated by a distance r scales as $1/r$. In contrast, the energy of a point charge interacting with a dipole scales as $1/r^2$ and the energy of two interacting dipoles scales as $1/r^3$. These physical properties determine the length scales for which truncation of the various interaction types is appropriate in order to achieve a given level of accuracy.

3.4.2 Molecular Mechanics Description of Charge

Common molecular mechanics force fields utilize Coulomb's law to describe the electrostatic interaction between atoms due to the uneven distribution of electron density across bonded atoms with different electronegativities. A carbonyl moiety, for example, is described with a partial negative charge located at the oxygen atom center and an equal but opposite charge located at the carbon atomic center. These fixed partial charges can be directly used in Coulomb's law to give molecular mechanics electrostatic energies and the associated forces.

For groups with net charge, such as a carboxyl moiety, the formal charge of the group must also be accounted for. This can be done by delocalizing the formal charge over the appropriate atoms and adding this delocalized formal charge to the partial charges due to the dipolar nature of bonds between different atoms. Thus, for the carboxyl moiety, one can think of the system as two dipoles, one across each C-O bond, as well as point charges of -0.5 located in the oxygen atoms to account for the net negative charge of the group.

Some force fields, such as MMFF, have been conceptualized using this scheme, and the charge parametrization for an arbitrary molecule is achieved by assigning dipoles across bonds as well as delocalized formal charges on atomic centers and summing up these values to obtain fixed partial charges on atomic centers. Other force fields, such as AMBER*, assign charges to arbitrary molecules by dividing a given molecule into small groups of atoms with unit charge, like amino acids, for which there exist fixed partial charge data. Even in the latter case, the fixed partial charges can be transformed into a system of bond dipoles and delocalized formal charges.

In addition to assigning partial charges based on the selection of a force field, MacroModel allows you to specify custom partial charges for a system in the input structure file (see the CHGF opcode in the *MacroModel Reference Manual*). MacroModel is capable of decomposing this list of partial charges into delocalized formal charges and bond dipoles. In the event that the partial charges so specified for a molecule do not sum to the net formal charge of that mole-

cule, the input structure partial charges will be modified to meet this constraint before the decomposition is performed.

3.4.3 Bond Dipole Cutoffs (BDCO)

Since a system of atoms can be described as consisting of bond dipoles and atom-centered delocalized formal charges, the energetics of the charges can be considered as due to charge-charge, charge-dipole, and dipole-dipole interactions. This conceptualization is at the heart of the Bond Dipole Cutoff (BDCO) method. By applying different cutoff distance criteria to these three different types of interactions based on how they scale with distance, energetic accuracy is preserved while retaining the benefits of having a short non-bonded pair list. The implementation of BDCO *does not* involve any change in the functional form of a molecular mechanics force field. Coulomb's law in its familiar form is preserved. Rather, the BDCO algorithm applies different, successively shorter, cutoff distances to charge-charge, charge-dipole, and dipole-dipole interactions, where the charges are delocalized formal charges and the dipoles are bond dipoles. The resultant list of interactions is then used to generate charge products for pairs of atoms; the "charge product" is defined as the numerator in Coulomb's law. The charge product for a given pair of atoms is used directly in Coulomb's law to give the electrostatic energy and associated force for that pair.

3.4.4 BDCO and Molecular Modeling

The reason BDCO works is that although long cutoffs are used for charge-dipole and charge-charge interactions, there are relatively few such interactions in proportion to the number of dipole-dipole interactions in systems commonly modeled using molecular mechanics. For example, of the 20 types of amino acids seen in proteins, only 5 usually have charged (Lys, Arg, His, Glu, Asp) side chains. For a uniform distribution of amino acid types, 1/4 of the amino acids in a sequence will be charged. Assuming approximately 15 atoms per amino acid, the delocalized formal charge for a charged amino acid will on average reside on 2 of the 15 atoms. $1/4 * 2/15 = 1/30$, so only 1 of every 30 atoms in an average protein is expected to have delocalized formal charge. The number of dipoles in the system is approximately equal to the total number of atoms. Since delocalized formal charges occur on 1/30 of the atoms, the proportion of charge-charge interactions is $1/30 * 1/30 = 1/900$ of the total number of charge-charge, charge-dipole, and dipole-dipole interactions. Likewise, the proportion of charge-dipole interactions is $1/30 * 1 = 1/30$. The number of charge-charge and charge-dipole interactions in such a system is very small compared to the number of dipole-dipole interactions, and assigning long cutoffs to the former two types of interactions will not impact the size of the non-bonded list. However, use of long cutoffs for such interactions will positively impact the accuracy of the calculation due to the long range nature of these interactions.

Refer to the BDCO opcode in the *MacroModel Reference Manual* for more information.

3.5 General Guidelines for Convergence

It is as important to fully explore the potential energy surface as it is to describe that surface correctly by means of a high-quality force field. We use the term “convergence” to describe this issue. A common problem in modeling is obtaining unconverged results, i.e., the calculation gives a significantly different answer if the calculation continues or is repeated from different initial conditions. The term “convergence” has somewhat different implications, depending on the nature of the procedure that is being carried out, and for each type of energy calculation, there are different convergence issues. These issues are discussed in the following sections.

3.5.1 Minimization

It is important that, during energy minimization, the energy is minimized to a low gradient norm. Energy minimization calculations generally converge readily, although there is no guarantee that final structures are low in energy relative to the actual global minimum. By default, we set a convergence criterion for minimization of a gradient to < 0.05 kJ/Å-mol, and this is usually satisfactory. However, the default setting for the maximum number of iterations for a minimization calculation may not be sufficient to meet the convergence criterion, particularly for large structures, or ones with poor starting geometries. The latter situation often arises during conformational searching. For this reason, we suggest that all structures located in a conformational search be reminimized with a large value (5000) for the maximum number of iterations and with all the tests for uniqueness in place. This will help ensure that the set of conformations obtained in the conformational search represent unique conformations and not structures that would be duplicates of other structures if completely minimized. The final report for a conformational search highlights with an asterisk any structure that did not converge to the criterion in force during minimization.

3.5.2 Conformational Searching

Conformational searches involving structures with up to a dozen rotatable bonds generally converge easily, but convergence can be problematic with more flexible structures. When performing conformational searches with these more flexible molecules, it is therefore wise to carry out multiple searches with different starting geometries to verify that the same final structures are found.

In addition to ensuring that the final set of conformers obtained from a conformational search is completely minimized, there is another issue of convergence which affects the stochastic searching methods (LMOD, LLMOD, MCMM, and SPMC). For these methods, it is important to establish that the structures obtained in a given search represent a thorough sampling of all

possibilities. There is no way to prove that the search was fully convergent, but there are several methods which can at least indicate that it was not.

As the search approaches convergence, the number of new unique conformations which are found should begin to approach zero. That is to say, after the search has been in progress for some time, many, many trials should elapse between the finding of previously unobserved unique structures. In practice, however, this may not be a good test unless all conformations found minimized to low gradients, as described in the previous section.

We generally recommend the following procedure to ensure that a stochastic search has been exhaustive:

- Conduct a search of approximately 1000-2000 MC steps for every variable degree of freedom.
- Ensure all the structures saved are thoroughly minimized as described above.
- Run a second search (possibly “seeded” with the selected structures from the of the first search—see the *MacroModel Reference Manual* for details on how to do this). Ensure that a SEED command is used so that a different sequence of random numbers will be used in the search.
- Combine the results of the two searches and reminimize keeping only the unique conformations.
- Repeat this procedure of searching and combining the results until the total number of accumulated conformations (especially those in the low energy region) becomes constant.

In any modeling study, it is up to you to demonstrate that results are converged and what the bounds of uncertainty are. *Unconverged results are meaningless.*

3.5.3 Molecular and Stochastic Dynamics

Obtaining convergence in dynamics simulations has traditionally been problematic, because of the slow frequency at which systems undergoing these processes cross barriers between the various minima on the potential energy surface. Stochastic dynamics may search conformational space more efficiently than does regular molecular dynamics, but neither method gives frequent crossing of barriers much larger than 3 kcal/mol (13 kJ/mol). A much higher rate of convergence can be obtained by the use of the mixed-mode Monte Carlo/Stochastic Dynamics (MCSD) procedure. We therefore recommend using this procedure whenever appropriate, such as for simulations of acyclic systems in which generation of the canonical ensemble is desired and no time-dependent information is required. For cyclic systems, the Jumping Between Wells (JBW) method can be used. However, even with these methods, it is important to test

that converged results have been obtained. Converged simulations will satisfy the criteria listed below. Satisfaction of these criteria, like those described in the last section, are necessary, but not sufficient, to prove convergence:

- Average quantities—for example, potential energy, or the fraction of times a hydrogen bond appears—should exhibit numerical stability. For example, once the simulation is converged, doubling the length of the simulation should not change the averages. It is possible to obtain the average values of geometric and hydrogen bond monitors periodically during the simulation.
- Dihedral angle distributions in achiral compounds should be symmetric. For example, dihedral angle distributions produced using the MDDA command should exhibit the same populations of +gauche and -gauche values for three-way torsions.
- Simulations starting from different starting geometries should give the same results.

3.6 Problem Size: Limits of Operation

MacroModel is distributed with a number of arbitrary fixed limits on its operation. Some of these restrictions have been lifted since the program has gone to dynamic memory allocation, but some still remain. For example, we distribute a version of MacroModel which can perform Monte Carlo conformational searches on up to 500 rotatable bonds; however, these limits usually far exceed the size of a calculation that can be reasonably performed. For example, it is highly unlikely that a converged conformational search could be performed in a reasonable amount of time on a system with even 100 degrees of rotational freedom. Because of the wide variety of platforms on which MacroModel is used—ranging from Linux-based PCs to super computers or a network of scores of workstations—and also because chemical problems differ greatly, we cannot make categorical pronouncements setting forth the “largest” problem one can reasonably address. Instead, we encourage you to consider the issues of convergence as described below, and use the principles described there as a guide.

3.7 Nonbonded Cutoffs

By default, MacroModel energy calculations are performed with cutoff distances in place for nonbonded interactions. If the distance between any nonbonded pair of atoms is greater than the cutoff distance, then the non-bonded term (electrostatic or van der Waals energies) for that pair will be ignored. These terms decrease in magnitude with increasing interatomic distance, and so terms for nonbonded atom pairs separated by a distance greater than a reasonable cutoff distance will not contribute greatly to the overall energy. The introduction of nonbonded cutoffs greatly reduces the time it takes to compute the energy for a large molecular system. Several implications of using cutoffs, however, deserve special mention.

- For systems with formally charged atoms, the default cutoff distances (7 Å for van der Waals interactions and 12 Å for electrostatics) are probably too small for accurate results. Better values might be 8 Å and 20 Å, respectively. These values can be set by selecting Extended from the Cutoff option menu in the Potential folders of the Maestro MacroModel panels.
- The original MM2 and MM3 programs do not use non-bonded cutoffs. When making comparisons with these programs, MacroModel calculations should be performed using non-bonded cutoff distances which are effectively infinite; that is, the values should be larger than the largest molecular dimension. See [Section 3.1 on page 31](#) for other considerations involving comparisons between MacroModel results and results from the standard implementations of other force fields.
- We recommend that calculations performed using the MacroModel solvation model be carried out with effectively infinite cutoffs.

3.8 Electrostatic Treatments

By default, MacroModel uses a constant dielectric constant during the evaluation of electrostatic interactions between atoms. In most situations the use of a constant dielectric is appropriate for such calculations. One such place is when making comparisons with gas phase ab initio molecular orbital calculations. Additionally, use of the MacroModel GB/SA continuum solvation model also requires the constant-dielectric electrostatic treatment. When this model is utilized, electrostatic interactions are automatically calculated using this method. Alternatively, a distance dependent dielectric electrostatic treatment may be specified by selecting Distance Dependant from the Electrostatic Treatment option menu in the Potential folder of the MacroModel panels.

3.9 File Conversion

Maestro reads Protein Data Bank (PDB), MDL SD (.mol) files, and sybyl mol2 formatted structure files. You can then use these structures in MacroModel calculations. Note, however, that Maestro may have some difficulty with some of these files. It may be necessary to examine the file in detail and to make repairs. If you have problems reading a PDB file, please read Maestro's "PDB Conversion" online help topic. This topic explains Maestro's color-coded warning system. File conversion can also be performed from the command line using the `pdbconvert`, `sdconvert`, `molmmmod`, and `mmmodmae` scripts in the `$SCHRODINGER/utilities` directory.

3.10 Modified Force Field, Solvent, or Atom Type Files

When a job is launched, MacroModel searches in the current working directory for a copy of the force field (*.fld), solvent(*.slv), and atom type(atom.typ) files. If it finds any of these files there, it uses them to perform the energy calculation. If MacroModel does not find copies of these files locally, it then looks in \$HOME/.schrodinger/macromodel. Only if MacroModel does not find the files in .schrodinger does it use the global settings in \$SCHRODINGER. If you do create locally modified files, manage them closely. The unintentional use of forgotten locally modified files can give surprising energetic results.

3.11 Hydrogen Treatment

3.11.1 Hydrogens on Hetero Atoms

While the AMBER* and OPLS force fields can be used without explicit hydrogen atoms on carbon atoms, all force fields require the use of explicit hydrogen atoms on hetero atoms. Note also that AMBER* requires explicit lone pairs on sulfurs. Original MM2 implementations utilize lone pairs on sp³ oxygens; however, we recommend removing them for best results in MM2*.

In Maestro, you can modify a structure's hydrogen treatment using settings in the Hydrogen Treatment panel, which you open from the Build menu.

3.11.2 Hydrogen Treatment in the AMBER* Force Field

The AMBER* force field file contains charge sets appropriate to either united-atom or all-atom representations; however, by default, the united-atom charges are used. This means that even those carbons having explicit hydrogens attached will, by default, get united-atom charges. In practice, this makes little difference in energetic results, but in some situations—for example, when making comparisons with the original AMBER program, and also when calculating absolute solvation free energies—it is important to use all-atom charges if such carbons exist. This can be done by selecting Force Field Defined from the Electrostatic Treatment option menu in the Potential folder of the MacroModel panels.

When using the united-atom models in AMBER* or OPLS*, the best results are obtained when explicit hydrogen atoms are used on sp² carbons. The command:

```
$SCHRODINGER/utilities/applyhtreat in.mae out.mae
```

where *in.mae* is the name of a file containing structures to which hydrogen atoms should be added produces a file called *out.mae* containing the corresponding structures containing the additional hydrogen atoms.

Running MacroModel From the Command Line

Running MacroModel from the command line gives you the ability to customize batch processing. In addition, some types of MacroModel calculations are not supported by the Maestro GUI. This chapter gives instructions on how to launch MacroModel jobs from the command line. For information on starting MacroModel jobs from the Maestro graphical user interface, see the chapter in this manual that corresponds to the particular type of job you want to submit.

4.1 Preparing for MacroModel Calculations

Before you can submit MacroModel jobs from the command line, you must set the `SCHRODINGER` environment variable and must have a valid input structure file and a command (*jobname.com*) file.

4.1.1 Environment Variables

The `SCHRODINGER` environment variable specifies the installation directory for your Schrödinger software. These files include the force-field files, the solvent files, and the `atom.typ` file. To run MacroModel as a stand-alone application, you must first define this environment variable. See [Section 2.2 on page 5](#) for more information on setting `$SCHRODINGER`. Whether you start MacroModel from the UNIX shell or from Maestro, local versions of the force-field, solvent, and `atom.typ` files override versions in the location specified by `$SCHRODINGER`.

You might also want to specify the default location for temporary (scratch) files. The default scratch directory is set in a `schrodinger.hosts` file, which should exist on each host on which you want to run MacroModel. This file can be in a local directory, in `$HOME/.schrodinger`, or in `$SCHRODINGER`. You can override the default by setting the `SCHRODINGER_TMPDIR` environment variable.

4.1.2 The Command File

Every time job files are written or a job is started from Maestro, a MacroModel command file and a structure file are created. The command file contains all necessary information for MacroModel to perform the job as specified. Most MacroModel tasks can be executed from

Maestro. However, there are a few tasks that still involve manually setting up and executing a command file.

In this chapter, examples of common command files are given that are set up through Maestro and provide some insight to the different commands included in these files. At the end of the chapter, more advanced command files are presented, including tasks that can not currently be set up from Maestro.

Many of the commands presented in the command files below have additional arguments available. For more information on any of the commands, refer to the [MacroModel Reference Manual](#).

4.1.3 The MacroModel Command File Format

To run MacroModel calculations, a molecular structure file and a command file are required. The molecular structure file contains the structures to be used as input in the calculation. The command file contains the name of the input structure file, the name of the output structure file, and an ordered list of operation codes, or *opcodes*, for the calculations.

A generalized form of a MacroModel command file is given below:

```
jobname.mae
jobname-out.mae
MMOD      0      1      0      0      0.0000      0.0000      0.0000      0.0000
FFLD      10     1      0      1      1.0000      0.0000      0.0000      0.0000
BDCO      0      0      0      0      20.0000 99999.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
CONV      2      0      0      0      0.0500      0.0000      0.0000      0.0000
MINI      1      0      500     0      0.0000      0.0000      0.0000      0.0000
```

This file must have a name in the form of *jobname.com*, where *jobname* is replaced with the actual name of the job. Thus, a job that was actually named “jobname” could be submitted from the command prompt with a command such as

```
$SCHRODINGER/bmin jobname
```

Maestro uses a similar command for computations prepared with the interface. The first line of the instruction file above is the name of the input structure file. The input structure file can be named with any valid UNIX filename. The resulting output structure file is given the name listed on the second line of the instruction file. The full path to the structure files may be given if the files are not in the current directory.

In addition, any input substructure (.sbc) or velocity (.vel) files should contain the same prefix as the input structure file. Similarly the output energy listing (.mmo), substructure, and dihedral drive (.grd) files have the same base name as the output structure file. The

jobname.log file contains text messages tracing the progress of the job. Experienced users of MacroModel should be aware that this mechanism has not changed from previous versions of MacroModel.

There are some default file naming changes, indicated in the sample file above, that are used by the Maestro interface due to the advent of Maestro formatted structure files. Although not mandated, to be consistent, all Maestro formatted structure files are given the suffix *.mae*. This is a different default behavior from the former MacroModel user interface and the older MacroModel structure format. Previously, input files were given the names *jobname.dat* and *jobname.out*. Now that all Maestro formatted structure files are, by default, named with the *.mae* extension, Maestro automatically names the input structure file *jobname.mae* and the output structure file *jobname-out.mae*.

The remaining lines in the command file provide the instructions to MacroModel concerning the type and order of calculations to be performed. The opcode lines must be of the following fixed format:

```
#OPCD 123456 123456 123456 123456 FFFFF.FFFF FFFFF.FFFF FFFFF.FFFF FFFFF.FFFF
```

Each opcode has four letters, and is preceded by a blank space. A specific opcode can be ignored (commented out) in a command file by placing a character other than blank space before the opcode. The eight fields after the opcode are referred to as *arguments*, and are often referred to simply as “arg”s. Arg1 through arg4 are integer arguments (Fortran I6 format), while arg5 through arg8 are floating point arguments (Fortran F10.4 format). The opcodes indicate individual energetic calculations or options and the arguments indicate additional options or quantify the parameters of the calculation. Many arguments have default values, which are indicated by a value of zero as the argument. Thus is it unnecessary to indicate explicitly all arguments in the instruction file if the default values are sufficient. The default values are included with the opcode descriptions in the *MacroModel Reference Manual*. It is important to strictly adhere to the format of the command file. We recommend using an existing command file as a template, rather than to build one from scratch. Tabs are not allowed in MacroModel instruction files.

4.2 Submitting Jobs From the Command Line

The *bmin* shell script located in the *\$SCHRODINGER* directory should be used to launch stand-alone MacroModel jobs from any directory, once the *SCHRODINGER* environment variable has been set. To launch a job, enter the following command in a terminal window:

```
$SCHRODINGER/bmin [options] jobname
```

replacing *jobname* with the stem of the input file name (not needed if the `-TEST` or `-HELP` options are specified) and including any of the following options.

Table 4.1. Options for the `bmin` command.

Option	Meaning
<code>-ARCH</code> <i>pattern</i>	The platform tag must match <i>pattern</i> (e.g., <i>mips4</i>).
<code>-COMPAT</code> <i>execdir</i>	Require compatibility with executables in <i>execdir</i> .
<code>-DEBUG</code>	Show details of operation of the top level script.
<code>-DEBUGGER</code> <i>debugger_name</i>	Run <code>bmin</code> in the foreground with a debugger called <i>debugger_name</i> .
<code>-DEBUG2</code>	Print out debugging information from the scripts used to start <code>bmin</code> .
<code>-HELP</code> <code>-help</code> <code>-h</code>	Print the <code>bmin</code> command syntax, options, and their definitions and exit.
<code>-HOST</code> <i>hostname</i>	Run job remotely on the indicated host.
<code>-HOSTFILE</code> <i>hostfilename</i>	The name of the hosts file to use for this run.
<code>-INTERVAL</code>	The maximum time in seconds for updating the monitoring files.
<code>-LOCAL</code>	Do not place files in a temporary directory. Keep files in the local directory.
<code>-NICE</code>	Run the job at reduced priority.
<code>-NO_REDIRECT</code>	Run in the foreground and send output to standard output.
<code>-PROJ</code> <i>name</i>	Assign job to a Maestro project <i>name</i> .
<code>-REL</code> <i>version</i>	The version number must match <i>version</i> (e.g., 80 or 80015).
<code>-TEST</code>	Run the standard test suite (<i>jobname</i> not needed).
<code>-TMPDIR</code> <i>directory</i>	The name of the directory used to store files temporarily during a job.
<code>-USER</code> <i>username</i>	Launch job as user <i>username</i> .
<code>-VER</code> <i>pattern</i>	Executable directory pathname must match <i>pattern</i> .
<code>-WAIT</code>	Do not return a prompt until the job finishes or is submitted to a batch queue.

Additional diagnostic options, which do not run the job but provide information are listed in Table 4.2.

Table 4.2. Options to the `bmin` command that provide information only.

Option	Meaning
-ALL	Ignore platform compatibility (list all platforms).
-ENTRY	List the relevant host entries from the hosts file.
-HOSTS	List the hosts available for remote jobs.
-LIST	List all the platform-compatible versions of the product.
-WHICH	Report the product executable directory and <code>MMSHARE_EXEC</code> .
-WHY	Show details of how the selected executable directory was chosen.

Launching a MacroModel job using the `bmin` script creates the output structure file and the `jobname.log` file. The `.log` file contains MacroModel job information, such as job number and C-shell timings, and any system errors, as well as output from the job.

4.2.1 Using an Executable Command File

If you have a command file that must be run repeatedly (such as the MacroModel test scripts), you can embed the contents of the command file in an executable shell script. Below is an example, though we generally recommend using shell scripts, as described in the previous section.

```
#!/bin/csh -f
cat <<EOC > a_run.com
jobname.mae
jobname-out.mae
FFLD
READ
ELST      2
EOC
$SCHRODINGER/bmin a_run
```

If the above contents were contained in an executable file called `mm2-ELST`, you could invoke MacroModel to run the included commands by simply naming the file on the command line, e.g.,

```
mm2-ELST &
```

4.3 Running Remote and Distributed Calculations

Remote MacroModel jobs run on a different host from the one on which the command is entered to start the job. Distributed MacroModel allows certain types of calculations to be run in a distributed fashion across a number of different machines. Remote MacroModel jobs can be run from Maestro, but at present, distributed MacroModel jobs can only be run from the command line.

Schrödinger's Job Control facility controls both remote and distributed MacroModel jobs. The Job Control facility enables remote and distributed MacroModel jobs to run reliably on a wide range of computer platforms. If you intend to run jobs on various hosts, you must set up the hosts for remote access and provide information on the hosts to the Job Control facility through a file named `schrodinger.hosts`. How to provide this information is described in the [Job Control Guide](#). Brief instructions for setting up the required information on remote hosts and for configuring batch queues is given in the [Installation Guide](#).

Any MacroModel job can be run remotely. The normal commands are used to run the job, with additional information specifying which host from the hosts file to use:

```
$SCHRODINGER/bmin -HOST remote_host_name jobname
```

MacroModel can divide some types of calculations into “tasks”. These tasks may be distributed over a number of processors to reduce the calculation time. Calculations that can be distributed using MacroModel include many of the conformational searching methods as well as Free Energy Perturbation (FEAV/FESA) and MBAE calculations. For a thorough description of the distributed MacroModel calculations, including the specific types of calculations supported, see the [MacroModel Reference Manual](#). Below are two simple examples that illustrate how to perform these calculations. To run these examples, first ensure that the hosts and accounts have been prepared and an appropriate `schrodinger.hosts` file is available.

For distributing calculations focused on a single input structure, internal distributed calculations are used. For example, to distribute a conformational search of a protein-ligand complex across two processors, add the line:

```
NPRC      2      20      10      1      0.0000 0.0000 0.0000 0.0000
```

just after the `MMOD` command in the example `.com` file given in [Section 10.3.3 on page 117](#). The job can then be started just like any other MacroModel job:

```
$SCHRODINGER/bmin jobname
```

Another type of job that can be distributed is MCMM serial searches in which each input structure is subjected to a separate search (see [Section 10.3.2 on page 116](#)). You can do this with the `para_bmin` command, which can be run directly on the unmodified non-distributed files:


```
$SCHRODINGER/utilities/para_bmin -NJOBS 5 -HOST "comp1:1 comp2:2" serial-lmcs
```

This command divides the job into five tasks, all of which are run on computers `comp1` or `comp2`. These computers must be described in the `schrodinger.hosts` file. At any one time, up to one task is running on `comp1` and up to two tasks are running on `comp2`.

For distributed MacroModel calculations that involve a random or accumulated aspect (e.g. MCMM conformational searches) the results may differ from those obtained from an equivalent non-distributed calculations. However, the different results from either type of calculation are valid.

4.4 Interacting With a Running MacroModel Job

The Schrödinger Job Control facility allows you to interact with running MacroModel jobs. Because all jobs submitted, either from Maestro or the command line, are controlled by the Job Control facility, all jobs can be monitored using the Monitor panel in Maestro or using the `jobcontrol` script from the command line. For specific instructions on using the Monitor panel or the `jobcontrol` script, see the *Job Control Guide*.

General Energetic Settings

5.1 General MacroModel Calculation Setup

All Maestro energetic calculation panels have the same basic setup process:

1. Specify the job input source from the Use structures from option menu (see [Section 5.1.1](#) below).
2. Specify settings in the tabbed folders:
 - For the Potential folder, see [Section 5.2 on page 54](#)
 - For the Constraints folder, see [Section 5.4 on page 59](#)
 - For the Substructure folder, see [Section 5.5 on page 61](#)
 - For folders specific to the energetic calculation, see the appropriate chapter
3. Click Start to set up and launch the job (see [Section 5.6 on page 64](#)) or click Write to write the files for future use (see [Section 5.7 on page 65](#)).

5.1.1 Specifying Job Input Source

For most MacroModel jobs run from the Maestro interface, you can use as input either the structures included in the Workspace or the entries selected in the Project Table. For these jobs, if you have a file containing the structures of compounds you want to evaluate, first import the structures into a project using the Import panel, and then select the desired structures in the Project Table. In addition, there are some MacroModel calculations, such as eMBRAcE Minimization, that allow you to use an input file directly, without having to import the structures into the Project Table.

The default setting for MacroModel job input is Workspace (included entries). If you want to perform MacroModel calculations only on the structures currently appearing in the Workspace, you do not need to change the Use structures from setting. However, if you want to use Project Table entries as input, change the Use structures from setting to Project Table (selected entry) or Project Table (selected entries), depending on which of these is present. When using this setting, make sure that you have selected the desired entries in the Project Table panel.

5.2 The Potential Folder

The Potential folder appears on all of the MacroModel energy panels. This folder contains tools for specifying a force field, a solvent treatment, and an electrostatic treatment for the calculation. The Potential folder is displayed by default when the energetic panel is opened.

Selecting a Force Field: To specify a force field for a MacroModel calculation, select the desired force field from the Force field option menu. The default force field is MMFFs. For information about MacroModel implementation of the supported force fields, see [Section 3.1 on page 31](#).

Selecting a Solvation Treatment: MacroModel calculations are carried out in the gas phase by default, but you may elect to carry out the calculation with a solvation treatment. MacroModel provides the GB/SA continuum solvation treatment, and water, octanol, or chloroform may be chosen. To change the solvent, select the name of the desired solvent from the Solvent option menu on the Potential folder of the energy panel.




Figure 5.1. The Potential folder.

Selecting an Electrostatic Treatment: The default electrostatic treatment for all built-in force fields is to use a dielectric constant of 1.0. You can also select a distance-dependent dielectric from the Electrostatic treatment option menu. You can set the dielectric constant; the default value is 1.0, which corresponds to Coulomb's law in a vacuum. Selecting the force field defined treatment is equivalent to choosing a constant dielectric with a value of 1.0. The MMFF and MMFFs force fields use a "buffered" constant dielectric treatment; see Halgren [9,10].

Using Charges From a Structure File: The charges used in the electrostatic portion of an energy calculation can either be assigned by the force field or obtained from the structure. Regardless of the source, charges are written to the structure file when a job is started. By default, the Force field option is used. To use charge information from the structure, select Structure file from the Charges from option menu.

Using an Existing Command File: To read a MacroModel command file and have the settings on the energetics panel update to reflect all potential energy settings, click Read Settings From Command File. Setting up MacroModel calculations in this way is helpful because you can easily reproduce your potential energy settings for multiple calculations.

Choosing a Nonbonded Cutoff Setting: The maximum distances over which hydrogen bonding, van der Waals, and electrostatic contributions to the molecular potential energy are computed may be controlled using the Cutoff option menu. Default cutoff distances are 4 Å for hydrogen bonding, 7 Å for van der Waals and 12 Å for electrostatics. The van der Waals and electrostatic cutoff distances are the center of a soft cutoff that starts at 1 Å smaller than the specified distance and ends at 1 Å larger than the specified distance.

5.3 Fixing and Freezing Atoms in Large Systems

Two important reasons for imposing constraints on molecular systems are to force the modelled system to meet desired geometric conditions, and to reduce the cost of the calculation by eliminating interactions that are expected to have little influence on the results. When modeling molecular systems, you might want to constrain some of the degrees of freedom, such as atom positions or distances, angles and dihedral angles. Constraints are useful in the following situations:

- When it is inherently hard to get the system to adopt the desired geometries in the absence of constraints
- When the potential functions are inadequate for the system
- When parts of the system that would normally lead to the appropriate geometries are not represented

When modelling a large system, it is often desirable to focus calculations on important regions and ignore any large portions of the system from which no significant influence is expected. For example, in the minimization of a protein-ligand complex, the region of interest is the ligand and the area surrounding the ligand. The large outer regions of the protein may not significantly influence the structure of the minimized ligand or its vicinity, and these regions may therefore be excluded from the calculation. Fixing or freezing remote portions of the system can help further reduce computational time and provide buffer zones between fully moving and ignored regions.

Frozen atoms cannot move, while fixed atoms are restrained to lie close to the requested positions. Moving atoms interact with all other moving, fixed and frozen atoms in the usual way. However, among the fixed or frozen atoms only stretch interactions involving fixed atoms are considered. Eliminating most of the interactions among the fixed or frozen atoms can save a great deal of computer time and is a useful approximation provided that the potentials restraining the fixed atoms are strong enough to prevent dramatic distortions of the molecular structures.

5.3.1 Methods for Freezing and Fixing Atoms

In MacroModel there are two ways of constraining parts of a molecular system:

- **Using constraints:** Manually fixes atoms (and also distances, angles, or torsions) through the use of `FXAT` (or `FXDI`, `FXBA`, or `FXTA`) commands, either using the Constraints panel in Maestro or by directly modifying a command file. Only the geometry of the items explicitly specified is affected. All other parts of the system interact and move normally.
- **Using substructures:** Uses Maestro's substructure facility to define, within a molecular system, flexible substructures and fixed and frozen shells. You can also directly modify a command file to create or modify substructures. Atoms that are not in the substructure or specified as fixed or frozen are ignored.

To constrain only a few atoms (the first method) using `FXAT` commands is best, and is also easily achieved through point-and-click operations in Maestro. To constrain a larger number of atoms, you will find it is much more efficient to use substructures.

Note: Constraints can be violated if you use automatic setup (see [Section 10.2.4 on page 103](#)). If so, MacroModel will fail.

With either method, a positive force constant can be specified for the restraining harmonic potential for fixed atoms. Free movement of the atoms within a specified range can be facilitated with the use of flat-bottomed potentials in which the two halves of the harmonic potential are separated by a specified distance within which the potential is zero.

5.3.2 Freezing and Fixing Individual Atoms

Individual atoms can be constrained either by picking the desired atoms using the tools in the Constraints folder of the MacroModel panels within Maestro, or by using the `FXAT` command in a `jobname.com` file.

Atoms are fixed or frozen in their original position by using tools in the Constraints folder within the various MacroModel energy panels in Maestro. The desired atom is selected through on-screen picking, and is then assigned values for how much it should be allowed to move (half-width, in Å) and the penalty for moving outside the defined limits (force constant, in kJ/mol-Å²). This process is repeated for all atoms to be constrained. See [Section 5.4 on page 59](#) for more information on specifying constraints using Maestro.

For additional information about using the `FXAT` command when manually constructing MacroModel command files (`jobname.com`), refer to the pages describing the `FXAT` opcode in the *MacroModel Reference Manual*.

All atoms not defined by a `FXAT` command remain fully flexible, free to sample the potential energy surface defined by the force field.

5.3.3 The Substructure Facility

You can define substructures by manually placing the appropriate lines in the MacroModel command file (`jobname.com`), or by using the Substructure folder in the MacroModel energy panels to define groups of flexible, fixed, and frozen atoms. The settings are saved in a substructure file consisting of `SUBS` commands for the flexible atoms and `FXAT` commands for fixed/frozen atoms.

Note that atoms not defined by either `SUBS` or `FXAT` commands are completely ignored (as if they had been deleted) during the calculation.

After defining the flexible part (typically the ligand and atoms in its near vicinity) of the system, you can define a combination of fixed/frozen atom shells (having different force constants) around the movable part. These shells are defined based on a user-selected radius (in Ångstroms) from the flexible part of the system. To avoid cutting amino acid residues in two by the radius definition, Maestro offers a Fill residues option that ensures that if any atom of a residue is within the defined radius, the full residue is included in the shell. Additionally, there is an option for adding isolated atoms to any shell. Each shell can be assigned a different force constant (for a fixed shell), or can be completely frozen (force constant takes a negative value).

Looking at a protein-ligand complex as an example, the flexible part is most easily defined by picking. By choosing Molecules from the Pick menu in the Atoms for substructure section, you can select an atom in the ligand to define the whole ligand as part of the substructure. Atoms

belonging to the substructure are identified by white markers on screen. If parts of the active site are to be included in the substructure as well, selected residues can easily be picked by choosing Residues from the Pick menu and clicking on an atom in each residue to be included. Alternatively, a shell of flexible residues can be defined using the Atom Specification Language (ASL) in the ASL text box or in the Atom Selection dialog box. Assuming that the ligand is molecule number one and the protein is molecule number two, the following command in the substructure definition field places the ligand and a 3.5 Å shell of filled out residues in the substructure part:

```
m. 1 or fillres within 3.5 m. 1
```

The different shells defined around the substructure have markers of different colors (orange, purple, etc.), in order to easily visualize the various selections on screen.

After the substructure and shell selections are made, the substructure (*filename.sbc*) file can be written to disk. The name of the substructure file must match the input file. The substructure file has the following format:

Command	arg1	arg2	arg3	arg4	arg5	arg6	arg7	arg8
SUBS	1	2	3	4	0.	0.	0.	0.
SUBS	5	6	7	8	0.	0.	0.	0.
[...]								
FXAT	123	0	0	0	200.	0.	0.	0.
FXAT	124	0	0	0	-1.	0.	0.	0.
[...]								

For the SUBS part of the substructure file:

arg1-4: Atom numbers of atoms included in the substructure (flexible) part of the system.

arg5-8: Not normally used.

For the FXAT part of the substructure file:

arg1: Atom number of atom to be included in fixed/frozen part of the system.

arg2-4: Not used.

arg5: Force constant (kJ/mol Å²) used for atom listed in arg1. If the force constant is negative, the atom is frozen.

arg6-8: Desired X, Y, Z coordinates of atom defined in arg1. If all are zero (the default), the starting coordinates are used.

Note that the SUBS command in the MacroModel command file takes on two meanings, depending on whether or not a substructure file is to be used. If arg1 of SUBS is zero in the command file, MacroModel looks for a substructure file to obtain information on flexible and fixed/frozen parts of the system. There are other possible uses of the SUBS command: in partic-

ular it may be useful to combine the use of substructures both in a *filename.sbc* file and separate SUBS lines in the MacroModel command file. Please refer to the SUBS command in the *MacroModel Reference Manual*.

For more information on using Maestro's substructure facility, see [Section 5.5 on page 61](#).

5.4 The Constraints Folder

In MacroModel, subsets of atoms can be constrained in almost all energetic calculations, and the constraints can be of two basic types. The atoms can be completely frozen and not allowed to move, or atoms can be fixed, which means that the atoms in question are tethered in place using a constraint. Fixed atoms have the ability to move, frozen atoms do not. The following two sections describe how to prepare MacroModel computations from Maestro that include fixed and frozen atoms. For more information, see the *MacroModel Reference Manual* sections on FXAT, FXDI, FXBA, FXTA, and SUBS.



Figure 5.2. The Constraints folder.

Making selections in the Constraints folder adds the appropriate opcode entries into the MacroModel command file (*jobname.com*) for the computation being prepared. When you make selections in the Substructure folder, however, the constraints are added to a separate input substructure file, usually named *jobname.sbc*. This mechanism allows large numbers of constraints to be shared between different energetic jobs for the same structural system. Complementing this feature, Maestro can read and write substructure files. For more information about the Substructure Facility, see [Section 5.3 on page 55](#). Structural constraints can be added in both the Constraints folder and the Substructure folder for a single computation.

Note: Constraints can be violated if you use automatic setup (see [Section 10.2.4 on page 103](#)). If so, MacroModel will fail.

5.4.1 Constraining Atoms, Distances, Angles, and Torsions

Individual atom positions, bond distances, bond angles, and torsion angles can all be defined from the Constraints folder. When you click the button corresponding to the structural element you want to constrain, a panel is displayed, which you can use to pick the desired elements and define the relevant settings. The panels for defining constrained atoms, distances, angles, and torsions are generally the same. Exceptions are noted below.

Specifying Atoms to Be Constrained: To constrain an atom, distance, angle, or torsion, display the appropriate constraining panel by clicking on the corresponding button in the Constraints folder: Atoms, Distances, Angles, or Torsions. Choose a structural unit from the Pick menu and click on the desired atom(s) in the Workspace. The picked atoms appear in the list in the upper portion of the panel. To constrain many atoms simultaneously, you can use an ASL expression. ASL expressions can be built using the Atom Selection dialog box, which you open using the Select button. For more information on ASL, see the online help or the [Maestro Command Reference Manual](#).

Specifying the Force Constant and Atom Flexibility: By default the constrained structural elements have harmonic potentials with a force constant of 100 kJ/mol-Å² with the current value of the element used for the potential minimum. The force constant can be modified. For fixed distances, angles, and torsional angles, the position of the potential minimum can be adjusted. A flat-bottom potential, essentially a harmonic potential with the two halves separated by a region of 0 potential, can be requested by specifying the half width for the flat region in the +/-: text box.

Removing Constraints: To remove one constraint, select it from the list at the top of the specific constraints panel and click Delete. To remove all constraints of a particular type, click Delete All in the corresponding constraints panel. To remove all constraints of all types, click Reset All in the Constraints folder.

5.4.2 Freezing Atoms

“Frozen” atoms have no forces acting on them and do not move from their initial position in the structure. Frozen atoms do influence other, non-frozen atoms in the structure, mostly through nonbonded interactions. Frozen atoms differ from constrained atoms in that constrained atoms may move slightly.

To freeze atoms, click **Atoms** in the **Freeze** section of the energy panel. The **Frozen Atoms** panel is displayed. To select the atoms, choose a structural unit from the **Pick** menu and pick an atom in the **Workspace**. You can constrain many atoms simultaneously using an ASL expression in the **Atom Selection** dialog box, which you open with the **Select** button. To remove constraints from an atom, select its name from the list and click **Delete**. To remove constraints from all currently frozen atoms, click **Delete All**.

5.5 The Substructure Folder

The **Substructure** folder is used to set the definition of a substructure and any shells of constrained or frozen atoms that surround it.

If any substructure or shell atoms are specified, an input substructure file, *jobname.sbc*, is created and used. All instructions relating to the fixed and frozen atoms are entered in the input substructure file, and read from this file when the calculation is launched.

5.5.1 Defining a Substructure

You can select the atoms you want to include in a substructure definition by picking the desired atoms, chains, residues, molecules, or entries from the **Workspace** structure. To do this, select the structural unit you want to use from the **Pick** menu in the **Atoms for substructure** section, and then pick atoms in the **Workspace**. If **Show markers** is selected, the atoms in the substructure are marked with white markers.

In addition to defining substructure atoms by picking them from the structure, you can specify many atoms simultaneously by entering an appropriate ASL expression in the ASL text box in the **Atoms for substructure** section. For example, the expression:

```
fillres within 6.0 mol.n 3
```

picks all the atoms within 6 Å of molecule number 3, and also adds to the selection the remainder of any residues partially selected using the 6 Å proximity criterion. You can also construct ASL expressions using the **Atom Selection** dialog box. To open this dialog box, click **Select**.




Figure 5.3. The Substructure folder.

Once you have picked atoms from the Workspace structure or have specified substructure atoms using ASL, you can expand the “shell” of atoms to a specified distance from the defined substructure atoms. To expand the shell, enter the desired radius in the Expand to atoms within radius of text box.

To ensure that the substructure consists of only complete residues, select Complete residues. The defined substructures then include all the atoms in the residues to which the defined atoms belong.

5.5.2 Creating a Shell of Atoms

Defined substructures indicate which atoms may freely move during the energy calculations. Often, subsequent sets of atoms surrounding the substructure are intended to be fixed and/or frozen during the computation. This is accomplished by defining shells around the substructure in the Shells portion of the Substructure folder. To create a shell around a defined substructure,

click New Shell in the Substructure folder of the relevant MacroModel energy panel. If Show markers is selected, the shell atoms are indicated in the Workspace by colored markers: orange, purple, yellow, green, and so on.

Once you have created a new shell, you can define and modify its contents in the following ways:

Specifying a Shell Radius: You can assign a radius by selecting a shell from the list and then entering a value into the Radius text box. The shell includes all atoms within the given radius of the previous shell, or the substructure if it is the first shell.

Picking Additional Shell Atoms: To add individual atoms to a shell, select the desired shell number from the Shells list, select a structural unit from the Pick menu, and pick Workspace atoms. You can also add shell atoms by entering an ASL expression in the Additional atoms for shell section, or by constructing an ASL expression in the Atom Selection dialog box, which you open with the Select button. Atoms may be added even if a substructure has not been defined.

Setting the Force Constant for a Shell: Shell atoms that are to be constrained rather than frozen are allowed to move a finite distance based on the force constant value you supply. Enter the desired force constant in the Force constant text box.

Freezing Atoms in a Shell: To freeze atoms in a shell so that the atoms cannot move at all, select the desired shell number in the Shells list and select Freeze atoms.

Filling Out Residues in a Shell: You can include in the shell all atoms in residues having any members in the original shell definition by clicking Complete residues.

Deleting a Shell: To delete an existing shell, select the shell number in the Shells list and click Delete Shell.

5.5.3 Reading and Writing Substructures

You may want to write your substructure definitions to a separate file so that they can be used in subsequent calculations. To write a substructure file, click Write .sbc File, specify the desired file name, and click OK. Include the .sbc suffix in the name. To read an existing substructure file, click Read .sbc File, navigate to the desired file and select it, then click OK. If you want, you can then edit the substructure to create a new substructure, using the controls in the Substructure folder.

5.6 Setting up the Job

When you are finished with all the settings for the calculation, you can click **Start** to open the Start dialog box (see [Figure 5.4](#)) and choose options for how and where the job is run.

Output Settings

- **Incorporate**—MacroModel job output structures can be incorporated into the project by appending the structural results as new entries (choose **Append new entries**), or by replacing the input structure entries with the output structures and their energetic properties (choose **Replace existing entries**). If you do not want to include the calculation results in the project at all, choose **Do not incorporate**.

Job Settings

- **Name**—By default, when the Current Energy Start dialog box is first opened, the job name `mmod_energy` appears in the Name text box. The name entered in this box is used as the base name for all input and output files associated with that job. To run a job with the default name, simply do not enter a new name in the Name text field. However, because Maestro tracks the files associated with all jobs, you cannot run subsequent jobs with the same name unless you remove the files associated with the first job. If you attempt to run a job with the same name as a job that already exists in the job database, a dialog will appear asking you whether you want to delete the existing job files or rename the new job.
- **Host**—This option menu displays all the hosts defined in the `$SCHRODINGER/schrodinger.hosts` file, with the number of processors on the host in parentheses. To run a MacroModel job on a remote computer, you must first set up the remote hosts according to the instructions in the *Job Control Guide*, then the remote hosts available for running MacroModel jobs appear in the Hosts option menu. Choose a host on which to run the job. The default, `localhost`, runs the job on the local host.




Figure 5.4. The Current Energy Start dialog box.

- **Username**—If you elect to run a job on a remote host, and your user name on that host is different than that on your local machine, enter the correct user name for the remote host.

5.7 Writing Job Files

You may want to set up a MacroModel calculation and save the generated input files without actually submitting the job. This may be the case, for example, if you want to use MacroModel from the command line. To create the structure file and the command (.com) file needed to launch the job at a later time, specify the desired values using the settings on the desired MacroModel panel, and then click the Write button, located in the lower right corner of the panel.

5.8 Running and Monitoring Jobs

To start a MacroModel calculation after you have selected the desired values from the panel settings, click the Start button, located in the lower left corner of the panel. The job is started on the host whose name you have selected from the Hosts option menu. For more information on setting up and selecting remote hosts, see the *Job Control Guide*.

After you click the Start button on a MacroModel panel, the Monitor panel is displayed, and the contents of the *jobname.log* file are displayed as it is written. The last job in the job database portion of the Monitor panel is automatically selected.

While jobs are running, the Detach, Pause, Resume, Stop, Kill, and Update buttons are active. When there are no jobs currently running, only the Monitor and Delete buttons are active. You can choose to monitor or delete the currently selected job by clicking the corresponding buttons. To monitor or delete another job in the database, first click on the job name in the list and then on the appropriate button. Jobs, once launched, appear in the job database portion of the Monitor panel until they are explicitly deleted.

By default, Show jobs from current project only is selected. To show other jobs, deselect this option.

Current Energy Calculations

Determining the current energy for a structure is a basic type of calculation in molecular modeling. Setting up a current energy calculation in MacroModel involves many of the same steps as other calculations within MacroModel, so be sure to read this chapter carefully.

In addition to calculating the total energy for the structure, MacroModel can be instructed to produce an `.mmo` file during a current energy calculation. An `.mmo` file contains detailed information on the various contributions to the potential energy and can lead to a number of insights. Maestro's Force Field Viewer is an excellent tool for interactively examining this information.

6.1 The Current Energy Panel

You can use the Maestro Current Energy panel to prepare and submit MacroModel Current Energy calculations. The Current Energy panel contains 4 folders, each of which contains tools and settings related to a different aspect of a Current Energy calculation—Potential, Constraints, Substructure, and ECalc. The first three folders were described in [Chapter 5](#), so only the ECalc folder is described in this chapter.

To open the Current Energy panel, choose Current Energy from the MacroModel submenu of the Applications menu in the main menu bar.

6.2 The ECalc Folder

The Energy Listing option menu allows you to select the amount of information in the energy listing for the current job. Selecting the **None** setting results in the total molecular mechanics energy being listed in the `jobname.log` file. Selecting **Minimal** prints in addition a minimal summary in the Energy Summary file (`jobname.mmo`). Selecting **Complete** lists the complete energy output to the Energy Summary file, including all internal coordinate components. A “complete” energy listing is required to use Maestro's force field viewing tool. You should note, however, that performing calculations with a complete energy listing can create large output files.




Figure 6.1. The ECalc folder of the Current Energy panel.

6.3 Command File Examples

To run MacroModel calculations, a molecular structure file and a command file are required. The molecular structure file contains the structures to be used as input in the calculation. The command file contains the name of the input structure file, the name of the output structure file, and an ordered list of operation codes (*opcodes*) for the calculations.

Once you set up a job in Maestro and click either **Start** or **Write**, Maestro writes out a molecular structure file and a command file. For many types of jobs, command files written this way are complete and adequate, but for some types of jobs you may need to adjust the Maestro-generated command file.

The command files and the log files for the examples given in this section can be found in `$SCHRODINGER/macromodel-vversion/samples/Examples`.

Below is an example of a command file for an energy calculation with solvation. A description of the opcodes and their arguments follows.

```
ecal.mae
ecal-out.mae
MMOD      0      1      0      0      0.0000      0.0000      0.0000      0.0000
FFLD      10      1      0      1      1.0000      0.0000      0.0000      0.0000
EXNB      0      0      0      0      0.0000      0.0000      0.0000      0.0000
BDCO      0      0      0      0      41.5692 99999.0000      0.0000      0.0000
SOLV      3      1      0      0      0.0000      0.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
ELST      1      0      0      0      0.0000      0.0000      0.0000      0.0000
```

MMOD: Creates and updates an intermediate structure file so that structures can be displayed in Maestro as the job progresses.

FFLD: Force field selection. Arg1 denotes the actual force field used in the calculation (in this case MMFF94). Arg2 defines the electrostatic treatment for the calculation. Default (arg=0) is to use the dielectric treatment encoded in the force field, however, in this case a constant dielectric is used due to the use of solvation model 3 (see SOLV below). Arg4 is MMFF94 specific. Arg4=1 defines the MMFF94s version of the force field, which ensures planarity around delocalized sp² nitrogens.

EXNB: Extended non-bonded interaction cut-off. This is set by default when solvation is used in a calculation. The default values for extended cut-off are: 8 Å vdW (arg5); 20 Å charge-charge electrostatic (arg6); 4 Å hydrogen bonding (arg7).

BDCO: Use the Bond Dipole CutOff (BDCO) method for truncating electrostatic interactions. Arg5 and arg6 are used to specify the cutoffs used for charge-dipole and charge-charge interactions, respectively.

SOLV: Specify the implicit solvation treatment to be used in the calculation. Arg1 defines the type of solvation model to be used. Arg1=3 means that the GB/SA solvation model will be used. Arg2 defines what type of solvent is used. Arg2=1 selects water as the solvent.

READ: Directs MacroModel to read the input file.

ELST: Calculates the single point energy of the input structure(s). Arg1 determines the extent of output listing, and to which files the output will be written. Arg1=1 gives a listing of the total molecular mechanics energy to the log file and a complete listing with all internal coordinates to the .mmo file. Arg2 allows for switching between kJ/mol and kcal/mol units for the energies listed.




Figure 6.2. The Force Field Viewer *panel*.

6.4 The Force Field Viewer

If you want to view details of the force field used in a MacroModel calculation, you can use the Force Field Viewer panel in Maestro. You first ensure that Complete is selected from the Energy Listing option menu before you run the job. The data displayed by the Force Field Viewer is located in .mmo files, which are created when energy operations are performed using Complete energy listings. To read an .mmo file, click Browse and navigate to the desired file.

The Force Field Viewer panel has a number of buttons that represent the various types of interactions that can be present in a calculation (stretch, bond angle etc.). After an .mmo file has been read, buttons that correspond to types of interactions present in that particular file become active. Each button, when selected, displays a corresponding panel. In each panel, the interactions are listed, and selection of an interaction in the list marks the interaction in the displayed structure with a magnifying glass icon and a line or an asterisk. By default the first item in the list is selected. These panels are described in detail in the online help.

The Wilson Angle button and the Improper Torsion buttons are never active simultaneously, since they perform similar operations, but are associated with different force fields. Wilson angles are used in MMFF and MMFFs force field calculations, and improper torsions are used in the others.

A brief summary of each button and its corresponding panel is given below.

Stretch

The Stretch panel is used to display the stretching interactions for a given calculation. The selected interaction is marked with a dashed yellow line.

Bond Angle

The Bend panel, opened by clicking Bond Angle, is used to display the bending interactions for a given calculation. The selected interaction is marked with a dashed red line.

Torsion Angle

The Torsion Angle panel displays the torsion interactions for a calculation. The selected interaction is marked with a pale green dashed line.

Improper Torsion

In the Improper Torsion panel, the improper torsion interactions associated with a particular calculation can be viewed. The selected improper torsion in the displayed structure is marked with a gold-colored asterisk.

Wilson Angle

The Wilson Angle panel displays the Wilson angle interactions for a given calculation. The selected Wilson angle is marked with a peach-colored asterisk.

GB Solvation

The GB Solvation panel can be used to view the Generalized Born (GB) portion of the solvation interactions for a given calculation. The selected atom in the displayed structure is marked with a yellow asterisk.

SA Solvation

The SA Solvation panel can be used to view the Surface Area (SA) solvation interaction components for a given calculation. The selected entry in the SA Solvation List is marked with a peach-colored line.

Van der Waals

The Van der Waals panel can be used to view all van der Waals interactions for a given calculation. The selected interaction in the displayed structure is marked with a dashed orange line.

Electrostatic

The Electrostatic panel can be used to view all electrostatic interactions for a given calculation. The selected interaction is marked with a dashed purple line.

6.5 Checking and Interpreting Results

The energy estimate for an energy calculation depends critically on the nature of the structure it is applied to. It is often crucial that the particular structure was produced using a calculation, such as a minimization or a dynamics simulation, employing the same force field. Otherwise, the energy may contain significant contributions due to strain resulting from differences in the equilibrium structures for the different force fields.

You should also bear in mind the points raised in [Section 3.2 on page 33](#).

Detailed energy listings can be useful, particularly in conjunction with Maestro's Force Field Viewer tool, for the following reasons:

- To gain insight into what interactions are key within the system
- To understand how important parameters of limited quality are for a particular calculation
- To examine problematic structures

Minimizations

Energy minimization is a key type of calculation in molecular modeling, in part because it gives a distinct structure that is often related to a subset of conformers found under thermal conditions. It also plays an important role in many compound calculations, such as conformational searching, and thus this chapter is a useful one to review carefully. MacroModel provides a number of well-tested and efficient minimization methods.

7.1 The Minimization Panel

The MacroModel Minimization panel is used to set up and submit minimization calculations from within Maestro. The Minimization panel consists of five parts. The first part, the upper portion of the panel, contains controls for general aspects of job set up, such as job name and job source. This portion of the panel also appears on the other MacroModel energy panels. The Minimization panel contains four tabbed folders, the first three of which also appear in the other MacroModel panels. The controls in these folders, and those in the upper portion of the panel, are discussed in detail in [Section 5.1 on page 53](#) through [Section 5.5 on page 61](#).

To open the Minimization panel, choose Minimization from the MacroModel submenu of the Applications menu in the main menu bar.

7.2 Performing a Minimization Calculation

To set up a minimization calculation:

1. Select the entry that you want to use as input, or display the structure you want to use in the Workspace.
2. Open the Minimization panel
3. Set the controls in the upper portion of the panel and in the first three folders—Potential, Constraints, and Substructure.
4. In the Mini folder, select a minimization method, the maximum number of iterations, and the convergence criterion and threshold.

The features of the Mini folder are described in detail in the next few sections.

7.2.1 Minimization Methods

From the Method option menu, select a minimization method for your calculation. Choose from the following supported methods:

- **PRCG (Polak-Ribiere Conjugate Gradient)**

This is a conjugate gradient minimization scheme that uses the Polak-Ribiere first derivative method with restarts every 3N iterations [17]. This is the best general method for energy minimization, but it should not be used to find transition states.

- **TNCG (Truncated Newton Conjugate Gradient)**

TNCG uses second derivatives and line searching and is highly efficient for producing very low gradient structures. It generally converges in one tenth the number of iterations necessary for a PRCG, but each iteration takes more time. Often FMNR re-minimization of TNCG structures gives the lowest final gradients [18].



Figure 7.1. The Mini folder of the Minimization panel.

- **OSVM (Oren-Spedicato Variable Metric)**

This is a variable metric minimization that uses the Oren-Spedicato modification [19] of the Fletcher-Powell quasi-Newton method. Convergence to saddle points is possible. Typical convergence occurs in $3N - 6N$ iterations, but note that iteration speed is relatively slow. OSVM is not recommended for structures with poor starting geometries.

- **SD (Steepest Descent)**

This is a steepest descent minimization method. The “SD” should not be used to find saddle points, and convergence is poor towards the end of minimization. This is a good method for starting geometries that are far from the minimum, but switching to another method is recommended when derivatives fall below 10 kJ/mol-Å or so. It is generally not optimal to use SD minimization. PRCG is usually a better choice.

- **FMNR (Full Matrix Newton Raphson)**

With this method, convergence to saddle points is not uncommon. Use FMNR with pre-minimized structures having RMS gradients of less than 0.1 kJ/mol-Å. Use line searching (select after prompt) for problematic cases, or if the RMS First Derivative is greater than 0.1. The preminimization requirement derives from the Newton Raphson assumption of a quadratic potential surface. The method works only if the assumption is valid. FMNR is the most effective method for fully converging structures, but the computational resources required are significant for large structures.

To find saddle point structures, you will need to start very close to the saddle point, and you will have to disable line searching.

FMNR has excellent convergence properties, and typically converges in two to ten iterations.

- **LBFGS (Low-memory Broyden-Fletcher-Goldfarb-Shanno)**

A method that performs well with large structures.

7.2.2 Convergence Parameters

You can specify convergence criteria for your calculation using the Converge on option menu in the Mini folder. Choose from Gradient, Energy, Movement, or Nothing.

Gradient is the default value. With the Energy setting, MacroModel determines when to halt the calculation based on the energy difference between iterations. With the Movement setting, MacroModel determines when to stop the calculation based on the maximum atomic movement at each iteration. When the Nothing option is selected, the calculation runs for the maximum number of iterations as specified in the Max # Iterations text field.

The value in the Maximum iterations text box determines when MacroModel should end the calculation if the specified convergence criterion hasn't been met. The default value is 50 for SD and FMNR minimization methods, and 500 for the remaining methods. To specify a different number of iterations, enter it into the Maximum iterations text box.

The value entered in the Convergence threshold text box determines the threshold applied to the Converge on method. The default convergence setting is convergence on Gradient with a threshold of 0.05. These settings will suffice in most cases.

7.3 Command File Examples

To run MacroModel calculations, a molecular structure file and a command file are required. The molecular structure file contains the structures to be used as input in the calculation. The command file contains the name of the input structure file, the name of the output structure file, and an ordered list of operation codes (*opcodes*) for the calculations.

Once you set up a job in Maestro and click either Start or Write, Maestro writes out a molecular structure file and a command file. For many types of jobs, command files written this way are complete and adequate, but for some types of jobs you may need to adjust the Maestro-generated command file.

The command files and the log files for the examples given in this section can be found in \$SCHRODINGER/macromodel-vversion/samples/Examples.

Below is an example command file for a minimization calculation. Explanations of the opcodes that appear in the file follow.

```
mini-constr.mae
mini-constr-out.mae
MMOD      0      1      0      0      0.0000      0.0000      0.0000      0.0000
FFLD      10      1      0      1      1.0000      0.0000      0.0000      0.0000
BDCC      0      0      0      0      41.5692 99999.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
FXTA      14      15      16      17      100.0000 -44.7270      5.0000      0.0000
CONV      2      0      0      0      0.0500      0.0000      0.0000      0.0000
MINI      9      0      500      0      0.0000      0.0000      0.0000      0.0000
```

MMOD: Creates and updates an intermediate structure file so that structures can be displayed in Maestro as the job progresses.

FFLD: Force field selection. Arg1 denotes the actual force field used in the calculation (in this case MMFF94). Arg2 defines the electrostatic treatment for the calculation. Default (arg=0) is to use the dielectric treatment encoded in the force field. However, in this case a constant dielectric is used due to the use of solvation model 3 (see SOLV below). Arg4 is MMFF94-

specific. Arg4=1 defines the MMFF94s version of the force field, ensuring planarity around delocalized sp² nitrogens.

BDCO: Use the Bond Dipole CutOff (BDCO) method for truncating electrostatic interactions. Arg5 and arg6 are used to specify the cutoffs used for charge-dipole and charge-charge interactions, respectively.

READ: Directs MacroModel to read the input file.

FXTA: Constrains a torsional angle. Arg1–arg4 give the atom numbers of the four atoms defining the dihedral angle that will be constrained. The angle is kept in place through a force constant, defined in arg5 (kJ/mol). Arg6 lists the desired value for the fixed torsion, and arg7 gives the half-width of the flat-bottom potential used (i.e., the “flexibility” of the constrained angle).

CONV: Defines convergence criteria. Arg1=2 signifies derivative convergence (if no CONV command is present, default criterion is 0.05 kJ/mol-Å; this value is set in arg5).

MINI: Starts the minimization. Arg1 defines the type of minimization algorithm to be used. Arg1=9 means that Truncated Newton-Raphson Conjugate Gradient will be used. In arg3, the number of minimization steps is defined. Arg3 can be set to a large number since the calculation stops automatically as soon as the convergence criterion has been reached.

7.4 Checking and Interpreting Results

Any individual minimization is intended to minimize the structure in the local minimum which may not be the global minimum of the system. This fact can make it hard to put the minimized system in the proper context without additional calculations (e.g., conformational searches). Typical output for a MacroModel minimization looks like this:

```
Starting conjugate gradient minimization.
Minimization converged; gradient = 0.387E-01 .LT. 0.500E-01
Iterations = 110 out of 500
Conf 1 E = -594.116 ( 0.039) kJ/mol
BatchMin normal termination
Total number of structures processed = 1
BatchMin: normal termination 20-Feb-2002
15:31:29
```

The line:

```
Minimization converged; gradient = 0.387E-01 .LT. 0.500E-01
```

indicates that the minimization met the convergence conditions. In this case the gradient attained, 0.0387 kJ/mol-Å, was less than the level required, 0.05 kJ/mol-Å.

If the calculation had not converged, the minimization converged line would be missing and the number of iterations would match the total specified.

Unconverged results typically have little or no value. If the minimization does not converge within the specified number of iterations, you may need to repeat or continue the minimization with more iterations.

Most MacroModel calculations consider minimizations converged when the RMS gradient of the energy is less than 0.05 kJ/mol-Å. While this is adequate for most types of calculations, this value may be modified or another property, such as atomic motion, may be used to detect when a minimization is sufficiently converged.

The results of a minimization are force-field-dependent: minimizing the same initial structure using different force fields usually results in similar but not identical structures.

Multiple Minimizations

Multiple minimization is a tool for automatically minimizing a collection of structures. The collection of structures can involve different molecules or conformers. If the structures are conformers, then it is possible to do the minimization in conjunction with the elimination of redundant conformers. MacroModel uses multiple minimization a number of ways, some of which are described in this chapter.

8.1 The Multiple Minimization Panel

The Multiple Minimization panel is used to set up and submit minimization jobs that use as input either a Maestro- or MacroModel-formatted file containing the structures to be minimized. If the structures are conformers, then redundant conformers are eliminated if comparison atoms are specified. The Multiple Minimization panel can also be used to perform partition coefficient estimation.

The Multiple Minimization panel consists of several parts, four of which—the upper portion, Potential, Constraints, and Substructure folders—are common to all MacroModel energy panels. These components are described in detail in [Section 5.1 on page 53](#) through [Section 5.5 on page 61](#). The Mini folder is common to all of the MacroModel panels except Current Energy. For an explanation of the controls in this folder, see [Section 7.1 on page 73](#). The Mult folder, which includes access to the Comparison Atoms panel, is unique to the Multiple Minimization panel.

To open the Multiple Minimization panel, choose Multiple Minimization from the MacroModel submenu of the Applications menu on the main menu bar.

8.2 Setting Up Multiple Minimization Calculations

In the Multiple Minimization panel, set the options in the upper portion of the panel and in the Potential, Constraints, and Substructure folders. See [Section 5.1](#) through [Section 5.5](#) for a detailed discussion of these settings. Then set the options in the Mini folder as described in [Section 7.2 on page 73](#).

With these portions of the job setup complete, open the Mult folder. The Mult folder contains the controls for setting the parameters of the multiple minimization calculation. These controls are described in the sections below. Once you have adjusted the settings in the Mult folder,




Figure 8.1. The Mult folder of the Multiple Minimization panel.

either click **Start** to set up and launch the job, or click **Write** to write the job commands to a file to be run or edited later.

Input file

In this text box you can enter the name of an input file. This must be a file in either Maestro or MacroModel format that contains one or more valid structures. Often these structures are the results from a previous conformational search, which are to be reminimized. If no file is entered, the structural input is taken from the source indicated under **Source of job input** in the upper portion of the panel. The **Open** button displays a file selector that can be used to browse the file system for the desired input file.

Energy window for saving structures

If a structure is minimized and found to be above the global minimum by more than the energy window value (in of kJ/mol), it will be rejected and not included in the output file. This value has no effect if the structures being minimized are not conformers.

Comparison Atoms

During multiple minimizations of conformers, redundant structures are eliminated if comparison atoms are specified. The structures are compared against other low-energy structures that have already been minimized. The comparison is performed by rigid superposition, comparing only those atoms specified as “comparison atoms” in the setup, taking into account the topological symmetry of the molecule.

To define comparison atoms, click the Comparison Atoms button and use one or more of the following options:

- Click Heavy Atoms + O-H,S-H

This adds all the non-hydrogen atoms and the hydrogen atoms attached to oxygen and sulfur to the list of comparison atoms.

- Click Heavy Atom

This adds only the non-hydrogen atoms to the list of comparison atoms.

- Pick atoms in the Workspace

Choose a structural unit from the Pick menu in the Define comparison atoms section, and pick atoms in the Workspace to add the atoms in the structural unit to the list of comparison atoms.

- Select atoms using the Atom Selection dialog box




Figure 8.2. The Comparison Atoms panel.

For more complex combinations of comparison atoms, you can select atoms using the Atom Selection dialog box to build an Atom Specification Language (ASL) expression for the required atoms. See [Section 5.3](#) of the *Maestro User Manual* for more information on the Atom Selection dialog box.

- Select all atoms

Click All to add all the atoms in a structure to the list of comparison atoms.

When comparison atoms are picked, Maestro places light green “=” labels on them. To distinguish the currently selected atom, the program colors its label turquoise. To hide these markers, clear the Show Markers button.

To delete a defined comparison atom, select it in the list of comparison atoms, then click Delete. To redefine a comparison atom, pick a new atom while the comparison atom is selected in the list. To delete all the defined comparison atoms, click Delete All.

Max distance between atoms in equal structures

The threshold for determining whether structures should be considered to be equivalent is set by default to 0.25 Å. When the structures are compared, the maximum distance between pairs of corresponding atoms must be less than this threshold for the structures to be considered equivalent. You can change the threshold in this text box.

8.3 Partition Coefficient Estimation

Multiple minimization calculations can also be used to obtain estimates of the partition coefficients of a set of solutes between two solvents. After you have set up other parameters for a multiple minimization, including a primary solvent, select Enable LogP in the Mult folder and choose the secondary solvent from the Secondary Solvent menu.

Note: Once you select Enable logP, the options above it are unavailable, because they refer to redundant conformer elimination, which is not relevant to partition coefficient estimation.

Each molecule in the input file is minimized twice, once in each solvent. The results of each minimization are used to evaluate the free energy difference for the two solvents. Hence the logarithm of the partition coefficient, $\log P_{\text{solvl},\text{solv2}}$, which is defined by the relationship

$$\log P_{\text{solvl},\text{solv2}} = (\Delta G_{\text{solv2}} - \Delta G_{\text{solvl}})/(2.30RT)$$

where $P_{\text{solv1,solv2}} = [\text{Solute}]_{\text{solv1}} / [\text{Solute}]_{\text{solv2}}$, ΔG_{solvN} is the free energy of solvation of the molecule in solvent N , R is the gas constant, and T is the temperature in kelvin. The results are collected in a table at the end of the log file.

The solvation models are parametrized for ambient conditions. If the temperature deviates significantly from these conditions the solvation energy estimates and hence the calculated $\log P_{\text{solv1,solv2}}$ become less reliable.

8.4 Automatic Setup (AUTO)

You can use automatic setup for multiple minimizations of a collection of conformers with redundant conformer elimination. Automatic setup is enabled by adding the [AUTO](#) opcode to the command file for a calculation. You should place the [AUTO](#) opcode, with appropriate arguments, just before the [MINI](#) opcode in the command file. You must also add a [MULT](#) opcode to the command file, as seen in the example in [Section 8.6.2 on page 86](#). [AUTO](#) is available from Maestro only for conformational searches. Automatic setup should not be used with minimizations of single structures or minimizations with multiple non-conformers as input.

[AUTO](#) selects comparison atoms ([COMP](#)), so that redundant conformers can be identified and removed without explicitly designating comparison atoms in the command file or from Maestro. In addition, chiral atoms ([CHIG](#)) and torsion checks ([TORC](#)) are automatically identified by [AUTO](#). [AUTO](#) applies the desired structural criteria to the minimized conformers as if the [CHIG](#) and [TORC](#) opcodes were present in the command file. [AUTO](#) has options to adjust or turn off these features. See the [AUTO](#) opcode in the *MacroModel Reference Manual* for more information.

In addition, [AUTO](#) uses the information in any substructure specified for the computation, including information on fixed or frozen atoms. The substructure must be specified in a substructure file and read in with a [SUBS](#) 0 command, and cannot be specified in the command file. Likewise, fixed or frozen atoms must be specified in the substructure and cannot be specified independently.

Finally, if only automatic comparison atom setup is required for a multi-conformer multiple minimization with redundant conformer elimination, then [COMP](#) arg1=0 may also be used in the command file. Examples are given in [Section 8.6.2 on page 86](#).

[AUTO](#) is not currently compatible with [LOOP](#), and must be used carefully with eMBrAcE. For more information on the parameters that [AUTO](#) is assembling, use [DEBUG](#) 520 or 521.

8.5 Preparation of Multi-Ligand Structure Files With premin

MacroModel's multiple minimization capabilities are used in a script named `premin`, which is designed to prepare multi-ligand structure files for use in Glide and other applications. Database files sometimes contain problematic ligand structures that are either chemically incorrect or have species that are not covered by the force field parameters. This script culls out such problematic structures and minimizes the energy of all other structures. The problematic structures and the minimized structures are saved in separate files. You can also filter out problematic structures without doing a preminimization.

The preminimization employs a $4r$ distance-dependent dielectric model and uses MacroModel's efficient truncated-Newton minimizer. As well as the recommended default force field, MMFFs, `premin` also supports the OPLS_2001 and OPLS_2005 force fields.

The syntax for `premin` is shown below:

```
$SCHRODINGER/utilities/premin [options] input.mae
```

where `input.mae` is the input file that contains the ligand structures. The options are described in [Table 8.1](#).

Table 8.1. Options for the premin script.

Option	Description
-h	Print help text and exit.
-lic ligprep	Use LigPrep license instead of MacroModel license
-s output	Output file for successfully processed structures. Default is <code>input-min.mae</code> for minimizations, and <code>input-filtered.mae</code> for filtering.
-u bad	Output file of structures that could not be minimized. Default is <code>input-bad.mae</code> .
-v	Print version number and exit.
-doc	Print a brief guide to using <code>premin</code> .
-f 11	Use OPLS_2001 instead of the default force field, MMFFs.
-f 14	Use OPLS_2005 instead of the default force field, MMFFs.
-custom	Use custom <code>filter.com</code> and <code>goodmin.com</code> files from the current working directory.
-filter -m filter	Filtering only: skip minimization.

By default, successfully minimized structures are saved in *input-min.mae*, and unsuccessful structures are saved in *input-bad.mae*. If you choose to filter out problematic structures without minimization (*-filter* or *-m filter*), successfully filtered structures are saved in the file *input-filtered.mae*, and problematic structures are again saved in *input-bad.mae*.

The script creates two MacroModel command files, *filter.com* and *goodmin.com*, and writes them to the current working directory. You can then customize these files and use them in place of the default files with the *-custom* option. If you perform filtering only, the file *goodmin.com* is not used.

This script uses the SPAT opcode to filter out problematic structures that contain metal ions or generalized atom types.

8.6 Command File Examples

To run MacroModel calculations, a molecular structure file and a command file are required. The molecular structure file contains the structures to be used as input in the calculation. The command file contains the name of the input structure file, the name of the output structure file, and an ordered list of operation codes (*opcodes*) for the calculations.

Once you set up a job in Maestro and click either **Start** or **Write**, Maestro writes out a molecular structure file and a command file. For many types of jobs, command files written this way are complete and adequate, but for some types of jobs you may need to adjust the Maestro-generated command file.

The command files and the log files for the examples given in this section can be found in `$SCHRODINGER/macromodel-vversion/samples/Examples`.

8.6.1 Energy Minimization of Multiple Non-Conformers

Below is an example of the command file for an energy minimization of multiple non-conformers and explanations of the opcodes used in the file.

```
mult-min.mae
mult-min-out.mae
MMOD      0      1      0      0      0.0000      0.0000      0.0000      0.0000
FFLD      10      1      0      1      1.0000      0.0000      0.0000      0.0000
BDCC      0      0      0      0      41.5692 99999.0000      0.0000      0.0000
DEMX      0      0      0      0      50.0000      0.0000      0.0000      0.0000
BGIN      0      0      0      0      0.0000      0.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
CONV      2      0      0      0      0.0500      0.0000      0.0000      0.0000
MINI      9      0      500      0      0.0000      0.0000      0.0000      0.0000
END        0      0      0      0      0.0000      0.0000      0.0000      0.0000
```

MMOD: Creates and updates an intermediate structure file so that structures can be displayed in Maestro as the job progresses.

FFLD: Force field selection. Arg1 denotes the actual force field used in the calculation (in this case, MMFF94). Arg2 defines the electrostatic treatment for the calculation. The default (arg=0) is to use the dielectric treatment encoded in the force field. Arg4 is MMFF94-specific. Arg4=1 defines the MMFF94s version of the force field, ensuring planarity around delocalized sp^2 nitrogens.

BDCO: Use the Bond Dipole CutOff (BDCO) method for truncating electrostatic interactions. Arg5 and arg6 specify the cutoffs used for charge-dipole and charge-charge interactions, respectively.

DEMX: This command prevents saving of high-energy conformers during the search. Arg5 defines the allowed energy window above the current global minimum. New conformers that are not within arg5 kJ/mol are discarded. Additionally, a preliminary energy test can be performed during the energy minimization to ensure that a reasonable structure has been found. Arg2 sets the number of energy iterations to perform before the preliminary test (a good value is approximately 1/3 of the total number of energy iterations), while arg6 defines the energy above which conformers are discarded (a value of about 1.5–2 times arg5 is recommended).

BGIN/END: Defines the start/end of a loop. All commands between the BGIN and END lines are performed for each structure in the input file.

READ: Directs MacroModel to read the input file.

CONV: Defines convergence criteria. Arg1=2 signifies derivative convergence. The default criterion, if no **CONV** command is present, is 0.05 kJ/mol-Å; this value is set in arg5.

MINI: Starts the minimization. Arg1 defines the type of minimization algorithm to be used. Arg1=9 means that Truncated Newton-Raphson Conjugate Gradient will be used. In arg3, the number of minimization steps is defined. Arg3 can be set to a large number since the calculation will automatically stop as soon as the convergence criterion is reached.

8.6.2 Multiple Conformer Minimization With Automatic Redundant Conformer Elimination

Below are two examples of command files for energy minimization of multiple conformers with redundant conformer elimination. The first uses the **AUTO** opcode; the second uses the **COMP** opcode with arg1=0. Explanations of the opcodes used in the file that are essential to this example are given.

```
mult_conformers.mae
```

```

mult_auto-out.mae
MMOD      0      1      0      0      0.0000      0.0000      0.0000      0.0000
FFLD     10      1      0      1      1.0000      0.0000      0.0000      0.0000
SOLV      3      1      0      0      0.0000      0.0000      0.0000      0.0000
EXNB      0      0      0      0      0.0000      0.0000      0.0000      0.0000
BDCCO     0      0      0      0      89.4427 99999.0000      0.0000      0.0000
CRMS      0      0      0      0      0.0000      0.2400      0.0000      0.0000
DEMX      0      0      0      0      51.0000      0.0000      0.0000      0.0000
MULT      0      0      0      0      0.0000      0.0000      0.0000      0.0000
BGIN      0      0      0      0      0.0000      0.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
CONV      2      0      0      0      0.0100      0.0000      0.0000      0.0000
MSYM      0      0      0      0      0.0000      0.0000      0.0000      0.0000
AUTO      1      0      0      0      -1.0000      0.0000      0.0000      0.0000
MINI      1      0      5000      0      0.0000      0.0000      0.0000      0.0000
END        0      0      0      0      0.0000      0.0000      0.0000      0.0000

```

MULT: Mandatory for this computation, due to the presence of the **AUTO** opcode.

AUTO: Automatic parameter assignment. Arg1 turns on **AUTO** only for the first conformer. Arg5 turns off torsional parameters, which are not needed for this minimization.

```

mult_conformers.mae
mult_comp-out.mae
MMOD      0      1      0      0      0.0000      0.0000      0.0000      0.0000
FFLD     10      1      0      1      1.0000      0.0000      0.0000      0.0000
SOLV      3      1      0      0      0.0000      0.0000      0.0000      0.0000
EXNB      0      0      0      0      0.0000      0.0000      0.0000      0.0000
BDCCO     0      0      0      0      89.4427 99999.0000      0.0000      0.0000
CRMS      0      0      0      0      0.0000      0.2400      0.0000      0.0000
DEMX      0      0      0      0      51.0000      0.0000      0.0000      0.0000
BGIN      0      0      0      0      0.0000      0.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
CONV      2      0      0      0      0.0100      0.0000      0.0000      0.0000
MSYM      0      0      0      0      0.0000      0.0000      0.0000      0.0000
COMP      0      0      0      0      0.0000      0.0000      0.0000      0.0000
MINI      1      0      5000      0      0.0000      0.0000      0.0000      0.0000
END        0      0      0      0      0.0000      0.0000      0.0000      0.0000

```

COMP: Arg1=0 specifies that all heavy atoms are to be used as the comparison atoms when redundant conformer elimination is performed.

8.6.3 Partition Coefficient Estimation

Below is an example of a command file for conducting $\log P_{\text{octanol,water}}$ calculations on a series of molecules. It is very similar to that used for multiple minimization of non-conformers.

```

logP.mae
logP-out.mae

```

FFLD	11	1	0	1	1.0000	0.0000	0.0000	0.0000
BDCO	0	0	0	0	41.5692	99999.0000	0.0000	0.0000
SOLV	3	9	0	0	0.0000	0.0000	0.0000	0.0000
LOGP	1	1	1	0	0.0000	0.0000	0.0000	0.0000
BGIN	0	0	0	0	0.0000	0.0000	0.0000	0.0000
READ	0	0	0	0	0.0000	0.0000	0.0000	0.0000
AUTO	0	-1	0	0	0.0000	1.0000	0.0000	0.0000
CONV	2	0	0	0	0.0500	0.0000	0.0000	0.0000
MINI	1	0	5000	0	0.0000	0.0000	0.0000	0.0000
END	0	0	0	0	0.0000	0.0000	0.0000	0.0000

SOLV: Specifies the first, or primary, solvent. Arg2 = 9 corresponds to octanol.

LOGP: Turns on partition coefficient estimation and specifies the second solvent. Arg2 = 1 corresponds to water.

AUTO: Necessary for the proper processing of LOGP calculations. Arg2 = -1 instructs **AUTO** not to make a list of comparison atoms, and arg6 = 1.0000 indicates that this is a serial calculation.

In the .log file produced when such a calculation is run, a report appears for each ligand processed:

```
=====
LOGP(octanol,water) Calculation for ligand # 1
Ligand Name: Case74C
LOGP(octanol,water) = 1.41024
=====
```

and a table appears at the end of the .log file:

```
=====
LOGP(octanol,water) Results

Ligand      Solvation Free Energies      LOGP(octanol,water)  Converged
              octanol          water
              kJ/mol.          kJ/mol.
      1      -26.99          -18.94          1.41          T
      2      -20.93          -12.48          1.48          T
      3      -20.02          -11.02          1.58          T
      4      -29.41          -16.19          2.32          T
      ...
     12      -16.84          -15.40          0.25          T
=====
```

Here, T or F in the Converged column signifies whether or not the minimizations (one in each solvent for each molecule) converged. If they did not converge, the results should not be considered reliable.

Geometry Scans

Geometry scans can be used to probe the energy of a molecule as a function of a small number of geometric parameters: bond lengths, bond angles, or dihedral angles. The process involves stepping through the specified parameters, restraining the parameters and performing a minimization on the remaining geometric parameters. For example dihedral angle scans are used to map out the potential energy for one or two dihedral angles, such as a Ramachandran ϕ - ψ plot for amino acids, and can be very useful in understanding conformational behavior. Support for geometry scans in Maestro is limited to dihedral angle scans (formerly known as dihedral drives). Distance and bond angle scans can be done from the command line—see [Section 3.7](#) of the *MacroModel Reference Manual* for more information.

9.1 Performing a Dihedral Scan in Maestro

You can use the MacroModel Dihedral Driving panel to set up and submit dihedral scan calculations on either one or two dihedral angles. The tools in the Drive folder allow you to define the dihedrals for the calculation, to set the rotation angle through which the dihedrals are to be driven, and to specify the number of degrees you want to move through with each increment.

All of the components of the Dihedral Driving panel, with the exception of the Drive folder, appear on other MacroModel energy panels. For information about the upper portion of the Dihedral Driving panel, the Potential, Constraints, and Substructure folders, see [Section 5.1 on page 53](#) through [Section 5.5 on page 61](#). For information about the controls in the Mini folder, see [Section 7.1 on page 73](#).

To open the Dihedral Driving panel, choose Dihedral Driving from the MacroModel submenu of the Applications menu on the main menu bar.

To set up the first dihedral, choose either Atom or Bond from the Pick menu, and then pick the atoms or bonds in the Workspace to define the dihedral angle. After you make the selection, the atom numbers appear in the list at the top of the folder. Repeat the picking operations to define a second dihedral.

To set up the initial and final angles and the angle increment for the driving calculation, select a dihedral entry from the list at the top of the folder, then enter values in the Start, Final, and Increment text boxes. The default start value is 0°, and the default final value is 360°.




Figure 9.1. The Drive folder of the Dihedral Driving panel.

The value entered in the Increment text field determines the number of degrees through which a dihedral is driven before a minimization is done. An increment can be negative, but only if the starting angle is larger than the final one. The minimum increment is 1°, but an increment of 10° is usually sufficient. The default value is 30°.

To delete a defined dihedral, select it in the list of dihedrals and click **Delete**. To delete all defined dihedrals, click **Delete All**.

9.2 Command File Examples

To run MacroModel calculations, a molecular structure file and a command file are required. The molecular structure file contains the structures to be used as input in the calculation. The command file contains the name of the input structure file, the name of the output structure file, and an ordered list of operation codes (*opcodes*) for the calculations.

Once you set up a job in Maestro and click either Start or Write, Maestro writes out a molecular structure file and a command file. For many types of jobs, command files written this way are complete and adequate, but for some types of jobs you may need to adjust the Maestro-generated command file.

The command files and the log files for the examples given in this section can be found in `$SCHRODINGER/macromodel-vversion/samples/Examples`.

Below is an example command file for a dihedral drive calculation and explanations of the opcodes used in the files.

```
ddrive.mae
ddrive-out.mae
MMOD      0      1      0      0      0.0000      0.0000      0.0000      0.0000
FFLD      10      1      0      1      1.0000      0.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
BGIN      0      0      0      0      0.0000      0.0000      0.0000      0.0000
DRIV      1      14     15     16     0.0000     360.0000     30.0000     0.0000
DRIV      6      25     26     27     0.0000     360.0000     30.0000     0.0000
CONV      2      0      0      0      0.0500      0.0000      0.0000      0.0000
MINI      9      0      500     0      0.0000      0.0000      0.0000      0.0000
END       0      0      0      0      0.0000      0.0000      0.0000      0.0000
```

MMOD: Creates and updates an intermediate structure file so that structures can be displayed in Maestro as the job progresses.

FFLD: Force field selection. Arg1 denotes the actual force field used in the calculation (in this case MMFF94). Arg2 defines the electrostatic treatment for the calculation. The default (arg2=0) is to use the dielectric treatment encoded in the force field. In this case a constant dielectric (arg2=1), with a dielectric constant of 1 explicitly requested (arg5). Arg4 is MMFF94-specific. Arg4=1 defines the MMFF94s version of the force field, ensuring planarity around delocalized sp^2 nitrogens.

READ: Directs MacroModel to read the input file.

BGIN/END: Defines the start/end of a loop. All commands between the BGIN and END lines will be performed for each structure in the input file.

DRIV: Directs MacroModel to carry out a dihedral drive using the indicated dihedral angles (arg5=0° to arg6=360°) and step sizes (arg7=30°).

CONV: Defines convergence criteria. Arg1=2 signifies derivative convergence. The default criterion, if no CONV command is present, is 0.05kJ/mol-Å; this value is set in arg5.

MINI: Starts the minimization. Arg1 defines the type of minimization algorithm to be used. Arg1=9 means that Truncated Newton-Raphson Conjugate Gradient will be used. In arg3, the

number of minimization steps is defined. Arg3 can be set to a large number since the calculation automatically stops as soon as the convergence criterion is reached.

9.3 Plotting Scan Results in Maestro

You can display the results of a geometry scan in the 1DPlot panel or the 2DPlot panel in Maestro, depending on how many parameters you chose. These panels were designed for dihedral scans, but can be used for any kind of scan. The panels and their use is described below. Since there are many common features of these panels, the general features are described first. The information given below is also available in the online help.

To open the 1D Plot panel or the 2D Plot panel, choose 1D Plot or 2D Plot from the Tools menu in the main window.

9.3.1 1D and 2D Plot Panel General Features

Each panel consists of a plot area on the right and a set of controls on the left. The common controls are described below; the controls specific to a given plot are described in later sections.

Minimum Energy, Maximum Energy

These text boxes allow you to set the energy range for the plot.

Full Scale

Restores a plot whose range has been altered to its full scale.

Energy Units

These options allow you to display data in either kJ/mol or kcal/mol energy units.

Energy Scale

These options allow you to display the energy as either an absolute or a relative value. The Relative option is useful for estimating rotational barriers.

Constrain to Square

This option forces the plot to be displayed as a square, regardless of how the panel is resized.

Angle, Energy

These text boxes (Angle 1 and Angle 2 for 2D plots) display the values of the angle or angles and the energy at the pointer position when you middle-click and hold in the plot area. In addi-

tion, the structure corresponding to the given angles is displayed in the Workspace. (For distance scans, the “angle” is the distance.)

PostScript

Opens the 1D Plot PostScript panel or 2D Plot PostScript panel, which allow you to save a PostScript image of the plot area. When you have entered the file name, and selected a paper size and orientation, click Write to write the file.

9.3.2 The 1D Plot Panel

After a scan calculation with a single geometric parameter has been successfully completed, the resulting .grd file can be displayed graphically using this panel. To open a .grd file, click the Open button, navigate to the file, select it, and click Open. The data from the .grd file is displayed in the plot area. The corresponding structures are read into Maestro as a scratch entry, and one of the structures is displayed in the Workspace.




Figure 9.2. The 1D Plot panel.

The appearance of a displayed plot can be changed using the controls on the left side of the panel and the controls in the 1D Data panel, which is opened using the Data Sets button.

- The ranges of the axes can be changed by entering values in the Minimum Angle, Maximum Angle, Minimum Energy and Maximum Energy text boxes.
- The symbol shape and color and the curve style, color, and width are all controlled from the 1D Data panel. To change these attributes, select the data set in the list at the top of the panel, then choose the attributes from the option menus.

You can display multiple plots from `.grd` files in the 1D Plot window. Each plot is added by clicking the Open button and selecting the `.grd` file. The data sets are listed in the 1D Data panel. You can delete a plot by selecting it in the list and clicking Delete.

The plot area itself is an interactive tool as well as a display area. If you middle-click in the plot area, the structure in the Workspace is updated to reflect the angle values that correspond to the pointer position on the plot, and the energy and angle values are displayed in the Energy and Angle text boxes. You can drag horizontally with the middle mouse button, and the structure changes as you move the mouse. Note that the displayed structure's movement is rigid-rotor. The original calculations used to obtain the data for the plot, however, are performed with all other degrees of freedom minimized. Even so, the molecule's movement should give some idea of the geometry at various points on the plot. If you have more than one plot in the plot area, only the last structure added is displayed and updated.

To create a tab-separated file that can be read by a spreadsheet, click Tabular Data. A dialog box opens, in which you can enter a file name for the file. The data for all available plots is written to this file.

9.3.3 2D Plot Panel

After a 2D scan calculation has been successfully completed, the resulting `.grd` file can be displayed as a contour plot using this panel. To open a `.grd` file, click the Browse button, navigate to the file, select it, and click Open. The data from the `.grd` file is displayed in the plot area as a contour map. The corresponding structures are read into Maestro as a scratch entry, and one of the structures is displayed in the Workspace. Only one `.grd` file can be displayed at any given time. Reading in a second file clears the data from the first.

The appearance of a displayed plot can be changed using the controls on the left side of the panel. You can set the number of contour lines in the Number of Contours text box, the line thickness in the Contour Width text box, and use dashed lines for negative contour values by selecting Negative Dashed. You cannot, however, change the range of angles displayed.

To determine the contour intervals, Maestro divides the range of energy values specified in the Minimum Energy and Maximum Energy text boxes (or the maximum and minimum values in




Figure 9.3. The 2D Plot panel.

the file) by $n+1$, where n is the specified number of contour lines. This means that the maximum and minimum energy values are not represented by contour lines.

The plot area itself is an interactive tool as well as a display area. If you middle-click in the plot area, the structure in the Workspace is updated to reflect the angle values that correspond to the pointer position on the plot, and the energy and angle values are displayed in the Energy and Angle text boxes.

9.4 Checking and Interpreting Results

Dihedral angle driving calculations involve minimizations, and you should check that minimizations for all angles converge.

Conformational Searches

A frequent question in molecular modeling studies is, What conformations are important in this system? One way to address this question is to perform a conformational search in order to find a set of low energy conformers. MacroModel excels at conformational searching with a number of state-of-the-art conformational searching methods, described in the sections below. Conformational searches typically cycle through the process of generating a new structure, minimizing it, and then determining if the structure should be retained. Structure retention can be based upon energy relative to that of the lowest energy conformer found so far in the search and redundancy with other generated structures.

10.1 Conformational Search Methods

A variety of conformational search methods is available. The methods available in MacroModel are described below.

Monte Carlo Multiple Minimum (MCMM)

The Monte Carlo method (MCMM) implemented in MacroModel is highly efficient in performing global searching, exploring close as well as distant areas of the potential energy surface (PES). Random changes are made in torsion angles during the search. Although there is no limit to the number of variable torsions allowed in the search, more than about 10 or 15 flexible torsions greatly increases the complexity of the search. You should therefore expect increased searching times, and possibly also non-convergence of the search. See References [20](#) and [21](#) for more information.

Systematic Pseudo-Monte Carlo (SPMC)

If the task is to locate all possible conformational minima on the surface, MacroModel's systematic method (SPMC) ensures the exploration of unusual torsional values. SPMC allows for high efficiency throughout the search, even at the end where other methods tend to revisit already found conformers. The SPMC searching method is most efficient for smaller molecules. See Reference [22](#) for more information.

Low-Mode Conformational Search Methods

If you have little or no prior knowledge of the system to be searched, the Low-Mode Conformational Search (LMCS) [[23](#)], based on the principles behind saddle point searching, allows

for automatic searching without the need to define parameters to be varied during the search. One important consequence of this is that LMCS allows for a serial search of multiple different molecules without user intervention.

A related method, Large Scale Low-Mode (LLMOD) [24, 25], has been developed for large-scale conformational searching, such as protein loop optimization, homology model refinement, and fully flexible docking for induced fit modeling. LLMOD is similar to LMOD, but computes low-mode eigenvectors of a Hessian matrix that is referenced only implicitly, through its product with a series of vectors. LLMOD is the first conformational search method that can be applied to fully flexible, unconstrained protein structures. Serial LLMOD calculations are not currently supported.

Mixed MCMM/Low-Mode Conformational Search Methods

A pure LMCS search is local in scope, but performs exceptionally well in global searches through hybridization with MCMM [26]. By defining key torsions to be varied through MCMM steps, the mixed MCMM/LMCS search has proved faster and more efficient than any other searching algorithm for a variety of systems. MCMM/LMCS has successfully been applied to difficult tasks, such as searching the conformational space of a ligand in the active site of a protein. Mixed MCMM/LMC2 searches are also supported.

10.2 Performing Conformational Searches

You can use the Conformational Search panel to set up and submit MacroModel conformational search calculations. The panel has six parts, five of which are common to other MacroModel panels. For more information about the controls in the Mini folder, see [Section 7.1 on page 73](#). For information about the upper portion of the Conformational Search panel and the Potential, Constraints, and Substructure folders, see [Section 5.1 on page 53](#) through [Section 5.5 on page 61](#).

The CSearch folder is unique to the Conformational Search panel. Controls in this folder allow you to choose a search method; define an energy window for saving structures; set up torsion, chiral center, and distance checks; and define other search parameters.

To open the panel, select Conformational Search from the MacroModel menu in the main menu bar. If the MacroModel menu is not in the main menu bar, choose MacroModel from the Applications menu.

To set up a conformational search, you must first select any entries that you want to use as input from the Project Table, or display the structure you want to use in the Workspace. In the Conformational Search panel, set the options in the upper portion of the panel and in the Potential, Constraints, and Substructure folders, then set the options in the Mini folder.

With these portions of the job setup complete, open the CSearch folder. The CSearch folder contains the controls for setting the parameters of the conformational search. These controls are described below. When you have finished adjusting the settings in the CSearch folder, either click **Start** to set up and launch the job, or click **Write** to write the job commands to a file to be run or manually edited later.

10.2.1 Conformational Search Method

You can select the conformational search method to be used during a calculation from the Method list. The conformational search method options are described below.

Torsional sampling (MCMM)

This is the recommended conformational search method. The input structure is modified by random changes in torsion angles and/or molecular position as specified in the panels opened




Figure 10.1. The CSearch folder of the Conformational Search panel.

by the Torsion Rotations or Molecule Trans/Rot buttons. Ordinarily, whether a single structure or multiple structures appear in the input file, they will first be read in, minimized, and treated as if already found by the MCMM procedure. This allows a new search to be initialized from the output of a previous search, by using the output file of the old search as input for the new one. See [page 97](#), and the **MCMM** opcode description in the *MacroModel Reference Manual*, for a detailed description of this method.

Serial torsional sampling (MCMM)

Serial torsional sampling (MCMM) is an automated procedure for performing a separate MCMM search on each structure in the input structure file. Perform automatic setup during calculation is always selected when Serial MCMM is selected. This option causes MacroModel to set up each search during the calculation automatically, and is described in [Section 10.2.4 on page 103](#).

Systematic torsional sampling (SPMC)

This method is similar to MCMM, but uses a systematic search instead of a random search. The search begins at low torsional resolution (120°), searches all angles without duplicating coverage, then doubles the resolution. This method has the advantage of not retracing its path and consequently converges the final stages of the conformational search more efficiently than MCMM. Like MCMM, the method is effectively open-ended: it will search conformational space until you stop it. See [page 97](#), and the **SPMC** opcode command description in the *Macro-Model Reference Manual*, for a detailed description of this method.

Low-mode sampling

This method, termed LMOD, is highly efficient and has the advantage that ring structures and variable torsion angles do not need to be specified. This conformational search method works by exploring the low-frequency eigenvectors of the system, which are expected to follow “soft” degrees of freedom, such as torsions. LMOD methods search conformational space aggressively enough to switch the chirality of atoms within the structures provided. Chirality checking should be used for chiral atoms for which such chirality switching is undesirable (see [Section 10.2.4 on page 103](#)). See [page 97](#), and the **LMCS** opcode description in the *Macro-Model Reference Manual*, for a detailed description of this method.

Serial low-mode sampling

Serial Low-Mode is an automated procedure for performing a separate low-mode search on each structure in the input structure file. Perform automatic setup during calculation is always activated when Serial MCMM/low-mode is selected.

Mixed torsional/Low-mode sampling

This method uses a combination of the random changes in torsion angles and/or molecular position from the MCMM method, together with the low-mode steps from the LMOD method used in pure low-mode. See [page 98](#), and the [MCMM](#) and [LMCS](#) opcode descriptions in the *MacroModel Reference Manual*, for a detailed description of this method.

Serial torsional/Low-mode sampling

Serial MCMM/Low-Mode is an automated procedure for performing a separate mixed MCMM/Low-Mode search on each structure in the input structure file. Perform automatic setup during calculation is always activated when Serial MCMM/Low-mode is selected.

Large scale low-mode sampling

Large Scale Low-Mode (LLMOD) is similar to Low-Mode (LMOD), except that it uses techniques to reduce the memory requirements so that it can be applied to much larger systems such as protein-ligand complexes. Like LMOD, LLMOD methods search conformational space aggressively enough to switch the chirality of atoms within the structures provided. Chirality checking should be used for chiral atoms for which such chirality switching is undesirable (see [Section 10.2.4 on page 103](#)). See the [LMC2](#) and [ARPK](#) opcode descriptions in the *MacroModel Reference Manual* for a detailed description of this method.

Mixed torsional/Large scale low-mode sampling

This method uses a combination of the random changes in torsion angles and/or molecular position from the MCMM method, together with the low-mode steps from the LLMOD method used in Large scale low-mode. See the [MCMM](#), [LMC2](#), and [ARPK](#) opcode descriptions in the *MacroModel Reference Manual* for detailed descriptions of this method.

10.2.2 Global Search Restrictions

At the top of the CSearch folder, you can set restrictions on the scope of the search and its output. In addition, you can specify the criterion for equivalent structures in the lower portion of the panel.

Maximum number of steps

For each search method, specify the number of steps to be performed in the Maximum number of Steps text box. When the number of generated trial structures matches the value in field, the conformational search is terminated.

Distinguish enantiomers

By default, the search algorithm compares each trial structure and its enantiomer with the list of unique structures to determine if the trial structure is unique. If you select this option, enantiomers are considered to be separate structures.

Number of structures to save for each search

Use this text box to specify the number of unique structures to save when the search is complete, counting from the structure that is lowest in energy. If the value is zero, all unique structures are saved.

Energy window for saving structures

The default value is set to 50 kJ/mol. The value in this text box is used to compare trial structures. Any new structures generated and minimized are saved as results of the search only if they are within the energy window value above the current global minimum. Lowering this value results in the search saving fewer structures.

Max distance between atoms in equal structures

If comparison atoms are selected or Perform automatic setup during calculation is selected, structures produced during a conformational search are compared to see if they are unique. One of the conditions for structures to be considered different is that all corresponding atoms in the two structures lie within the value entered in this field. The default is 0.25 Å.

10.2.3 Low-Mode Parameters

In the lower portion of the Conformational Search panel are three text areas that are relevant only to low-mode searches. These are active only when any of the four methods involving Low-Mode Conformational Searching is the current search method.

Probability of a torsion rotation/molecule translation

This control is available only with the mixed methods. This text box allows you to set a probability that any defined torsion rotations and molecule translations will be made during the search. The value should be a number from 0.0 to 1.0.

Minimum distance for low-mode move

Maximum distance for low-mode move

Used for setting the minimum and maximum distance for a low-mode move. During a search, the fastest moving atom is moved at randomly generated distances that are between the minimum and maximum values specified in these text boxes.

10.2.4 Automatic Setup of Conformational Search Variables

Before you start a conformational search, the values for certain variables need to be set, such as which torsions are to be rotated and which rings are to be opened. There are three ways to set up conformational search variables from the CSearch folder: using Perform Automatic Setup, using Perform automatic setup during calculation, or setting variables manually. Setting conformational search variables manually is described in the next section.

The automatic setup procedures select comparison atoms and other conformational search parameters for you. For MCMM and mixed MCMM/LMOD calculations, comparison atoms, chiral atoms, torsion checks, molecule moves, variable torsions, and ring closures are automatically set up. For low-mode and large-scale low-mode calculations, comparison atoms, chiral atoms, and torsion checks are automatically set up.

Perform Automatic Setup

Use the Perform Automatic Setup button to automatically generate all the necessary variables for the search. The settings chosen by the automatic setup procedure can then be reviewed using the buttons in the center section of the panel. Clicking on one of these buttons opens the corresponding panel, which you can use to define and edit the conformational search variables. These panels are described in [Section 10.2.5 on page 104](#). Automatic setup may not be available for all conformational search methods. If it is not available, the Perform Automatic Setup button is dimmed. Perform Automatic Setup applies only to the structures in the Workspace. To set up conformational search variables for other structures, use Perform automatic setup during calculation.

If you perform an automatic setup on an entry that contains a substructure, only the atoms in the substructure are automatically set up with MCMM parameters and comparison atoms. Atoms in fixed and frozen regions are not included in automatic setup.

To clear any previously set conformational search variables before performing an automatic setup, click Reset All Variables.

Perform automatic setup during calculation

Perform automatic setup during calculation performs the same tasks as Perform Automatic Setup, but it permits you to set up parameters for all the structures that are processed when the search is performed. This option disables the controls used with Perform Automatic Setup. You can use this option with any of the conformational search methods, but you must use it with serial MCMM searches. This option inserts AUTO opcodes into the command file. See [Section 8.4 on page 83](#) and the *MacroModel Reference Manual* for more information on AUTO.

Note: Automatic setup is not intended for use with constraints. If you set constraints and use automatic setup, MacroModel can fail because the automatic setup can request movement of atoms that have constraints.

Torsion sampling options

You can use the options on this menu to determine which of the torsions around amide and ester linkages and other planar groups (azo and imido groups) are selected for sampling during Automatic Setup. The options are as follows:

- **Restricted**—Do not sample amide and ester derivatives or other planar groups. Torsional constraints are applied to ensure that the relative conformation of these groups is maintained during the MCMM and minimization process. (Sets AUTO arg8=0.)
- **Intermediate**—Sample the torsions of amide or ester linkages of non-standard groups like anhydrides, carbamates, hydrazones, and so on. Normal ester and amide linkages are not sampled, nor are azo and imido linkages. (Sets AUTO arg8=1.)
- **Enhanced**—Sample all amide-like and ester-like linkages, including standard amides and esters. Azo and imido linkages are not sampled. (Sets AUTO arg8=2.)
- **Extended**—Sample rotations around C=N and N=N bonds, in addition to all amide and ester derivatives. (Sets AUTO arg8=3.)

Use *N* steps per rotatable bond

Use this text box to specify the number of unique structures to save per rotatable bond. The total number of structures is still limited by the Maximum number of steps setting. This option is useful for serial searches, where the structures might have different numbers of rotatable bonds.

10.2.5 Setting Conformational Search Variables Manually

To set conformational search variables manually, or to examine them, use the buttons in the center of the CSearch folder. Each button opens a panel that can be used to set a particular kind of variable for a conformational search.

A description of each of these panels is given below. Each panel has some common features: a text box at the top that lists the defined features; a Show markers button, which controls display of the defined features; a Delete button, which deletes the selected definition from the list; and a Delete All button, which deletes all the defined features in the list. In any of the subpanels, you can edit a defined feature by selecting it in the list and picking Workspace atoms.

Markers can be displayed in the Workspace to identify the conformational search variables. Most of the panels have a Show markers option to control the display of markers for the vari-

ables set in that panel. When you close the panel, the markers are undisplayed. You can also display or undisplay markers for all variables by clicking Display All Markers or Undisplay All Markers.

To clear any previously set conformational search variables, click Reset All Variables.

10.2.5.1 Ring Closures

During a Monte Carlo-based conformational search (MCMM or SPMC), rings must be opened before varying their torsions. Then, after random torsion variations have been evaluated, the ring must be re-closed. The appropriate location for ring closures must be defined before a calculation can be started. “Ring Closures” consist of four atoms that define the bond to be broken (and then re-closed) during an MCMM or SPMC step.

The simplest way to create ring closures is to use the Perform Automatic Setup button. MacroModel makes reasonable decisions about where a ring should be opened. The list of atoms in the ring that will be opened appears in the text box at the top of the Ring Closures panel. Each ring closure is defined by four atoms. The opening is made between the second and third atoms. If Show markers is selected, the ring closure is marked by a light green line, which is broken between the second and third atoms. A lightning bolt is also placed between the second and third atoms to show where the opening will be made.

The functions of the Ring Closures panel are described below.




Figure 10.2. *The conformational search Ring Closures panel.*

Minimum distance

You must specify the minimum acceptable distance between the second and third ring closure atoms. If the distance between the ring closure atoms is less than this minimum distance, then the ring is reopened and a different set of random variations is performed. The default value for this is 0.5 Å. This value should suffice for most searches.

Maximum distance

The maximum acceptable distance between the second and third atoms of a re-closed ring can be specified. If the distance between the ring closure atoms is greater than this maximum distance, the ring is reopened and a different set of random variations is performed. The default value for this is 2.5 Å. This value should suffice for most searches.

Define ring closures

To define the location of a ring closure manually, choose **Atom** or **Bond** from the **Pick** menu and pick four atoms or three bonds from the ring. Picked atoms are marked by purple boxes. The bond is broken between the second and third picked atoms. Once four atoms have been picked, a new entry is displayed in the list and the ring closure is marked by the green line and lightning bolt, as described above.

10.2.5.2 Torsion Rotations

During a Monte Carlo conformational search, random changes are performed on the structure of the molecule. Structure energies are then evaluated. To set up a Monte Carlo-based search using either the MCMM or SPMC method, you must specify which torsions in the molecule can be rotated. These torsions are defined by pairs of atoms (defining bonds around which the structures can be rotated) using the **Torsion Rotations** panel.

The simplest way to define torsion rotations is to use the **Perform Automatic Setup** button. MacroModel detects torsions that need to be rotated and generates a list of torsion rotation definitions. This list appears in the text box at the top of the **Torsion Rotations** panel.

If **Show markers** is selected, the bonds defined as torsion rotations are indicated in the **Work-space** by a purple line and an icon (a purple horizontal line encircled by an arrow).

The functions of the **Torsion Rotations** panel are described below.




Figure 10.3. *The conformational search* Torsion Rotations *panel*.

Minimum rotation

The minimum acceptable rotation of a torsion must be specified. Each time that a torsion must be rotated, an increment for rotation larger than this value is selected. The default minimum rotation is 0°.

If a torsion rotation is added for a double or amide bond (i.e., to search both E- and Z- isomers around these bonds), the minimum value for rotations should be set to 90°. This maximizes the probability that interconversion between the two isomers will occur.

Maximum rotation

This text field allows you to specify a value for the maximum acceptable torsion rotation increment. Each time that a torsion must be rotated, an increment for rotation smaller than this value is selected. The default maximum rotation is 180°.

Define torsions to rotate

To define torsion rotations manually, choose Atom or Bond from the Pick menu, and pick two atoms or one bond in the Workspace structure to define the torsion rotation. When you pick the first of a pair of torsion rotation atoms, Maestro places a crimson box around it. Once both the atoms have been picked, a new entry is displayed in the list at the top of the panel, and the bond is marked as described above if Show markers is selected.

10.2.5.3 Molecule Trans/Rot

During a Monte Carlo conformational search, random changes are performed on the structure of the molecule. Structure energies are then evaluated to find the lowest-energy structure possible. To set up a Monte Carlo-based search using either the MCMM or SPMC method, you must specify which molecules will be translated and rotated. The Molecule Trans/Rot panel

facilitates the specification of molecules to be rotated and translated during a Monte Carlo (MCMM or SPMC) conformational search.

The simplest way to define molecules for translation and rotation is to use the **Perform Automatic Setup** button. MacroModel locates molecules that need to be translated and rotated, and generates a list of molecule translation and rotation definitions. This list appears in the text box at the top of the **Molecule Trans/Rot** panel.

Maestro labels the selected molecules with peach-colored vertical arrows encircled by “rotation” arrows if **Show markers** is selected. The currently selected molecule is distinguished by a turquoise label.

The controls of the **Molecule Trans/Rot** panel are described below.

Minimum rotation

The minimum acceptable rotation of a molecule must be specified. Each time a molecule must be rotated, an increment for rotation that is larger than this value is selected. The default minimum rotation is 0°.

Maximum rotation

The **Maximum rotation** text field allows you to specify a value for the maximum acceptable molecule rotation increment. Each time a molecule must be rotated, an increment for rotation that is smaller than this value is selected. The default maximum rotation is 180°.




Figure 10.4. The conformational search **Molecule Trans/Rot panel**.

Minimum translation

The value in the Minimum translation text box sets the lower limit for molecule translation. Each time a molecule must be translated, a random increment larger than this value is used. The default value is 0.0.

Maximum translation

This text box sets the upper limit for molecule translation. Each time a molecule must be translated, a random increment smaller than this value is used. The default value is 1.0.

Define molecules to translate/rotate

To define molecules to be translated or rotated manually, click on an atom in the Workspace structure. (Define molecules to translate/rotate is selected by default.) A new entry is displayed in the list at the top of the panel. The atom is colored peach and the vertical arrow icon is placed by the atom in the Workspace. Although only the picked atom is marked, the entire molecule is selected. You should not pick more than one atom in any molecule, even though this operation is allowed.

10.2.5.4 Comparison Atoms

During conformational searches, new structures are generated and minimized. The structures are compared against other low-energy structures that have already been found in the search. The comparison is performed by rigid superposition, comparing only those atoms specified as “comparison atoms” in the setup.

The simplest way to define comparison atoms is to click Perform Automatic Setup. MacroModel locates comparison atoms and generates a list of the atoms. This list appears in the text box at the top of the Comparison Atoms panel.




Figure 10.5. The Comparison Atoms panel.

The features of the Comparison Atoms panel work in the same way as for multiple minimization (see [page 81](#)).

To define comparison atoms, use one or more of the following options:

- Click Heavy Atoms + O-H,S-H.

This adds all non-hydrogen atoms and the hydrogen atoms attached to oxygen and sulfur to the list of comparison atoms.

- Click Heavy Atom.

This adds only the non-hydrogen atoms to the list of comparison atoms.

- Pick atoms in the Workspace.

Choose an object from the Pick menu, and click atoms in the Workspace to add the atoms to the list of comparison atoms.

- Select atoms using the Atom Selection dialog box.

For more complex combinations of comparison atoms, click Select to select atoms using the Atom Selection dialog box.

- Select all atoms.

Click All to add all atoms to the list of comparison atoms.

Comparison atoms are marked in light green in the Workspace with an “=” icon beside them. The currently selected comparison atom is marked in aquamarine.

10.2.5.5 Chiral Atoms

Because Monte Carlo conformational searches can generate and then minimize highly strained structures, chiral atoms in a molecule might be inverted. To prevent structures with inverted centers from appearing on the final list of optimized structures, you must identify a molecule's chiral atoms before beginning a search. Once defined, the chirality of each center is compared against that in the starting structure. If inversion has occurred, the search result is rejected.

The simplest way to define chiral atoms is to click Perform Automatic Setup. MacroModel locates chiral atoms and generates a list of the atoms. This list appears in the text box at the top of the Chiral Atoms panel.

If it is necessary to define chiral atoms manually, check that Define chiral atoms is selected (the default), and pick the atoms in the Workspace to add them to the list.

Chiral atoms are marked in peach with “RIS” labels beside them. The currently selected atom is colored turquoise. To hide these markers, deselect Show markers.




Figure 10.6. The conformational search Chiral Atoms panel.

10.2.5.6 Distance Check

Occasionally you might want to restrict the scope of a conformational search to generate only structures that are consistent with certain geometric constraints—for example, when experimentally obtained results such as NOE constraints are available. Distance checks are used to reject structures in which certain distances do not meet the specified criteria.

Distance checks cannot be set using the automatic setup process. They must be added by picking atoms from the structure in the Workspace.

To define a distance for checking, choose either Atom or Bond from the Pick menu, then pick two atoms or one bond in the Workspace. A new entry is added to the list at the top of the Distance Check panel. Maestro marks the defined distance check with a purple dotted line and




Figure 10.7. The conformational search Distance Check panel.

a “check” icon. The currently selected distance check is distinguished by two solid lines on either side of the dotted line.

For each pair of atoms that you select, you must define the desired distance and the allowed variation of that distance. The value entered in the **Desired distance** text box defines the target distance between two atoms. Structures with distances that differ from this value by more than the value in the **Allowed variation** text box are rejected.

10.2.5.7 Torsion Check

Because highly strained structures may be generated during Monte Carlo conformational searches, the geometry around double or amide bonds may be changed (from E- to Z- isomers, for example). You might want to restrict the scope of the conformational search to structures that retain the original torsions around these bonds. Torsion checks can be used to reject structures that do not retain the original geometry for these special cases.

The simplest way to define torsion checks is to click **Perform Automatic Setup**. MacroModel generates a list of the amides and double bonds in a structure. This list appears in the text box at the top of the Torsion Check panel.

Torsion checks can be added by picking atoms in the Workspace. To define a torsion for checking, choose either **Atom** or **Bond** from the **Pick** menu, then pick four atoms or three bonds in the Workspace. A new entry is added to the list at the top of the Torsion Check panel, and the torsion check is marked with a solid red line and a check mark. The currently selected torsion check is marked by solid lines on either side of the red line.




Figure 10.8. The Torsion Check panel.




Figure 10.9. *The conformational search* Ligand Bonds *panel.*

For each torsion check, the minimum and maximum torsional angle must be defined. Any search structures in which the checked torsional angle is not between the minimum and maximum angle is rejected and not included in search results.

10.2.5.8 Ligand Bonds

Conformational searches of inorganic complexes produce structures in which ligand positions vary with respect to the metal center. Maestro permits the specific definition of the ligand bonds around which the reorientation takes place.

To use this feature, choose Atom or Bond from the Pick menu and click on the desired bonds or atom pairs that define these bonds in the Workspace. A new entry appears in the list at the top of the panel, and Maestro marks the bond with a purple dotted line and a scissor icon. The selected bond is distinguished by purple lines on either side of the dotted line.

After the bonds have been selected, click Perform Automatic Setup in the Conformational Search panel. Maestro generates molecular rotation and translation commands for each fragment created by the defined ligand bonds. In addition, it produces torsion rotations, torsion checks, and chiral atom definitions. To view or edit the automatically generated settings, open the Torsion Rotations, Torsion Check, or Chiral Atoms panels.

10.3 Command File Examples

To run MacroModel calculations, a molecular structure file and a command file are required. The molecular structure file contains the structures to be used as input in the calculation. The command file contains the name of the input structure file, the name of the output structure file, and an ordered list of operation codes (*opcodes*) for the calculations.

Once you set up a job in Maestro and click either **Start** or **Write**, Maestro writes out a molecular structure file and a command file. For many types of jobs, command files written this way are complete and adequate, but for some types of jobs, you may need to adjust the Maestro-generated command file.

The command files and the log files for the examples given in this section can be found in `$SCHRODINGER/macromodel-vversion/samples/Examples`.

10.3.1 Conformational Search Using MCMM

An example command file appears below for a conformational search calculation that uses the MCMM search method. Descriptions of the opcodes in the file follow.

```
mcmn.mae
mcmn-out.mae
MMOD      0      1      0      0      0.0000      0.0000      0.0000      0.0000
FFLD      10      1      0      1      1.0000      0.0000      0.0000      0.0000
BDCC      0      0      0      0      41.5692 99999.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
MCMM      100     0      0      0      0.0000      0.0000      0.0000      0.0000
MCNV      2      4      0      0      0.0000      0.0000      0.0000      0.0000
MCSS      2      0      0      0      50.0000      0.0000      0.0000      0.0000
MCOP      1      0      0      0      0.0000      0.0000      0.0000      0.0000
DEMX      0      166     0      0      50.0000      100.0000      0.0000      0.0000
MSYM      0      0      0      0      0.0000      0.0000      0.0000      0.0000
COMP      1      2      3      4      0.0000      0.0000      0.0000      0.0000
COMP      5      6      7      9      0.0000      0.0000      0.0000      0.0000
COMP      10     14     15     16     0.0000      0.0000      0.0000      0.0000
COMP      17     25     26     27     0.0000      0.0000      0.0000      0.0000
COMP      28     29     30     31     0.0000      0.0000      0.0000      0.0000
TORS      1      14      0      0      0.0000      180.0000      0.0000      0.0000
TORS      6      25      0      0      0.0000      180.0000      0.0000      0.0000
TORS      14     15      0      0      0.0000      180.0000      0.0000      0.0000
TORS      15     16      0      0      0.0000      180.0000      0.0000      0.0000
TORS      25     26      0      0      0.0000      180.0000      0.0000      0.0000
CONV      2      0      0      0      0.0500      0.0000      0.0000      0.0000
MINI      9      0      500     0      0.0000      0.0000      0.0000      0.0000
```

MMOD: Creates and updates an intermediate structure file so that structures can be displayed in Maestro as the job progresses.

FFLD: Force field selection. Arg1 defines the force field used in the calculation (in this case MMFF94). Arg2 defines the electrostatic treatment for the calculation. In this case a constant dielectric is used. Arg4 is MMFF94 specific: arg4=1 defines the MMFF94s version of the force field, ensuring planarity around delocalized sp^2 nitrogens.

BDCO: Use the Bond Dipole CutOff (BDCO) method for truncating electrostatic interactions. Arg5 and arg6 are used to specify the cutoffs used for charge-dipole and charge-charge interactions, respectively.

READ: Read the input file.

MCMM: Use Monte Carlo Multiple Minimum searching. Arg1 defines the number of MC steps for the search.

MCNV: Sets the number of degrees of freedom to be varied at each MC step. If arg1 and arg2 differ, the search varies a random number of degrees of freedom between the numbers defined in arg1 and arg2. We recommend setting arg1=2 and arg2=maximum number of degrees of freedom.

MCSS: Select starting structure for the search steps. Arg1=2 defines use-directed selection of starting structures, where the least used structures will be used as starting geometries, as long as they are low enough in energy (as defined in arg5). This is more efficient in exploring new areas of the potential energy surface than, for instance, a random-walk starting geometry scheme. Arg5 gives the energy window for selecting a new starting structure. The new starting structure must be within arg5 kJ/mol of the lowest energy conformer found in the search.

MCOP: Monte Carlo options determine what and how often data is written to the log file. Arg1=1 ensures printing to the log file at every search step.

DEMXX: This command is used to prevent saving of high-energy conformers during the search. Arg5 defines the allowed energy window above the currently found global minimum. New conformers that are not within arg5 kJ/mol will be discarded. Additionally, a preliminary energy test can be performed during the energy minimization, to ensure that a reasonable structure has been found. Arg2 sets the number of energy iterations to be performed before the preliminary test (a good value is approximately 1/3 of the total number of energy iterations), while arg6 defines the energy above which conformers will be discarded (a value of about 1.5–2 times arg5 is recommended).

MSYM: Invokes the numbering symmetry library `mmsym`, which automatically and more generally identifies a suitable numbering order for use in comparing molecular conformations.

COMP: Arg1-arg4 list the atom numbers of atoms to be used in structural comparison with all previously found conformers. A maximum of 200 atoms can be used in the comparisons. For arg1=0, all heavy atoms are compared. By default, structures that have equivalent atoms separated by more than 0.25 Å upon superposition are considered different.

TORS: Defines the variable torsions in the molecule. Arg1 and arg2 are the atom numbers of the two central atoms defining a variable torsion. Arg5 and arg6 define the minimum and maximum dihedral angle variations (in both directions).

CONV: Defines convergence criteria. Arg1=2 specifies derivative convergence (default criterion is 0.05 kJ/mol-Å, and this value is set in arg5).

MINI: Starts the minimization. Arg1 defines the type of minimization algorithm to be used. Arg1=9 means that Truncated Newton-Raphson Conjugate Gradient will be used. In arg3, the number of minimization steps is defined. Arg3 can be set to a large number since the calculation automatically stops as soon as the convergence criterion has been reached.

10.3.2 Multi-structure Conformational Search Using LMOD

An example command file appears below for a conformational search calculation using low-mode conformational searching. Descriptions of the opcodes in the file follow. However, opcodes that are also in the MCMC command file example are not repeated. See the explanations in [Section 10.3.1](#).

```
serial-lmcs.mae
serial-lmcs-out.mae
FFLD      10      1      0      1      1.0000      0.0000      0.0000      0.0000
BDCC      0      0      0      0      41.5692 99999.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
LMCS      100     0      0      0      0.0000      0.0000      3.0000      6.0000
MCSS      2      0      0      0      50.0000      0.0000      0.0000      0.0000
MCOP      1      0      0      1      0.0000      0.0000      0.0000      0.0000
DEMX      0      166     0      0      50.0000      100.0000      0.0000      0.0000
MSYM      0      0      0      0      0.0000      0.0000      0.0000      0.0000
COMP      0      0      0      0      0.0000      0.0000      0.0000      0.0000
CONV      2      0      0      0      0.0500      0.0000      0.0000      0.0000
MINI      1      0      500     0      0.0000      0.0000      0.0000      0.0000
```

LMCS: Use the Low-Mode Conformational Search method. Arg1=100 means that 100 Monte Carlo steps will be carried out before the calculation stops.

MCSS: Select starting geometries for Monte-Carlo search steps. Arg1=2 tells MacroModel to use as starting geometries structures whose energies are allowed by arg2 and arg5, and are used the fewest times as starting structures.

MCOP: Arg1=1 prints search results every step. Arg4=1 signifies that this is a low-mode serial search. A separate conformational search will be performed on all molecules in the input file. Arg5=0 implies a pure low-mode conformational search, rather than a mixed-mode, as in a previous example.

COMP: Setting arg1=0 allows for all heavy atoms to be compared for each structure in the input file. This removes the need to use atom numbers for the individual structures.

10.3.3 Mixed MCMM/Low-Mode Search Using a Substructure File

This example uses the CDK2 structure 1e1v with the co-crystallized ligand cmg. Alternate binding modes of the ligand are sampled using a mixture of torsional moves, low-mode moves, and rotation and translation of the ligand, in a conformational search. This example was prepared with Maestro. A substructure is used, in which the freely moving region includes the ligand and any residues with atoms within 3.0 Å of the ligand. A fixed and frozen region was also set up in the Substructure folder. The conformational search parameters were initially set by using the Perform Automatic Setup button in the CSearch folder. The parameters were then modified using the individual parameter panels in the CSearch folder to arrive at the sample instruction file below. Only torsions in the ligand are being explicitly sampled in this example.

Descriptions of the opcodes in the file follow. However, opcodes in the previous command files in the chapter are not repeated. See the previous examples for those explanations.

```
subs-mcmm-lmod.mae
subs-mcmm-lmod-out.mae
MMOD      0      1      0      0      0.0000      0.0000      0.0000      0.0000
FFLD      11      1      0      0      1.0000      0.0000      0.0000      0.0000
EXNB      0      0      0      0      0.0000      0.0000      0.0000      0.0000
BDCC      0      0      0      0      89.4427 99999.0000      0.0000      0.0000
SUBS      0      0      0      0      0.0000      0.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
CRMS      0      0      0      0      10.0000      0.5000      0.0000      0.0000
LMCS      10      0      0      0      0.0000      0.0000      3.0000      6.0000
MCNV      2      4      0      0      0.0000      0.0000      0.0000      0.0000
MCSS      2      0      0      0      50.0000      0.0000      0.0000      0.0000
MCOP      1      0      0      0      0.5000      0.0000      0.0000      0.0000
DEMX      0      666      0      0      50.0000      100.0000      0.0000      0.0000
COMP      117     121     128     133      0.0000      0.0000      0.0000      0.0000
COMP      136     140     158     162      0.0000      0.0000      0.0000      0.0000
COMP      165     169     174     176      0.0000      0.0000      0.0000      0.0000
COMP      223     225     231     236      0.0000      0.0000      0.0000      0.0000
COMP      439     445     466     470      0.0000      0.0000      0.0000      0.0000
COMP      473     476     482     487      0.0000      0.0000      0.0000      0.0000
COMP      1291    1296    1299    1303      0.0000      0.0000      0.0000      0.0000
COMP      1307    1309    1313    1318      0.0000      0.0000      0.0000      0.0000
COMP      1326    1362    1366    1370      0.0000      0.0000      0.0000      0.0000
COMP      1372    2105    2109    2115      0.0000      0.0000      0.0000      0.0000
COMP      2120    2155    2159    2164      0.0000      0.0000      0.0000      0.0000
COMP      2172    2323    2327    2331      0.0000      0.0000      0.0000      0.0000
COMP      2333    4599    4600    4601      0.0000      0.0000      0.0000      0.0000
COMP      4602    4603    4604    4605      0.0000      0.0000      0.0000      0.0000
COMP      4606    4607    4608    4609      0.0000      0.0000      0.0000      0.0000
COMP      4610    4611    4612    4613      0.0000      0.0000      0.0000      0.0000
COMP      4614    4615    4616      0      0.0000      0.0000      0.0000      0.0000
MSYM      0      0      0      0      0.0000      0.0000      0.0000      0.0000
```

CHIG	117	119	165	167	0.0000	0.0000	0.0000	0.0000
CHIG	223	439	468	1291	0.0000	0.0000	0.0000	0.0000
CHIG	1311	1364	2107	2157	0.0000	0.0000	0.0000	0.0000
CHIG	2325	0	0	0	0.0000	0.0000	0.0000	0.0000
TORS	4600	4610	0	0	0.0000	180.0000	0.0000	0.0000
TORS	4604	4605	0	0	0.0000	180.0000	0.0000	0.0000
TORS	4605	4606	0	0	0.0000	180.0000	0.0000	0.0000
TORS	4606	4611	0	0	0.0000	180.0000	0.0000	0.0000
MOLS	4633	0	0	0	0.0000	180.0000	0.0000	1.0000
TORC	135	134	132	133	90.0000	180.0000	0.0000	0.0000
TORC	164	163	161	162	90.0000	180.0000	0.0000	0.0000
TORC	1310	1309	1307	1308	90.0000	180.0000	0.0000	0.0000
TORC	2118	2117	2115	2116	90.0000	180.0000	0.0000	0.0000
CONV	2	0	0	0	0.1000	0.0000	0.0000	0.0000
MINI	9	0	1000	0	0.0000	0.0000	0.0000	0.0000

SUBS: With `arg1=0` look for an `.sbc` file that contains information on how the substructure is set up in the system.

MCOP: `Arg1=1` print search results every step. `Arg5=0.5` specifies that this mixed Monte Carlo/Low-Mode search will attempt Monte Carlo torsional and/or molecular translation moves half the time. The rest of the attempted moves will be generated using low-mode moves.

MOLS: `Arg1-4` specifies one atom from each molecule to be moved using molecular translation/rotation moves. `Arg5` and `arg6` specify the minimum and maximum rotation angles while `arg7` and `arg8` specify the minimum and maximum translation distances.

10.4 Checking and Interpreting Results

Like other methods that involve minimizations, it is important to check that all of the minimizations for the conformers generated are converged. If some of the results are not converged, they can be minimized further using multiple minimization. Often conformational searches yield large numbers of conformers. Clustering the conformers into families of similar structures can lead to useful insights. XCluster is an excellent tool for performing clustering analysis on collections of conformers obtained from conformational searches. See [Chapter 18](#) and the *MacroModel XCluster Manual* for more information.

Ligand Torsional Searches

When generating collections of conformations for ligands, it is often sufficient to generate broad rather than exhaustive coverage of conformational space. With growing interest in generating conformations for many ligands, speed is becoming an important consideration. For this purpose, MacroModel uses the ConfGen conformational generation utility, originally developed for rapidly and systematically exploring ligand conformations in Glide.

ConfGen carefully examines the structure of the ligand to understand where to expect local minima as a function of rotations about rotatable bonds. It then systematically generates the conformations that arise from various combinations of these local minima. Thus it provides a broad and fairly uniform coverage of the available conformational space. This systematic approach used in ConfGen avoids the enormous amount of resampling of conformations that occurs in most conformational searching methods designed for exhaustive sampling. In addition to sampling rotatable bonds ConfGen also samples ring conformations, chiral nitrogen atom inversions and, amide bond conformations. For a more detailed description of the conformational search, see [Section 3.2](#) of the *Glide User Manual*.

ConfGen's speed and ability to generate compact collections of quality candidate conformations forms a powerful combination with MacroModel's force fields, GB/SA solvation model, minimization procedures and redundant conformer elimination facilities to provide very good coverage of conformational space. This facility is focused on the efficient generation of ligand conformations and may not be applicable or suitable to other types of problems.

ConfGen functions differently from the conformational searching methods supported in MacroModel: it is specifically focused on extremely efficient generation of ligand conformations similar to those found in protein-ligand complexes. As a result, it is implemented somewhat differently and must be described separately from MacroModel's conformational searching procedures (see [Chapter 10](#)).

This chapter describes how to use ConfGen from the Ligand Torsional Search panel in Maestro, and provides an example command file with comments. A detailed description of ConfGen is provided under the [CGEN](#), [CGOP](#), and [CHYD](#) opcodes in the *MacroModel Reference Manual*.

ConfGen is a separate module and requires additional licensing. See [Obtaining a License](#) in the *Installation Guide* for information on obtaining licenses.

11.1 The Ligand Torsion Search Panel

To conduct ConfGen conformational searches from Maestro, choose Ligand Torsion Search from the MacroModel submenu of the Applications menu. This panel contains three folders: Potential, Mini, and LTSearch, which are described in detail in the following sections.

11.1.1 The Potential Folder

This folder is a modified version of the standard Potential folder described in [Section 5.2 on page 54](#). The differences are as follows:

- The force field list is limited to the MMFF, MMFFs, and OPLS_200X force fields.
- There is no text box for the H-bond cutoff distance since the supported force fields do not have a distinct representation for hydrogen bonds.
- There is an option labeled Suppress hydrogen bonding electrostatics.

Select this option to suppress electrostatic interactions between the ligand donor and acceptor groups. This option can aid in generating conformations closer to those found in protein-ligand complexes, where ligand donor and acceptor groups are likely to be hydrogen-bonded to the protein rather than to each other.




Figure 11.1. The Potential *folder of the* Ligand Torsion Search *panel*.

Deselect this option to model the interactions between donor and acceptor groups according to the force field and solvation model.

For more information, see the description of the [CHYD](#) opcode in the *MacroModel Reference Manual*.

It is strongly recommended you use either a distance-dependent dielectric or GB/SA solvation for these calculations.

11.1.2 The Mini Folder

This folder is very similar to the standard Mini folder described in [Section 7.2 on page 73](#). The only difference is that in this panel there are two settings for the maximum number of iterations instead of only one.

Pre-minimization of input structures

In this text box, you specify the maximum number of iterations for the input structures.

In most procedures, the same maximum number of iterations are used in minimizing the structures read in from the input structure file and those generated in subsequent processing. However, since ConfGen calculations are sensitive to the bond lengths and angles in the seed structures, it is a good idea to thoroughly minimize the structures from the input structure file.

Post-minimization of generated structures

In this text box, you specify the maximum number of minimization steps for the ConfGen generated structures.

Unlike other search methods employed by MacroModel, ConfGen attempts to generate conformations that are close to a local potential minimum. As a result, for most applications, it is not necessary to carefully minimize these generated structures.

11.1.3 The LTSearch Folder

This folder controls how the ConfGen search is conducted in terms of the number and types of conformations generated, and the elimination of high-energy or redundant conformers.

Number of search moves

Specify the target number of ConfGen generated conformations to be processed by MacroModel. If ConfGen finds more conformers than this for a given ligand, it uniformly selects conformations from this collection. If ConfGen finds fewer conformers than this for a given ligand, only that number of conformers are processed.




Figure 11.2. The LTSearch folder of the Ligand Torsion Search panel.

Save at most *N* conformations per ligand

Normally, all acceptable structures produced by MacroModel are written to the output structure file. With this option, you can enter an upper bound on the number of conformers saved for each ligand.

Distinguish enantiomers

Search mode

The search mode controls what combinations of rotations about terminal rotatable bonds (TRB) are sampled. In Rapid mode, the minima for each TRB is sampled separately with all other TRB's in their lowest energy orientation. In Thorough mode, all combinations of TRB minima are sampled.

Amide bonds

Select one of these options to control how amide bonds that do not lie in a ring are sampled by ConfGen:

- Vary amide bond conformation
- Retain original amide bond orientation
- Set amide bond conformation to trans

Sample rings

Select this option to enable ConfGen to seek specific templates for rings systems present in each ligand and to activate Maximum ring conformations. Deselect this option to prevent ConfGen from sampling ring systems.

Maximum ring conformations

Enter the maximum number of combinations of ring conformations to be sampled per ligand. This option is activated by the Sample rings option above.

Energy window for saving conformations

The lower limit of the energy window is the MacroModel force field energy of the lowest energy conformer for that ligand. The upper limit of the energy window is the lower limit plus the amount in this text box. Any conformer with a MacroModel force field energy greater than the upper limit is not retained in the output structure file.

Compare conformers by

Redundant conformers may be eliminated using geometric criteria. Redundant conformer elimination involves comparing each pair of conformers to see if they are different enough to justify retaining both of them.

- None: no redundant conformers are eliminated
- Heavy atoms: conformer comparisons are preformed using the positions of non-hydrogen atoms (activates Max distance between atoms in equal structures).
- Heavy atoms plus polar hydrogens: conformer comparisons are performed using the positions of non-hydrogen atoms and the torsional angles for polar hydrogens (activates both Max distance between atoms in equal structures and the Max torsional angle difference for polar H's).

Max distance between atoms in equal structures

Conformers are considered distinct if at least one atom in the superimposed conformers is different in position by a distance greater than the value in this text box. The default is 0.5 Å. (This text box is only activated if Compare conformers by is set to either Heavy atoms or Heavy atoms plus polar hydrogens).

Max torsional angle difference for polar H's

Conformers are considered distinct if a dihedral angle involving a polar hydrogen in the two conformers differs by more than the value in this text box. The default is 60°. (This text box is only activated if Compare conformers by is set to Heavy atoms plus polar hydrogens).

11.2 Setting up ConfGen Calculations

To set up a ConfGen calculation:

1. Choose Ligand Torsion Search from the MacroModel submenu of the Applications menu.
2. Choose Workspace (included entries) or Project Table (selected entries) from the Use structures from option menu.
3. Select the appropriate settings in the Potential, Mini, and LTSearch folders.

It is strongly recommended that you use either a distance-dependent dielectric or GB/SA solvation.

4. Click Start to set up and launch the job (see [Section 5.6 on page 64](#)), or click Write to write the jobs files so they can be run manually from the command line or edited for additional customization.

11.3 Command File Example

To run MacroModel calculations, a molecular structure file and a command file are required. The molecular structure file contains the structures to be used as input in the calculation. The command file contains the name of the input structure file, the name of the output structure file, and an ordered list of operation codes (*opcodes*) for the calculations.

Once you set up a job in Maestro and click either Start or Write, Maestro writes out a molecular structure file and a command file. For many types of jobs, command files written this way are complete and adequate, but for some types of jobs, you may need to adjust the Maestro-generated command file.

The command files and the log files for the examples given in this section can be found in \$SCHRODINGER/macromodel-*vversion*/samples/Examples.

Below is an example of a .com file for a ConfGen calculation. Descriptions of the opcodes in the file follow. More information on these opcodes can be found in [Chapter 3](#) of the *MacroModel Reference Manual*.

```
mmod_confgen.mae
mmod_confgen-out.mae
MMOD      0      1      0      0      0.0000      0.0000      0.0000      0.0000
FFLD      10      1      0      1      1.0000      0.0000      0.0000      0.0000
SOLV       3      1      0      0      0.0000      0.0000      0.0000      0.0000
BDCC       0      0      0      0      41.5692 99999.0000      0.0000      0.0000
READ       0      0      0      0      0.0000      0.0000      0.0000      0.0000
CRMS       0      0      0      0      0.0000      0.5000      60.0000      0.0000
CGEN      1000      0      0      2      0.0000      0.0000      0.0000      0.0000
```

CGOP	0	2	50	1	48.0000	8.0000	2.0000	0.0000
CHYD	1	0	0	0	0.0000	0.0000	0.0000	0.0000
MCOP	1	0	0	0	0.0000	0.0000	0.0000	0.0000
DEMX	0	0	0	0	25.0000	0.0000	0.0000	0.0000
MSYM	0	0	0	0	0.0000	0.0000	0.0000	0.0000
NANT	0	0	0	0	0.0000	0.0000	0.0000	0.0000
AUTO	0	4	1	2	-1.0000	1.0000	0.0000	0.0000
CONV	2	0	0	0	0.0500	0.0000	0.0000	0.0000
MINI	9	0	100	0	0.0000	0.0000	0.0000	0.0000

MMOD: Creates and updates an intermediate structure file so that structures can be displayed in Maestro as the job progresses.

FFLD: Specifies the force field and aspects of the electrostatic treatment to use in MacroModel. Arg1 defines the force field used in the calculations. 10 means MMFF or MMFFs will be used. ConfGen calculations are compatible with the MMFF, MMFFs, OPLS_2001 and OPLS_2005 force fields. Arg2 defines the electrostatic treatment for the calculation. In this case a constant dielectric is used due to the use of solvation model 3 (see SOLV below). Arg4 selects either MMFF (0) and MMFFs (1). Arg5 is used to indicate that the dielectric constant should be set to 1 and is appropriate for solvation model 3.

SOLV: Controls which solvation model to use in MacroModel. Arg1 specifies the type of solvation treatment to use. 1 means that GB/SA is used. Arg2 selects the solvent and a value of 3 corresponds to water.

BDCO: Use the Bond Dipole CutOff (BDCO) method for truncating electrostatic interactions. Arg5 and arg6 are used to specify the cutoffs used for charge-dipole and charge-charge interactions, respectively.

READ: Read the input file.

CRMS: Sets parameters for redundant conformer elimination. Arg6 controls the maximum distance between atoms in equal conformers and is set to 0.5 Å. Arg7 controls the maximum difference in dihedral angles for polar hydrogens and is set to 60°.

CGEN: Use ConfGen to generate conformers. Arg1 defines the maximum number of conformers MacroModel should process during the search. Arg4 controls what combinations of rotations of terminal rotatable bonds are sampled. A value of 2 means that all combinations of minima for terminal rotatable bonds are sampled.

CGOP: Selects options for ConfGen. Arg1 is used to control sampling of certain types of compact symmetric terminal groups e.g. $-\text{CF}_3$ or $-\text{SO}_3^-$. A value of 0 means do not conduct rotations of these groups. Arg2 controls how rings are sampled. Arg2=2 means use specific ring system templates when sampling ring flexibility. Arg3 controls minimization of conformers produced by ConfGen. If it is a negative value, the structures are not minimized. If

it is a positive value, the structures are minimized using this number of steps. Arg3=0 means the number of steps listed in MINI Arg3 are used. Arg4 controls the sampling of non-ring amide bonds. Arg4=1 specifies that both cis and trans conformations should be generated. Arg5 controls the maximum relative energy between conformers inside ConfGen. Conformers with estimated energies higher than 50 kJ/mol above the lowest energy conformer are eliminated. Arg6 controls the maximum number of ring conformations to sample for the ligand as a whole. Arg7 controls the maximum number of ring conformations to sample for each ring system.

CHYD: Suppress hydrogen bond electrostatic interactions.

MCOP: Options to control what and how often data is written to the .log file during a conformation generation calculation. Arg1=1 updates the log file for every search step.

DEMX: This command is used to prevent saving of high-energy conformers during the search. Arg5 defines the allowed energy window above the currently found global minimum. New conformers that are not within arg5 kJ/mol will be discarded. Additionally, a preliminary energy test can be performed during the energy minimization, to ensure that a reasonable structure has been found.

MSYM: Invokes the numbering symmetry library `mmsym`, which automatically and more generally identifies a suitable numbering order for use in comparing molecular conformations.

NANT: Regards enantiomers as distinct and retains both forms in the output structure file.

AUTO: Run MacroModel's automatic setup for conformational generation. Arg2 controls what atoms and types of comparisons amongst the atoms are used to identify redundant conformers. A value of 4 indicates that the positions of non-hydrogen atoms and dihedral angles for polar hydrogens should be used. Arg3 controls chiral atom checking during the calculation. Arg3=1 means identify and enforce chirality conservation for chiral atoms during the search. Arg4 controls the identification of torsional constraints during conformational searches. Arg4=2 means that checks are carried out to ensure that C=C bond geometries (i.e. cis or trans) match those in the input structure file. Arg5 controls the identification of torsions to sample within MacroModel and a value of -1 turns this off since this is handled within ConfGen itself. Arg6 is used to indicate whether this is a serial calculation (i.e. perform a separate conformational search on each structure in the input structure file) or not. Arg6=1 means that this is a serial calculation. All ConfGen/MMOD jobs should be run as serial calculations.

CONV: Defines convergence criteria. Arg1=2 specifies derivative convergence (default criterion is 0.05 kJ/mol-Å, and this value is set in arg5).

MINI: Starts the minimization. Arg1 defines the type of minimization algorithm to be used. Arg1=9 means that Truncated Newton-Raphson Conjugate Gradient will be used. Arg3 defines the number of minimization steps used in the structures read in from the input structure file

prior to passing them to ConfGen. Arg3 also defines the number of minimization steps used on structures generated by ConfGen if CGOP arg3=0.

Dynamics Calculations

Molecular dynamics simulations in MacroModel use classical mechanics (Newton's equations of motion) to mimic how the system would behave as a function of time, typically at or close to the temperature of interest. You can use molecular dynamics studies to learn about the thermal variations within a system or to permit it to relax out of a local minimum structure into related but more probable structures.

12.1 The Dynamics Panel

Using the Dynamics panel, you can set up dynamics calculations and either submit the calculation or write the job files for later use. The Dynamics panel consists of the upper, general portion that is common to other MacroModel panels. This section is discussed in [Section 5.1 on page 53](#). The panel contains six tabbed folders: Potential, Constraints, Substructure, Mini, Monitor, and Dynamics. The Monitor and Dynamics folders are unique to the Dynamics panel. The controls in these folders are explained below. For information on the other folders, see [Section 5.2 on page 54](#) through [Section 5.5 on page 61](#).

To open the Dynamics panel, select Dynamics from the MacroModel submenu of the Applications menu on the main menu bar.

12.2 Performing a Dynamics Calculation

To set up parameters for a dynamics calculation, first set values in the upper part of the Dynamics panel and in the Potential, Constraints, Substructure, and Mini folders, then set values in the Monitor and Dynamics folders, which are described below.

12.2.1 The Monitor Folder

During a molecular dynamics simulation, you can monitor a number of geometrical parameters. These are angles, inter-atomic distances, dihedral angles, the surface areas of individual atoms, and the population of hydrogen-bonds. The Monitor folder of the Dynamics panel contains a series of buttons that open panels that allow you to define these parameters. The results of the monitoring appear in the *jobname.mmo* file. Information is also recorded in the *jobname.log* and the *jobname-out.mae* files. The latter is used to convey the property values to Maestro's project table.




Figure 12.1. The Monitor *folder of the* Dynamics *panel*.

You can also sample structures during a simulation. To do this, simply enter the total number of structures you want to sample in the Number of structures to sample: text box. This number of structures is selected at regular intervals throughout the simulation and written to the output (*jobname-out.mae*) file.

Five of the buttons in the Monitor folder open panels that allow you to choose what to monitor. Each of these panels is described in a section below. Each panel has a text area in which the atoms that define each geometrical parameter are listed, an atom selection tool, and Delete and Delete All buttons to delete one or all items from the list. You can redefine an item by selecting it and then picking atoms in the workspace. When parameters are defined, the atoms defining them are marked in the Workspace if Show Markers is selected. When you close the panel, the markers are automatically cleared.

The remaining button, Reset All, removes all entries from all dynamics monitors.




Figure 12.2. The Monitored Atoms Surface Area panel.

12.2.1.1 Surface Areas

Use the Monitored Atom Surface Areas panel to specify the atoms whose surface areas are to be monitored during the simulation. To monitor the surface area of specific atoms, use the selection tool to pick atoms. The selection tool allows you to specify atoms in the following ways:

- Choose a structural unit from a Pick menu and pick atoms in the Workspace belonging to structural units of the selected type.
- Choose all atoms by clicking the All button.
- Choose atoms using the Atom Selection dialog box, which creates ASL expressions that define the selected atoms. Open the Atom Selection dialog box by clicking the Select button.

After each atom is picked, a new entry appears in the list of atoms in the Monitored Atoms List located at the top of the panel. When an atom is picked, Maestro marks it and places an eye icon beside it. The currently selected atom is colored teal, the other selected atoms are colored green.

12.2.1.2 Distances

Use the Monitored Distances panel to specify interatomic distances for monitoring. To define a distance to be monitored, choose Atom or Bond from the Pick menu, and click on two atoms or one bond in the Workspace. When the distance is defined, the two atoms appear in the list at the top of the panel. The defined distances are marked with a purple dotted line and eye icon. The currently selected distance is distinguished by a solid line on either side of the dotted line.




Figure 12.3. The Monitored Distances *panel*.

12.2.1.3 Angles

Use the Monitored Angles panel to specify angles for monitoring. To define an angle to be monitored, choose Atom or Bond from the Pick menu, and click on three atoms or two bonds in the Workspace that define an angle. When the angle is defined, the three atoms appear in the list at the top of the panel. The defined angles are marked with a green solid line, a dotted line through the angle, and an eye icon. The currently selected angle is distinguished by a solid line on either side of the dotted line.

12.2.1.4 Torsions

Use the Monitored Torsions panel to select dihedral angles for monitoring. To define a dihedral angle to be monitored, choose Atom or Bond from the Pick menu, and click on four atoms or three bonds in the Workspace that define a dihedral angle. When the dihedral angle is defined, the four atoms appear in the list at the top of the panel. The defined dihedral angles are marked




Figure 12.4. The Monitored Angles *panel*.




Figure 12.5. The Monitored Torsions *panel*.

with a red solid line, a dotted line through the dihedral angle, and an eye icon. The currently selected angle has a solid line on either side of the dotted line.

12.2.1.5 H-Bonds

During a dynamics simulation, MacroModel periodically examines the geometry around each monitored H-bond. If the bond meets the three user-specified H-bond criteria, it is counted in the H-bond population survey.

To select H-bonds from the structure in the Workspace, click the four atoms that define the bond and the associated angles and distances. The picking order of the atoms is important. Start with the heavy atom of the donor pair (designated X), followed by the hydrogen atom




Figure 12.6. The Monitored H-Bonds *panel*.

(designated H), the acceptor (designated Y), and an atom attached to the acceptor (designated Z) to define an angle. For example, in a typical N-H to O=C hydrogen bond, the picking order should be N, H, O, then C. After the fourth atom has been selected, a new entry appears in list at the top of the panel. The H-bond is marked with solid yellow lines between X and H and between Y and Z, and a dotted yellow line between H and Y. The currently selected H-bond has solid lines on either side of the dotted line.

Once you select an entry, you can specify the monitoring criteria for that entry by entering a maximum H...Y distance, a minimum X-H...Y angle, and a minimum H...Y-Z angle in the appropriate text boxes. If you do not specify values, the default values are used. The default value for the H...Y distance is 2.5 Å, the minimum X-H...Y angle has a default value of 120 degrees, and the minimum H...Y-Z angle is 90 degrees. If the distance is greater than the maximum or either of the angles is smaller than the minimum, the H-bond is not counted in the current population survey.

12.2.2 The Dynamics Folder

The Dynamics folder of the Dynamics panel contains settings for defining the dynamics calculation method and the calculation conditions.

The controls in the Dynamics folder are described below.

Method

This menu contains two options:

- **Stochastic dynamics (default):** The most common method because it includes random forces that simulate the buffeting of a system by solvent molecules. In combination with damping forces, the random forces provide a simple, robust method for controlling temperature in a simulation. The random forces can also assist in sampling the potential energy surface. See Reference [27](#) for more information.
- **Molecular dynamics:** Uses a standard constant temperature velocity-Verlet algorithm.

SHAKE

The SHAKE procedure [[28](#)] constrains selected bond lengths to their original values. This procedure allows the use of larger time steps than unconstrained simulations. For most simulations, the Bonds to Hydrogens choice from the SHAKE option menu should be sufficient. This choice allows time steps up to 2 fs.

Simulation temperature (K)

The temperature at which a simulation is carried out. The default value is 300.0 K.

Time step (fs)

Determines the time step used in the integration of the equations of motion during the simulation. Smaller values lead to more accurate but computationally intensive simulations for a given amount of time simulated. The default value is 1.5 fs.

Equilibration time (ps)

Determines the length of the “settling down” period at the start of the simulation. The equilibration period is needed to allow initial velocities to stabilize before any monitoring data or sampled structures are written to output files.

Simulation time (ps)

The total time allowed for a given simulation. The default is 10 ps. However, if you want converged results, specify a much longer time, even for small to medium-sized systems.




Figure 12.7. The Dynamics folder of the Dynamics panel.

12.3 Command File Examples

To run MacroModel calculations, a molecular structure file and a command file are required. The molecular structure file contains the structures to be used as input in the calculation. The command file contains the name of the input structure file, the name of the output structure file, and an ordered list of operation codes (*opcodes*) for the calculations.

Once you set up a job in Maestro and click either Start or Write, Maestro writes out a molecular structure file and a command file. For many types of jobs, command files written this way are complete and adequate, but for some types of jobs you may need to adjust the Maestro-generated command file.

The command files and the log files for the examples given in this section can be found in `$SCHRODINGER/macromodel-vversion/samples/Examples`.

12.3.1 Stochastic Dynamics

Below is an example of the command file for a stochastic dynamics simulation and explanations of the opcodes that appear in the file.

```
sdyn.mae
sdyn-out.mae
MMOD      0      1      0      0      0.0000      0.0000      0.0000      0.0000
FFLD      10      1      0      1      1.0000      0.0000      0.0000      0.0000
BDCCO      0      0      0      0      41.5692 99999.0000      0.0000      0.0000
READ       0      0      0      0      0.0000      0.0000      0.0000      0.0000
CONV       2      0      0      0      0.0500      0.0000      0.0000      0.0000
MINI       9      0      500    0      0.0000      0.0000      0.0000      0.0000
MDIT       0      0      0      0      300.0000      0.0000      0.0000      0.0000
MDYN       0      1      1      0      1.5000      1.0000      300.0000      0.0000
MDSA      10      0      0      0      0.0000      0.0000      1.0000      0.0000
MDDA      14      15      16      17      0.0000      0.0000      0.0000      0.0000
MDYN       1      1      1      0      1.5000      10.0000      300.0000      0.0000
WRIT       0      0      0      0      0.0000      0.0000      0.0000      0.0000
```

MMOD: Creates and updates an intermediate structure file so that structures can be displayed in Maestro as the job progresses.

FFLD: Force field selection. Arg1 denotes the actual force field used in the calculation (in this case MMFF94). Arg2 defines the electrostatic treatment for the calculation. The default (arg=0) is to use the dielectric treatment encoded in the force field, however, in this case a constant dielectric is used. Arg4 is MMFF94-specific. Arg4=1 defines the MMFF94s version of the force field, ensuring planarity around delocalized sp^2 nitrogens.

BDCO: Use the Bond Dipole CutOff (BDCO) method for truncating electrostatic interactions. Arg5 and arg6 are used to specify the cutoffs used for charge-dipole and charge-charge interactions, respectively.

READ: Read the input file.

CONV: Defines convergence criteria. Arg1=2 signifies derivative convergence (default, if no CONV command is present, criterion is 0.05 kJ/mol-Å; this value is set in arg5).

MINI: An energy minimization precedes the dynamics calculation in order to eliminate excess potential energy. Arg1=9 indicates that the TNCG minimization method should be used for arg3 minimization iterations.

MDIT: Apply random initial velocities corresponding to 300K (arg5) to all atoms.

MDYN: Perform the dynamics simulation. Arg2=1 selects the use of the SHAKE protocol to constrain hydrogen bonds to their natural values. Arg3=1 sets up a stochastic dynamics run. The time step (in fs) is set in arg5, the length of the simulation (in ps) in arg6, and the temperature (in K) of the simulation in arg7.

Two MDYN command lines appear in this command file. The first performs a short equilibration run, while the second is where the actual sampling occurs. Arg1 defines printing of energy listing. In the first MDYN line, arg1=0, meaning that a summary of the energies is printed to the log file only. In the second MDYN line, arg1=1, signifying that energies as well as monitoring results are written to the .mmo file.

MDSA: Perform structure sampling during the stochastic dynamics simulation. Intermediate structures generated in the dynamics simulation are saved to the output structure file at regular intervals during the calculation. Arg1 defines the total number of structures to sample (number of “snapshots” to take) during the simulation. Structures can also be written out based on time-intervals (in ps) using arg5. Arg7=1 forces deletion of all structures stored in the -out.mae file prior to each simulation.

MDDA: Monitor a dihedral angle during the simulation. Arg1-arg4 defines the atom numbers of the angle to be monitored. The default arg5 divides the reported results into 10-degree increments. Average values are reported in the .log file while more detailed information (including the distribution of angles) is reported in the .mmo file.

WRIT: Write the final structure to the output file.

12.3.2 Simulated Annealing

Simulated annealing is often used to relax structures into a lower energy state a dynamics technique by systematically lowering the temperature used in the simulations. Below is an example

command file for a simulated annealing calculation, where the target temperature for the system is changed in a step-wise manner. The temperature can also be changed in a continuous manner. See the description of the [MDFT](#) opcode in the *MacroModel Reference Manual* for more information.

```
sim-ann.mae
sim-ann-out.mae
FFLD      10      1      0      1      1.0000      0.0000      0.0000      0.0000
BDCO       0      0      0      0      41.5692 99999.0000      0.0000      0.0000
READ       0      0      0      0      0.0000      0.0000      0.0000      0.0000
CONV       2      0      0      0      0.0500      0.0000      0.0000      0.0000
MINI       9      0     500      0      0.0000      0.0000      0.0000      0.0000
MDIT       0      0      0      0     300.0000      0.0000      0.0000      0.0000
MDYN       0      1      1      0      1.5000     10.0000     300.0000      0.0000
MDYN       0      1      1      0      1.5000     20.0000     150.0000      5.0000
MDYN       0      1      1      0      2.0000     20.0000      50.0000      5.0000
MINI       9      0     500      0      0.0000      0.0000      0.0000      0.0000
```

MDYN: Note that the three separate MDYN command lines perform different tasks. First, a 10 ps equilibrium run at 300 K (initialized by the [MDIT](#) line) is performed, followed by a 20 ps simulation coupled to a thermal bath of 150 K. Arg8 gives the time constant of the bath (in ps). This second MDYN command slowly cools the system to 150 K, while the last MDYN command cools the system to 50 K.

A WRIT command is not necessary because the final MINI command ensures that the minimized structure is written to the output file.

See [Section 12.3.1 on page 136](#) for a description of the other opcodes in this file.

12.4 Checking and Interpreting Results

It is often difficult to know when you have sampled enough in a molecular dynamics simulation. It usually is helpful to have some prior knowledge of the time-scales for crucial processes within the system to know if convergence may have been achieved.

While molecular dynamics simulations provide a wealth of structural and temporal information, the information is often hard to interpret. If the goal is to thermally sample the local conformational minimum of a small molecule, then you might not need to examine the results carefully provided that you have simulated sufficiently long (50 to 100 ps may be enough). However, if the system is large and the conformational variation is large, simulations likely will not adequately sample the conformations available, and careful, problem-specific consideration of the results may be needed to learn from such studies. In such circumstances, it almost always helps to examine the trajectory visually with a tool such as Maestro's ePlayer. In addition, clustering tools such as XCluster may help you identify when key events occurred during the simulation.

MC/SD Calculations

Monte Carlo/Stochastic Dynamics (MC/SD) performs constant temperature calculations that take advantage of the strengths of Monte Carlo methods for quickly introducing large changes in a few degrees of freedom, and stochastic dynamics for its effective local sampling of collective motions. MC/SD is a good general choice for studying a system at constant temperature.

13.1 The MC/SD Panel

You can prepare, write job files for, and submit a Monte Carlo/Stochastic Dynamics calculation from the Maestro MC/SD panel. The MC/SD panel has a general setting section and a Potential folder, a Constraints folder, and a Substructure folder like other MacroModel energy panels. These portions of the panels are described in detail in [Section 5.1 on page 53](#) through [Section 5.5 on page 61](#). The MC/SD panel also has Monitor and Dynamics folders. These folders also appear on the Dynamics panel, and are discussed in [Section 12.2 on page 129](#). The MCSD folder is unique to the MC/SD panel. This folder contains reaction condition settings.

To open the MC/SD panel, choose MC/SD from the MacroModel submenu of the Applications menu in the main menu bar.

13.2 Setting Up MC/SD Calculations

The MC/SD procedure [\[29\]](#) differs from a normal dynamics simulation in that it uses a mixture of Metropolis Monte Carlo and dynamics steps in order to greatly increase the rate at which a simulation explores conformational space. For MC/SD simulations, you must specify torsions to be rotated and, if there is more than one molecule in the system, molecules to be translated and rotated.

The controls for MC/SD simulations are in the MCSD folder, with the exception of the temperature setting, which is in the Dynamics folder.

To define the torsions to be rotated and the molecule translation and rotation selections automatically, click the Perform Automatic Setup button. To view and edit the automatically generated settings, or to specify the settings manually, open the Torsion Rotations and Molecule Trans/Rot panels by clicking the corresponding buttons. These panels are the same as for conformational searches. For a more detailed description of these panels, see [Section 10.2.5.2 on page 106](#) and [Section 10.2.5.3 on page 107](#).




Figure 13.1. The MCSD folder of the MC/SD panel.

The value entered in the Ratio of SD to MC steps text box determines how many Monte Carlo trials will be performed for each dynamics time step. The default value is 1, but the number can be increased to give more dynamics steps per Monte Carlo step.

You can specify the minimum and maximum number of torsions to vary in an MC step in the corresponding text boxes. The number of torsions varied is randomly selected from the range defined by the minimum and maximum values. The default is to vary only one torsion.

We recommend against the use of SHAKE during MC/SD simulations, as it can result in less than optimal temperature control.

13.3 Command File Examples

To run MacroModel calculations, a molecular structure file and a command file are required. The molecular structure file contains the structures to be used as input in the calculation. The

command file contains the name of the input structure file, the name of the output structure file, and an ordered list of operation codes (*opcodes*) for the calculations.

Once you set up a job in Maestro and click either **Start** or **Write**, Maestro writes out a molecular structure file and a command file. For many types of jobs, command files written this way are complete and adequate, but for some types of jobs you may need to adjust the Maestro-generated command file.

This section contains an example of an MC/SD computation using a small organic molecule. The command files and the log files for the examples given in this section can be found in `$SCHRODINGER/macromodel-vversion/samples/Examples`.

```
mcsd.mae
mcsd-out.mae
MMOD      0      1      0      0      0.0000      0.0000      0.0000      0.0000
FFLD     10      1      0      1      1.0000      0.0000      0.0000      0.0000
BDCC      0      0      0      0      41.5692 99999.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
CONV      2      0      0      0      0.0500      0.0000      0.0000      0.0000
MINI      1      0      500    0      0.0000      0.0000      0.0000      0.0000
MCNV      1      3      0      0      0.0000      0.0000      0.0000      0.0000
TORS      1     14      0      0      0.0000     180.0000      0.0000      0.0000
TORS      6     26      0      0      0.0000     180.0000      0.0000      0.0000
TORS     14     15      0      0      0.0000     180.0000      0.0000      0.0000
TORS     15     16      0      0      0.0000     180.0000      0.0000      0.0000
TORS     16     17      0      0      0.0000     180.0000      0.0000      0.0000
TORS     17     18      0      0      0.0000     180.0000      0.0000      0.0000
TORS     26     27      0      0      0.0000     180.0000      0.0000      0.0000
MCSD      1      0      0      0      0.0000      0.0000     300.0000      0.0000
MDIT      0      0      0      0     300.0000      0.0000      0.0000      0.0000
MDYN      0      0      1      0      1.5000      1.0000     300.0000      0.0000
MDSA     20      0      0      0      0.0000      0.0000      1.0000      0.0000
MDYN      1      0      1      0      1.5000     10.0000     300.0000      0.0000
WRIT      0      0      0      0      0.0000      0.0000      0.0000      0.0000
```

Several opcodes have been discussed in previous chapters. Only the opcodes relevant to the MC/SD procedure are described below.

MINI: Minimize the structure before the dynamics simulation.

MCNV: The number of torsion angles to vary during the MC portion of the simulation. Here, between one and three torsion variations are specified.

TORS: Specify the torsion angles to vary by atom number.

MCSD: Specify MC/SD sampling at 300 K with an MC to SD step ratio of 1:1.

MDIT: Specify dynamics simulation at a temperature of 300 K.

MDYN: Perform equilibration dynamics run for 1 ps.

MDSA: Sample 20 structures during the following full dynamics run.

MDYN: Perform full MC/SD dynamics run for 10 ps at 1.5 fs intervals.

13.4 Checking and Interpreting Results

As with normal dynamics simulations, it is often difficult to know when you have sampled enough in a MC/SD simulation. It is usually helpful to have some prior knowledge of the time scales for crucial processes within the system to know if convergence may have been achieved. While MC/SD simulations provide a wealth of structural and temporal information, this information is often hard to interpret. If the goal is to thermally sample the local conformational minimum of a small molecule, then you may not need to examine the results carefully provided that you have simulated sufficiently long (50 to 100 ps may be enough). However, if the system is large and the conformational variation is large, simulations likely will not adequately sample the conformations available, and careful, problem-specific consideration of the results may be needed to learn from such studies. In such circumstances, it almost always helps to examine the monitored trajectory with a tool such as Maestro's ePlayer. In addition, clustering tools such as XCluster may help identify when key events occurred during the simulation.

Minta Calculations

MINTA can be used for fast computation of the conformational free energy of small- and medium-sized molecular systems *in vacuo* and in the presence of a continuum solvent model. MINTA is an excellent tool for calculating the binding free energy of molecular complexes composed of substrate molecules bound to small receptors used in the molecular recognition field or enzyme receptor models used in pharmaceutical research. Unlike available free energy simulation programs, MINTA calculations are user-friendly and are a simple tool for medicinal chemists familiar with conformational analysis. For more information on MINTA calculations, see [Appendix H](#) of the *MacroModel Reference Manual*.

14.1 The MINTA Panel

You can run the MINTA program from the MINTA panel using the Workspace contents, project table entries, or structures in a separate file as input. The MINTA panel has a general setting portion, which is described in [Section 5.1 on page 53](#). The panel also contains four tabbed folders: Potential, Substructure, Mini, and MINTA. For information on settings in the Potential folder, see [Section 5.2 on page 54](#). For substructure-related material, see [Section 5.5 on page 61](#). For the Mini folder description, see [Section 7.2 on page 73](#).

To open the MINTA panel, choose MINTA from the MacroModel submenu of the Applications menu in the main menu bar. The controls in the MINTA folder are described below.

Input File

Use the first part of the panel to enter the name of the input file. This must be a Maestro-formatted file that contains one or more valid conformations. These structures are usually the results of a previous conformational search. The Open button opens a file selector for locating the desired input file. If no file is entered, the structural input is taken from the source indicated under Source of job input.

Number of MINTA iterations

The MINTA numerical integrals are calculated in statistical blocks to achieve better convergence. The number of blocks used is referred to as the number of MINTA iterations. The default value is 5, and the minimum value is 1.




Figure 14.1. The MINTA folder of the MINTA panel.

Number of energy evaluations per MINTA iteration

MINTA integration is based on single point energy evaluations. The default number of energy evaluations per MINTA iteration is 2000, and the minimum value is 1.

Temperature (K)

The value entered in this text box sets the temperature for the MINTA calculations. The default value is 300 K, and the minimum value is 0 K.

Hard limit for sampling along normal modes (Å)

This value determines a distance to which sampling is limited from the equilibrium geometry of the structure along any of the normal mode directions in $3N-6(5)$ dimensional normal mode space, where N is the number of atoms. The default value is 1.0 Å, the minimum value is 0.0 Å, and the maximum value is 3.0 Å.

Soft limit for sampling along normal modes (st. dev)

This value determines the units of standard deviation for sampling along normal modes. Sampling is limited to different distances from the equilibrium geometry along different normal mode directions. For a particular mode i sampling is limited to a particular distance equal to this value times the standard deviation of the multidimensional Gaussian function, along the particular normal mode direction i .

14.2 Command File Examples

To run MacroModel calculations, a molecular structure file and a command file are required. The molecular structure file contains the structures to be used as input in the calculation. The command file contains the name of the input structure file, the name of the output structure file, and an ordered list of operation codes (*opcodes*) for the calculations.

Once you set up a job in Maestro and click either **Start** or **Write**, Maestro writes out a molecular structure file and a command file. For many types of jobs, command files written this way are complete and adequate, but for some types of jobs you may need to adjust the Maestro-generated command file.

The command file and the log file for the example given in this section can be found in `$SCHRODINGER/macromodel-vversion/samples/Examples`.

MintaMae.mae

MintaMae-out.mae

MMOD	0	1	0	0	0.0000	0.0000	0.0000	0.0000
FFLD	10	1	0	1	1.0000	0.0000	0.0000	0.0000
BDCO	0	0	0	0	41.5692	99999.0000	0.0000	0.0000
BGIN	0	0	0	0	0.0000	0.0000	0.0000	0.0000
READ	0	0	0	0	0.0000	0.0000	0.0000	0.0000
CONV	2	0	0	0	0.0500	0.0000	0.0000	0.0000
MINI	9	0	500	0	0.0000	0.0000	0.0000	0.0000
MNTA	5	2000	0	0	300.0000	1.0000	3.0000	0.0000
END	0	0	0	0	0.0000	0.0000	0.0000	0.0000

MMOD: Creates and updates an intermediate structure file so that structures can be displayed in Maestro as the job progresses.

FFLD: Force field selection. Arg1 denotes the actual force field used in the calculation (in this case MMFF94). Arg2 defines the electrostatic treatment for the calculation. The default (arg2=0) is to use the dielectric treatment encoded in the force field. However, in this case a constant dielectric is used. Arg4 is MMFF94-specific. Arg4=1 defines the MMFF94s version of the force field, ensuring planarity around delocalized sp² nitrogens.

BDCO: Use the Bond Dipole CutOff (BDCO) method for truncating electrostatic interactions. Arg5 and arg6 are used to specify the cutoffs used for charge-dipole and charge-charge interactions, respectively.

BGIN/END: The loop reads structures from a preceding conformational search. The input file should contain only conformers of the same molecule or molecular complex.

READ: Read the input file.

CONV: Defines convergence criteria. Arg1=2 signifies derivative convergence (default criterion, if no CONV command is present, is 0.05kJ/mol-Å. This value is set in arg5).

MINI: Starts the minimization. Arg1 defines the type of minimization algorithm to be used. Arg1=9 means that Truncated Newton-Raphson Conjugate Gradient will be used. In arg3, the number of minimization steps is defined. Arg3 can be set to a large number since the calculation will automatically stop as soon as the convergence criterion is reached.

MNTA: A Minta free energy calculation will be performed. Minta numerical integrations are performed in blocks in order to achieve better convergence, and arg1 defines the number of such blocks. Arg2 gives the number of energy evaluations per block, hence the total number of energy evaluations per structure is $\text{arg1} * \text{arg2}$. Minta can be run in an adaptive manner, which is slightly more accurate than the non-adaptive mode. Arg3=0 specifies the use of non-adaptive mode integration. Set arg3 not equal to 0 and arg4 to values from 1 to 10 for adaptive integration. Arg4 defines the number of “soft” (or “low-frequency” vibrational modes) degrees of freedom for which numerical integration should be applied. The default is to apply numerical integration to all degrees of freedom. However, this is recommended only for very small molecules (about 20 atoms or fewer). It is strongly recommended not to use values of arg4 greater than 50. Arg5 sets the simulation temperature, while arg6 and arg7 define hard and soft limits for sampling along normal modes.

14.3 Checking and Interpreting Results

For many MINTA calculations, the input set of conformations must be generated in a separate calculation such as a conformational search. The set of conformations needed for accurate MINTA calculations should be fairly extensive and include all low-lying conformers.

Protein Loop Construction

Before analyzing how a given protein-ligand pair will interact, you may want to refine the protein structure. This commonly involves examining the structure of loops—short sequences of amino acids which typically occur on the surface of the protein, but in the middle of the protein sequence.

The reasons for refining a protein structure include:

- The structure of this portion of the protein was not well-resolved.
- The conformation of the loop changed upon solvation.

In addition, when working with enzymes, you may want to examine how the structure of the loop changes as you change the sequence in the loop.

The MacroModel LOOP tool is a conformational search method for protein loops. It generates a variety of structures that are both diverse and geometrically appropriate for the protein loop, and minimizes the energy of these structures. LOOP functions in much the same manner as MCMM, and is used in conjunction with the [MINI](#), [MCOP](#), [MSYM](#), and [LPOP](#) opcodes. LOOP is best used in conjunction with other methods, like LLMOD, to sample the loop structures more thoroughly. For more information on LOOP, see the [MacroModel Reference Manual](#).

15.1 Performing LOOP Calculations

When performing a LOOP calculation, you can indicate the loop sequence you want to examine in one of two ways:

- Indicate that you want to use the loop sequence is that in the original protein.
- Create an auxiliary *filename.lsq* file, which lists the desired loop sequence.

15.1.1 Input Restrictions

For the current implementation of LOOP, input must meet the following criteria:

- Only one loop can be examined at a time, although successive calculations may focus on different loops.
- The protein loop and the residues it is immediately attached to must consist of amino acids containing the backbone atoms Np-C α -Cp, where Np and Cp are nitrogen and carbon atoms involved in peptide bonds.

- An all-atom representation of the protein must be used.
- No disulfide bonds are permitted within the loop or between the loop and the rest of the protein.
- No atoms in the loop can be frozen.
- LOOP will not work from a collection of input protein structures. It will work only from one such structure.

15.1.2 Structure Handling Considerations

When using LOOP, you should be aware of the following:

- Rings, such as the one in proline, are treated as rigid.
- If comparison atoms ([COMP](#)) are specified explicitly or implicitly (see [arg3](#)), then [MSYM](#) must be used.
- All LOOP runs renumber the atoms within the structure such that the loop atoms become the highest numbered atoms in the structure. During a run, MacroModel tracks the shifted atoms numbers for the atoms in the system. LOOP automatically produces a substructure file, *filename-out.sbc*, containing the shifted atom numbers for use in follow-up studies.

15.2 Example .lsq Input File

Consider the following contents of an example auxiliary file, *filename.lsq*.

```
GLY
ALA
SER
THR
ARG
GLU
SER
```

Based on the file contents, LOOP creates a new loop with the seven specified amino acids, starting with GLY at the N-terminus end of the loop.

The format and contents of the file meet the following criteria:

- Only alpha amino acids are used.
- The amino acids used all appear in the Maestro fragment tables, found in the Maestro Build panel.

- Standard three-letter abbreviations are used.
- Nonstandard (D) alpha amino acids from the Maestro fragment tables can also be used, but these require the four-letter abbreviation, e.g., DALA.

15.3 Command File Examples

To run MacroModel calculations, a molecular structure file and a command file are required. The molecular structure file contains the structures to be used as input in the calculation. The command file contains the name of the input structure file, the name of the output structure file, and an ordered list of operation codes (*opcodes*) for the calculations.

This section includes example command files for using LOOP, when using a LOOP sequence from an input protein structure and using the LOOP sequence from a file. The command files and the log files for the examples given in this section can be found in `$SCHRODINGER/macromodel-vversion/samples/Examples`.

15.3.1 LOOP Job Using the Input Structure Sequence

An example command file for a LOOP sequence from an input protein structure appears below. Descriptions of the opcodes used in the file follow.

```
loop_input.mae
loop_input-out.mae
SOLV      3      1      0      0      0.0000      0.0000      0.0000      0.0000
FFLD      11      0      0      0      0.0000      0.0000      0.0000      0.0000
EXNB      0      0      0      0      0.0000      0.0000      0.0000      0.0000
SUBS      0      0      0      0      0.0000      0.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
LOOP      423      467      0      4      0.0      0.26      0.0000      0.0000
DEMX      0      0      0      0      1000.0      0.0000      0.0000      0.0000
MSYM      0      1      0      0      0.0000      0.0000      0.0000      0.0000
MCOP      1      0      0      0      0.0000      0.0000      0.0000      0.0000
CONV      2      1      0      0      0.0000      0.0000      0.0000      0.0000
MINI      9      0      2000      0      0.0000      0.0000      0.0000      0.0000
```

SOLV: The GB/SA effective water model is being used.

FFLD: This example uses OPLS_2001 with constant dielectric electrostatics, which is appropriate if GB/SA solvation is being used.

EXNB: Extended non-bonded cut-offs should be used with GB/SA solvation.

SUBS: Read in a substructure from an `.sbc` file. Currently, best results are obtained if all the non-loop atoms in the system are frozen. A shell containing complete residues for all atoms within 6 Å is usually sufficient.

READ: Read in the protein structure.

LOOP: Turn on LOOP.

arg1: N-terminus atom number for the loop. This is the atom number for the peptide Nitrogen atom that joins the loop to the protein at the N-terminus of the loop. Atom # 423 for this example.

arg2: C-terminus atom number for the loop. This is the atom number for the peptide Carbon atom that joins the loop to the protein at the C-terminus of the loop. Atom # 467 for this example.

arg3=0: Use the sequence for the loop from the protein structure supplied.

arg4=4: Generate four candidate loop structures. Typically you would want to generate many more than this.

arg5=0.0: Save up to the default number of structures (10,000).

arg6=0.26: If a given pair of atoms is closer than arg6 times the sum of their van der Waals radii after loop generation just prior to minimization then the structure is rejected and another loop structure is generated. Default: 0.25.

arg7=0.0: Automatically add all heavy atoms in the loop to the comparison list used in identifying conformations that are identical (COMP) entries. Also automatically add all chiral atoms in the loop to the list of chiral atoms that need to be checked (CHIG entries) before accepting a generated conformer.

DEMEX: arg5=1000.0: Keep conformers that are up to 1000 kJ/mol higher in energy than the lowest energy conformers. Such a large value might be appropriate if a follow-up study using a conformational search method like LLMOD were going to be conducted.

MSYM: arg2=1: Use mmsym to compare conformers. The comparison is done in-place, that is, translations and rotations of the protein as a whole are not permitted when comparing structures.

MCOP: Print messages to the log file concerning every conformation generated.

CONV: Minimizations are converged when the RMS gradient is less than 0.1 kJ/mol-Å.

MINI: Minimize the structure using the TNCG minimizer.

This LOOP run also writes out COMP and CHIG command lines containing the shifted atom numbers to the .log file for explicitly and implicitly specified COMP atoms and CHIG atoms in the system. These lines could be inserted into the command files for subsequent studies.

15.3.2 LOOP Job Using the Sequence From a File

This example is identical to the previous command file except that `arg3` of `LOOP` is 1, indicating that a *filename.lsq* is used to specify the sequence for the loop.

15.4 Example LOOP Job Output

15.4.1 Atom Renumbering Using Input Protein Structure

The excerpt below shows the atom renumbering formation from the `.log` file generated during a run in which the loop sequence from the actual protein structure was used as input.

```
LOOP being generated using sequence from input structure
Note Substructure information may have changed.
Recording Substructure information in: MCPC603-out.sbc
```

Beginning of new list of COMP and CHIG commands

COMP	6607	6608	6609	6610	0.0000	0.0000	0.0000	0.0000
COMP	6611	6612	6613	6614	0.0000	0.0000	0.0000	0.0000
COMP	6615	6616	6617	6618	0.0000	0.0000	0.0000	0.0000
COMP	6619	6620	6621	6622	0.0000	0.0000	0.0000	0.0000
COMP	6623	6624	6625	6626	0.0000	0.0000	0.0000	0.0000
COMP	6627	6628	6629	6630	0.0000	0.0000	0.0000	0.0000
COMP	6631	6632	6633	6634	0.0000	0.0000	0.0000	0.0000
COMP	6635	6636	6637	6638	0.0000	0.0000	0.0000	0.0000
COMP	6639	6640	6641	6642	0.0000	0.0000	0.0000	0.0000
COMP	6643	6644	6645	6646	0.0000	0.0000	0.0000	0.0000
COMP	6647	6648	6649	6650	0.0000	0.0000	0.0000	0.0000
COMP	6651	6652	6653	6654	0.0000	0.0000	0.0000	0.0000
CHIG	6612	6617	6623	6626	0.0000	0.0000	0.0000	0.0000
CHIG	6630	6641	6650	0	0.0000	0.0000	0.0000	0.0000

End of new list of COMP and CHIG commands

15.4.2 Conformational Search Using Input Protein Structure

The excerpt below is from the same `.log` file as the excerpt above. This `.log` file was generated using the loop sequence from the actual protein structure.

```
Step      1 New global minimum.  E (kJ/mol) =    -3830.85
Conf     1 E = -3830.85 (  0.057) is unique and stored as structure 1
Search initialized with 1 structures from the .dat file

Conf     2 E = -3668.40 (  0.077) is unique and stored as structure 2
Conf     3 E = -3652.38 (  0.043) is unique and stored as structure 3
Conf     4 E = -3630.91 (  0.096) is unique and stored as structure 4
```

```
Final report:
  4 unique conformations found
  4 minimized with good convergence
Found      1 confs within  1.00 kcal/mol ( 4.18 kJ/mol) of glob. min.
Global minimum E = -3830.85 found      1 times.

BatchMin normal termination
Total number of structures processed = 4
Conformations with poor convergence marked with a *
Conformation      1 (  -3830.854      kJ/mol) was found      1 times
Conformation      2 (  -3668.398      kJ/mol) was found      1 times
Conformation      3 (  -3652.378      kJ/mol) was found      1 times
Conformation      4 (  -3630.911      kJ/mol) was found      1 times
      *** MC Statistics ***
Percent of minimized structures within energetic window: 100.0000000
Average number of duplicates:  1.000000000
Duplication standard deviation: 0.0000000000E+00
5  structures generated
0  rejected by ring closure
2  rejected by van der Waals
0  duplicate minimised structures
      Time in Monte Carlo generation loop:      5.8 CPU sec
      Time in energy minimizations:      153.4 CPU sec
      Time in geometry optimisation:      0.0 CPU sec
```

15.4.3 LOOP Run Using an .lsq File as Input

The .log file for LOOP runs that use auxiliary *filename.lsqa* files looks essentially the same as those in the previous two examples, except that when an .lsqa file is used, the following statement appears in the .log file output:

```
LOOP being generated using sequence in .lsqa file
```

eMBrAcE

You can obtain a set of ligands that have been pre-positioned with respect to a receptor from various sources, including Schrödinger's docking program, Glide. To study the association of the ligands with the receptor further, you can use the automated mechanism of Multi-Ligand Bimolecular Association with Energetics (eMBrAcE). With eMBrAcE, complexes can be studied using simple minimizations or conformational searches.

16.1 Minimizations With eMBrAcE

An eMBrAcE minimization is a type of multiple minimization in which each of the specified pre-positioned ligands is minimized, in turn, with the receptor. You can perform eMBrAcE minimization calculations in two modes: Interaction Mode, in which the interaction between each ligand and the substrate is studied, and Energy Difference Mode, in which energy changes upon association are estimated.

The eMBrAcE Minimization panel is used to set up and submit eMBrAcE minimization jobs. To open this panel, choose eMBrAcE Minimization from the MacroModel submenu of the Applications menu in the main menu bar. The upper and lower parts of the panel and the Potential and Substructure folders are common to all MacroModel panels. These components are described in detail in [Section 5.1 on page 53](#) through [Section 5.5 on page 61](#). The Mini folder is common to many of the MacroModel panels. For an explanation of the controls in this folder, see [Section 7.1 on page 73](#). The controls for the eMBrAcE minimization settings are located in the eMBrAcE folder.

To perform an eMBrAcE minimization, first configure the general job settings in the upper portion of the eMBrAcE Minimization panel and the Potential and Substructure folders as described in [Chapter 5](#). Because the receptor is usually large, you should consider defining substructures and fixing or freezing atoms that are far from the active site, to speed up the calculations. Then, configure the eMBrAcE settings in the eMBrAcE folder, discussed below.

Source of ligands

You can select ligands for an eMBrAcE calculation from the entries in the Project Table, or you can read the ligands from a file. To specify the file, you can enter the path to the file in the text box, or you can click Browse and navigate to the file.




Figure 16.1. The eMBrAcE Minimization *panel showing the eMBrAcE folder.*

Receptor

For an eMBrAcE calculation to execute correctly, it must correctly identify which structure is to be treated as the receptor.

If you are reading the ligands from file, you can nominate the first structure in the file as the receptor, or you can specify it as an entry in the current project. If the receptor is not in the file, you must specify it as an entry in the current project. You can either type the entry name in the text box or click Choose and select an entry in the Choose Entry dialog box that is displayed. In this case, the receptor is written to a separate file for the job.

If the ligands are entries in the current project, you can nominate the first selected entry as the receptor, or you can choose another entry for the receptor. To choose an entry, type the entry name in the text box or click Choose and select an entry in the entry selection dialog box that is displayed.

Association energy mode

The eMBrAcE calculation can function in two modes. The first mode is Energy difference mode. In this mode, the calculation is performed first on the receptor, then on the ligand, and finally on the complex. The energy difference is then calculated using the equation:

$$\Delta E = E_{complex} - E_{ligand} - E_{protein}$$

The full effects of relaxation and solvation are included in this mode.

The second is Interaction energy mode. In this mode, the atoms in the ligand and the receptor are separated into two sets, and the interaction energy between the two sets is calculated. Interaction energy mode deals with terms that can be considered pair-wise additive, so the surface energy term in the solvation energies is not included in the interaction energy. Interaction energy mode does not include the relaxation or the change in solvation of the ligand on binding.

To specify a preference for calculation of the association energy, select either Energy Difference Mode or Interaction Energy Mode.

Structures saved

The three options for how output structures are written to the output structure file are: Complexes only (minimized complexes); Ligands only (ligand structures extracted from minimized complexes); and Receptor, ligands, and complexes, in which the receptor minimized without a ligand and ligands minimized without the receptor are written. For interaction mode, the last option is equivalent to Complexes only.

16.2 Conformational Searches With eMBrAcE

eMBrAcE can perform conformational searches in addition to minimizations. Energy difference mode is the only mode supported for conformational searches with eMBrAcE. Searches are conducted on the receptor, each ligand, and each ligand-receptor complex. The energies used in the energy difference equation

$$\Delta E = E_{complex} - E_{receptor} - E_{ligand}$$

for the receptor and the ligand are the values from the lowest energy conformations found for those systems. However, multiple complex conformations may be retained and a separate $E_{complex}$ is used in the energy difference equation for each one (MCOP arg6 specifies the number of such conformations to keep).




Figure 16.2. The eMBrAcE Conformational Search *panel showing the CSearch folder.*

The eMBrAcE Conformational Search panel is used to set up and submit eMBrAcE conformational search jobs. To open this panel, choose eMBrAcE Conformational Search from the MacroModel submenu of the Applications menu in the main menu bar. The upper and lower parts of the panel and the Potential and Substructure folders are common to all MacroModel panels. These components are described in detail in [Section 5.1 on page 53](#) through [Section 5.5 on page 61](#). The Mini folder is common to many of the MacroModel panels. For an explanation of the controls in this folder, see [Section 7.1 on page 73](#). The controls for the eMBrAcE settings are located in the eMBrAcE folder, which is described in [Section 16.1 on page 153](#). The conformational search parameters are set up in the CSearch folder. This folder contains a subset of the controls found in the CSearch folder for regular conformational searches and is described below.

eMBrAcE conformation searches automatically employ the AUTO Automatic Setup mechanism. AUTO selects the MCMM and comparison atom parameters for each individual ligand-receptor complex. Computations prepared from Maestro have the AUTO opcode added to the job's command file (*jobname.com*) automatically. In addition, AUTO is substructure-aware.

Parameters are only indicated for the freely moving receptor and ligand atoms, but not for fixed or frozen regions of the receptor. Thus, only a proper receptor substructure needs to be indicated to prepare the appropriate MCMM parameters. Examples can be found in [Section 16.3 on page 158](#). In addition, only the substructure facility can be used to indicate fixed or frozen atoms in an eMBrAcE calculation, and ligand atoms cannot be fixed or frozen. Other than by adjusting the substructure, it is not possible to adjust the AUTO parameters for individual complexes in an eMBrAcE automated calculation.

Conformational searches on protein-ligand complexes are computationally intensive, and the use of substructures is strongly recommended to reduce the CPU time and memory required. Even with fairly small substructures and very short searches, such searches require hours or days for each ligand processed.

To perform an eMBrAcE conformational search, first configure the general job settings in the upper portion of the eMBrAcE Minimization panel and the Potential and Substructure folders. Then configure the eMBrAcE settings in the eMBrAcE folder and the conformational search settings in the CSearch folder, which are described below.

Method

There are only three methods available for performing an eMBrAcE conformational search: Torsional sampling (MCMM), Low-mode sampling, and Mixed torsional/Low-mode sampling. Large scale low-mode (LMC2) searches are not supported at this time. Low-mode calculations are limited to a few hundred atoms and should be used for systems with smaller ligand/substructure combinations. However, for such systems, low-mode searches can be comparatively efficient. MCMM conformational searches are less memory-intensive, and thus may be used on larger portions of the system provided that the number of degrees of freedom changed at any one time is limited and the substructure carefully selected.

Number of steps

The number of steps that will be performed in any search is determined by the value entered in this text box. When the number of generated trial structures matches the value given, the conformational search is terminated.

Number of input conformations available to seed each search

The value in this text box specifies the number of conformers available for each ligand to seed the search. The file containing the ligands must have exactly the number of conformers specified for each ligand, otherwise the search will fail when it reaches the ligand that does not have the specified number. Multiple conformers could, for example, come from the use of COPY/ALGN, from a separate conformational search in which a specified number of conformations is kept, or from Glide poses.

Number of structures to save for each search

Specify the number of structures to save for each search, counting from the lowest in energy. A zero value means “save all structures.”

Energy window for saving structures

This is the threshold value for comparison of trial structures. Any new structures generated and minimized are kept only if their energy is less than this value above the current global minimum. Lowering this value results in fewer structures saved. The default value is set to 500 kJ/mol.

Max distance between atoms in equal structures

This text box specifies the threshold for determining whether structures should be considered to be equivalent. When the structures are superimposed, the distance between all pairs of corresponding atoms must be less than this threshold for the structures to be considered equivalent.

Probability of a torsion rotation/molecule translation

Minimum distance for low-mode move

Maximum distance for low-mode move

These three text boxes are relevant only to the low-mode searches and are active only when a method involving low-mode conformational searching is selected.

The first text box is used only with the Mixed torsional/Low-mode sampling method and allows the setting of a probability that any defined torsion rotations and molecule translations are made at each step during the search. This should be a number from 0.0 to 1.0.

The other two text boxes are used for setting the minimum and maximum distance for a low-mode move. During a search, the fastest moving atom is moved randomly generated distances that are between these two limits.

16.3 Specifying a Substructure for eMBrAcE

The use of substructures can dramatically speed up eMBrAcE calculations. Specifying substructures is described in [Section 5.5 on page 61](#). Substructures used in eMBrAcE must meet an additional requirement: the receptor atoms must be numbered starting from 1 when constructing the substructure. Below are two simple recipes for setting up substructures. In the first, all atoms in the receptor are either fixed or frozen. In the second, you can define shells of constrained atoms near the association site in the receptor. If you already have a suitable substructure, you can simply read it in by clicking Read .sbc File near the bottom of the Substructure folder.

Creating a Substructure With All Receptor Atoms Fixed or Frozen:

1. Include the receptor, and only the receptor, in the Workspace.
2. Click New Shell.
3. In the Additional atoms for shell section, click All. By default, all atoms are fixed, that is, constrained to their current positions using a harmonic potential.
4. To freeze the atoms in place, select Freeze atoms.

Using a Ligand to Assist in Creating a Receptor Substructure:

In this example, we will use a ligand from an entry called `our_lig` to set up a substructure for the receptor. The substructure will contain a shell of all receptor atoms from complete residues that lie within 6.0 Å of the ligand. Two additional shells, one fixed and one frozen, are also specified.

1. Include the receptor in the Workspace and then include the ligand.

It is important to include the receptor first to ensure that its atoms are assigned atom numbers starting from 1. You can now create a substructure in the usual manner. Maestro's support for the eMBrAcE utility automatically ignores atoms within the substructure that are numbered higher than the atoms in the receptor so the inclusion of ligand atoms is not a problem when you use the Maestro interface. For each ligand that eMBrAcE is applied to, all atoms in the ligand are automatically added to the substructure without constraints.

2. In the Atoms for substructure section, select Molecule from the Pick menu and pick a ligand atom.
3. Select Complete residues and enter 6.0 in the Expand to atoms within radius of text box.

Alternatively, enter the following ASL expression in the ASL text box:

```
fillres within 6.0 entry.name our_lig
```

Now add a shell of fixed receptor atoms from complete residues within 5.0 Å of the substructure atoms:

4. Click New Shell to create a new shell.
5. In the Radius text box, enter 5.0 and select Complete Residues.

To add another shell containing frozen receptor atoms from complete residues within 4.0 Å of the previously defined shell, follow the same procedure and select Freeze atoms.

16.4 Command File Examples

To run MacroModel calculations, a molecular structure file and a command file are required. The molecular structure file contains the structures to be used as input in the calculation. The command file contains the name of the input structure file, the name of the output structure file, and an ordered list of operation codes (*opcodes*) for the calculations.

Once you set up a job in Maestro and click either **Start** or **Write**, Maestro writes out a molecular structure file and a command file. For many types of jobs, command files written this way are complete and adequate. For some types of jobs, however, including eMBrAcE conformational searches, you will need to adjust the Maestro-generated command file.

The command files and the log files for the examples given in this section can be found in `$SCHRODINGER/macromodel-vversion/samples/Examples`.

16.4.1 eMBrAcE Minimization Calculations

You can set up the files for eMBrAcE minimization calculations using Maestro. This is the recommended mechanism for setting up such calculations. This section explains the structures of eMBrAcE minimization .com files if you want to customize them.

16.4.1.1 Interaction Mode

Interaction mode uses the ASET mechanism to define atom sets and to calculate interaction energies within and among all such sets. All results for sets are saved in the log file for the run. Only interaction energies between sets 1 and 2 are saved as project properties within the output structure file. Note that the surface energy contribution in GB/SA calculations is associated with the interaction of a set with itself.

Below is an example command file for an eMBrAcE job run in interaction energy mode. The input file, `MBAE_Interaction.mae`, follows the Glide pose-viewer pattern, that is, the first structure in the file must be the receptor and the remaining structures should be ligands previously positioned appropriately relative to the receptor. A description of the opcodes used in the file follows, then an excerpt from the output.

```
MBAE_Interaction.mae
MBAE_Interaction-out.mae
SOLV      3      1      0      0      0.0000      0.0000      0.0000      0.0000
FFLD     11      1      0      0      0.0000      0.0000      0.0000      0.0000
EXNB      0      0      0      0      0.0000      0.0000      0.0000      0.0000
BDCE      0      0      0      0      89.4427 99999.0000      0.0000      0.0000
MBAE      0      0      0      0      0.0000      0.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
SUBS      0      0      0      0      0.0000      0.0000      0.0000      0.0000
BGIN      0      0      0      0      0.0000      0.0000      0.0000      0.0000
```

READ	0	0	0	0	0.0000	0.0000	0.0000	0.0000
ASET	0	0	0	0	0.0000	0.0000	0.0000	0.0000
ASET	0	0	0	0	2.0000	0.0000	0.0000	0.0000
ASET	1	4729	0	0	2.0000	-2.0000	0.0000	0.0000
ASET	1	4729	0	0	1.0000	2.0000	0.0000	0.0000
MINI	1	0	2000	0	0.0000	0.0000	0.0000	0.0000
ELST	-1	0	0	0	0.0000	0.0000	0.0000	0.0000
END	0	0	0	0	0.0000	0.0000	0.0000	0.0000
MBAE	-1	0	0	0	0.0000	0.0000	0.0000	0.0000

SOLV: The GB/SA effective water model is being used.

FFLD: Use OPLS_2001 with constant dielectric electrostatics, which is appropriate if GB/SA solvation is being used.

EXNB: Extended non-bonded cut-offs should be used with GB/SA solvation.

BDCO: Use the Bond Dipole CutOff (BDCO) method for truncating electrostatic interactions. Arg5 and arg6 are used to specify the cutoffs used for charge-dipole and charge-charge interactions, respectively.

MBAE: Turn on eMBrAcE.

arg1 – 0: Using Interaction energy mode.

arg2 – 0: MINI calculations being performed.

arg3 – 0: Minimized complex structures are written out.

If arg3 were 1, only the ligand structures extracted from the minimized complex structure would be recorded.

READ: The first READ command obtains the structure of the receptor.

SUBS: Substructures may be used. All opcodes related to the substructure (**SUBS**, **FXAT**, **FXDI**, **FXBT**, and **FXTA**) should be defined in an .sbc file. The .sbc file should contain only lines referring to the receptor with the receptor atoms numbered starting from 1.

BGIN: Start the loop that processes each substructure in turn.

READ: Read in a ligand. This also produces a combined structure consisting of the receptor and the ligand with the receptor atoms being numbered lower than the ligand atoms. All ligand atoms are automatically added to the SUBS list at this stage.

ASET: This series of ASET opcodes places the substrate atoms in set 1 and all other atoms in set 2. If the ligand structure has multiple molecules in it, then this case would place all of the molecules for that structure in set 2. The first ASET line removes all of the atoms from the existing set (by placing them in set 0, which is a dummy set). The second ASET line places all

atoms in set 2. The third and fourth ASET commands remove the receptor atoms (numbered from 1 to 4729) from set 2 and place them in set 1. Other ASET combinations may be useful depending on the system (e.g., you could leave out particular molecules, such as water molecules from the substrate CT). Note that while the interaction energies from all sets are recorded in the .log file, only those for sets 1 and 2 are recorded in the output structure file for eventual inclusion in Maestro's Project Table.

MINI: Minimize the energy of the structure using the PRCG minimizer.

ELST: Invokes the ASET mechanism for calculating interaction energies. Arg1 = -1 causes limited information to be stored in the .log file only. Other values may produce very large .mmo files.

END: End the loop over ligands.

MBAE: When arg1 = -1, MBAE is turned off. This command is not needed in this example.

Below is the data generated by eMBrAcE for one ligand. The complete output file contains one such set for each ligand in the analysis set. Note that Conf 2 refers to the first ligand in the file, which is the second structure in the input structure file. Note also that only the energetic results between sets 1 and 2 are recorded in the output structure file.

A table is given at the end of the log file containing a summary of the results for all the ligands.

```
Conf      2 E=-14227.220 ( 0.048) kJ/mol
Using numerical surfaces and analytical Born radii
Solvation GB set energies do not include constant
contributions exclusively involving fixed/frozen atoms.
Energetic Interactions Within Atom Sets (with no. of interactions):
Atom set   1:
    Total Energy (kJ/mol) =-0.5221442D+05 ( 2513327 )
        Stretch = 0.1725957D+03 ( 3133 )
        Bend = 0.1620882D+03 ( 659 )
    Proper Torsion = 0.2545627D+03 ( 779 )
    Out-of-Plane = 0.6995831D+01 ( 90 )
    Electrostatic =-0.6687509D+04 ( 444935 ) (Part of nonbonded)
    Van der Waals =-0.6857149D+03 ( 61959 ) (Part of nonbonded)
    Solvation SA = 0.1683423D+03 ( 637 )
    Solvation GB =-0.4576560D+05 ( 1999662 )
    Fixed-atom constraint = 0.1598148D+03 ( 1473 )
    Non-bonded =-0.7373224D+04 ( 506894 ) (Elect + Hbnd + vdW)
Atom set   2:
    Total Energy (kJ/mol) =-0.1137470D+03 ( 3040 )
        Stretch = 0.3625181D+01 ( 42 )
        Bend = 0.4390021D+02 ( 70 )
    Proper Torsion = 0.3172877D+02 ( 88 )
    Out-of-Plane = 0.2464711D+01 ( 14 )
    Electrostatic =-0.2183937D+03 ( 668 ) (Part of nonbonded)
```



```

Van der Waals = 0.5127575D+02 (      543 ) (Part of nonbonded)
Solvation SA = 0.1506673D+01 (      15 )
Solvation GB =-0.2985466D+02 (     1600 )
Fixed-atom constraint = 0.0000000D+00 (      0 )
Non-bonded =-0.1671179D+03 (     1211 ) (Elect + Hbnd + vdW)
Energetic Interactions Between Atom Sets (with no. of interactions):
Atom sets  1 and  2:
Use these for eMBrAcE interaction energies for ligand: 4erk
Total Energy (kJ/mol) =-0.2295194D+03 (    212457 )
Stretch = 0.0000000D+00 (      0 )
Proper Torsion = 0.0000000D+00 (      0 )
Out-of-Plane = 0.0000000D+00 (      0 )
Electrostatic =-0.9417975D+02 (     65023 ) (Part of nonbonded)
Van der Waals =-0.1651001D+03 (     9236 ) (Part of nonbonded)
Solvation GB = 0.2976043D+02 (    138198 )
Non-bonded =-0.2592798D+03 (     74259 ) (Elect + Hbnd + vdW)

```

16.4.1.2 Energy Difference Mode

Energy difference mode first applies the type of calculation selected (currently only minimization is supported) to the receptor by itself. For each ligand in turn, it then applies the calculation to the ligand by itself and the complex of the ligand with the receptor, and reports the energy difference.

Below is an example command file for an eMBrAcE job run in energy difference mode. A description of the opcodes used in the file follows.

```

MBAE_eDiff.mae
MBAE_eDiff-out.mae
SOLV      3      1      0      0      0.0000      0.0000      0.0000      0.0000
FFLD     11      0      0      0      0.0000      0.0000      0.0000      0.0000
EXNB      0      0      0      0      0.0000      0.0000      0.0000      0.0000
BDCE      0      0      0      0      89.4427 99999.0000      0.0000      0.0000
MBAE      1      0      1      0      0.0000      0.0000      0.0000      0.0000
SUBS      0      0      0      0      0.0000      0.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
MINI      1      0    2000      0      0.0000      0.0000      0.0000      0.0000
BGIN      0      0      0      0      0.0000      0.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
MINI      1      0    2000      0      0.0000      0.0000      0.0000      0.0000
END       0      0      0      0      0.0000      0.0000      0.0000      0.0000
MBAE     -1      0      0      0      0.0000      0.0000      0.0000      0.0000

```

SOLV: The GB/SA effective water model is being used.

FFLD: This example uses OPLS_2001 with constant dielectric electrostatics, which is appropriate if GB/SA solvation is being used.

EXNB: Extended non-bonded cut-offs should be used with GB/SA solvation.

MBAE: Turn on eMBrAcE.

arg1 = 1: Using Energy Difference mode.

arg2 = 0: MINI calculations being performed.

arg3 = 1: Ligand structures extracted from the minimized complex are written to the output structure file.

If arg3 were set to 0, the minimized complex structure would be recorded.

SUBS: Substructures may be used. All opcodes related to the substructure should be defined in an .sbc file. The .sbc file should contain only lines referring to the receptor with the receptor atoms numbered starting from 1.

READ: The first READ command obtains the structure of the receptor.

MINI: Minimize the energy of the receptor. Arg1=1 uses the PRCG minimization technique for up to arg3 (2000) iterations.

BGIN: Start the loop which processes each ligand in turn.

READ: Sequentially read the ligand structures.

MINI: Minimize first the current ligand structure, then minimize the receptor-ligand complex. Record the energy differences to the *.log file and to the output structure file.

END: End the loop for ligands.

MBAE: When arg1 = -1, MBAE is turned off. In this example, this command is not needed.

Below is the eMBrAcE output from the .log file for one ligand in a job run in energy difference mode, showing the energy differences upon complexation.

```
=====
MBAE Energy Differences for Ligand # 1
Ligand Name: 4erk
  Total Energy(constraints) =-0.8012573E+02 kJ/mol
  Total Energy(no constraints) =-0.8160452E+02 kJ/mol
    Valence Energy = 0.1018281E+02 kJ/mol
      VDW =-0.1592237E+03 kJ/mol
        Electrostatic =-0.9869821E+02 kJ/mol
          Explicit Hydrogen Bonds = 0.0000000E+00 kJ/mol
            Solvation = 0.1661343E+03 kJ/mol
              Constraint = 0.1478790E+01 kJ/mol
                Converged =      T
```

In this example, `Converged = T` means that all minimization converged. If the value were `F` instead of `T`, then the energy differences would be suspect and the calculation may need to be repeated with more minimization steps.

At the end of the log file is a table summarizing the results for all the ligands.

16.4.2 Conformational Searches With eMBrAcE

Below are example `.com` files for an MCMM and an LMOD conformational search with eMBrAcE. After these an example is provided using `COPY` and `ALGN` in conjunction with an MCMM/MBAE conformational search.

Like eMBrAcE minimization calculations, eMBrAcE conformational searches include a table of results at the end of the `.log` file. Energetic properties are also included in this output structure file.

eMBrAcE conformational searches can require very large amounts of computer time and memory. The memory requirements are similar to the equivalent non-eMBrAcE search of the complex of the protein with the largest ligand. Because of these requirements, for nearly all receptors it will likely be necessary to specify a substructure to reduce the required resources to a practical level. If you want to perform only the eMBrAcE search on the ligand within an essentially static receptor, then the substructure used should not include `SUBS` lines but only list the receptor atoms as fixed or frozen (see the `FXAT` opcode description for more information) and use only one step for the first MCMM or LMOD line (the one for the search of the receptor).

16.4.2.1 MCMM Conformational Search

This example uses Monte Carlo Multiple Minimum searching. A copy of this file is available at:

```
$SCHRODINGER/macromodel-vversion/samples/Examples/MBAE_MCMM.com
SEARCH_MCMM.mae
SEARCH_MCMM-out.mae
FFLD      11      1      0      0      1.0000      0.0000      0.0000      0.0000
SOLV       3      1      0      0      0.0000      0.0000      0.0000      0.0000
EXNB       0      0      0      0      0.0000      0.0000      0.0000      0.0000
BDCC       0      0      0      0      89.4427 99999.0000      0.0000      0.0000
MSYM       0      0      0      0      0.0000      0.0000      0.0000      0.0000
MBAE       1      1      0      0      0.0000      0.0000      0.0000      0.0000
MCMM      50      0      0      0      0.0000      0.0000      0.0000      0.0000
MCNV       1      4      0      0      0.0000      0.0000      0.0000      0.0000
MCSS       2      0      0      0      500.0000      0.0000      0.0000      0.0000
MCOP       1      0 10000      0      0.0000      1.0000      0.0000      0.0000
DEMX       0      0      0      0      500.0000      0.0000      0.0000      0.0000
```

SUBS	0	0	0	0	0.0000	0.0000	0.0000	0.0000
READ	0	0	0	0	0.0000	0.0000	0.0000	0.0000
AUTO	0	0	0	0	0.0000	0.0000	0.0000	0.0000
MINI	1	0	2000	0	0.0000	0.0000	0.0000	0.0000
MCMM	50	0	0	0	0.0000	0.0000	0.0000	0.0000
AUTO	0	0	0	0	0.0000	1.0000	0.0000	0.0000
BGIN	0	0	0	0	0.0000	0.0000	0.0000	0.0000
MCOP	1	0	10000	0	0.0000	2.0000	0.0000	0.0000
READ	0	0	0	0	0.0000	0.0000	0.0000	0.0000
MINI	1	0	2000	0	0.0000	0.0000	0.0000	0.0000
END	0	0	0	0	0.0000	0.0000	0.0000	0.0000
MBAE	-1	0	0	0	0.0000	0.0000	0.0000	0.0000

FFLD: Use OPLS_2001 with constant dielectric electrostatics, which is appropriate if GB/SA solvation is being used.

SOLV: The GB/SA effective water model is being used.

EXNB: Extended non-bonded cut-offs should be used with GB/SA solvation.

BDCO: Use the Bond Dipole CutOff (BDCO) method for truncating electrostatic interactions. Arg5 and arg6 specify the cutoffs used for charge-dipole and charge-charge interactions, respectively.

MSYM: Invokes the numbering symmetry library mmsym, which automatically and more generally identifies a suitable numbering order for using in comparing conformations.

MBAE: Turn on eMBrAcE

arg1 = 1: Using Energy Difference mode.

arg2 = 1: Conformation searches will be performed in eMBrAcE calculations.

arg3 = 0: The lowest energy structures of the ligand-receptor complex are written to the output structure file.

MCMM: Use Monte Carlo Multiple Minimum searching. Arg1 defines the number of steps to use for the search of the receptor.

MCNV: Sets the number of degrees of freedom to be varied at each MC step. With different arg1 and arg2 values, the search varies a random number of degrees of freedom between the numbers defined in arg1 and arg2. For eMBrAcE calculations, we recommend setting arg1 to a small number and arg2 to no more than 10.

MCSS: MC structure selection allows for the setting of selection of starting structure for the search steps. Arg1=2 defines use-directed selection of starting structures, where the least used structures will be used as starting geometries as long as they are low enough in energy (as

defined in arg5). This is more efficient in exploring new areas of the potential energy surface than, for instance, a random-walk starting geometry scheme. Arg5 gives the energy window for selecting a new starting structure, which must be within arg5 kJ/mol of the lowest energy conformer found in the search.

MCOP: Monte Carlo options that determine what and how often data is written to the log file and how many structures are saved.

arg1 = 1: Print information for every search step to the log file.

arg3 = 10000: Use a large number here to avoid spurious lines in the eMBrAcE summary table at the end of the run.

arg6 = 1.0: Save only the lowest energy conformation of the receptor (this value must be specified as 1.0).

DEMX: arg5 = 500.0 means retain only the conformations within 500 kJ/mol of the minimum energy conformation found so far.

SUBS: Substructures may be used. All opcodes related to the substructure (SUBS, FXAT, FXDI, FXBA, and FXTA) should be defined in an .sbc file. The .sbc file should contain only lines referring to the receptor with the receptor atoms numbered starting from 1.

READ: The first READ command obtains the structure of the receptor.

AUTO: The first AUTO automatically sets up the MCMM conformational search of the receptor. Here arg6 must be absent or 0.0 since the search of the receptor is not regarded as part of a serial search.

MINI: Minimize the energy of each conformation generated in the search of the receptor using the PRCG minimization technique, arg1=1, for up to arg3 (2000) iterations.

MCMM: The second MCMM causes Monte Carlo Multiple Minimum searching to be used in the search of each isolated ligand and each ligand-receptor complex. Arg1 defines the number of steps to use for each search. This need not match that used in the search of the receptor.

AUTO: The second AUTO sets up the MCMM search for each isolated ligand and each ligand-receptor complex. Here arg6 should be set to 1 because the searches of the isolated ligands and the complexes are regarded as part of a serial search.

BGIN: Start the loop which processes each ligand in turn.

MCOP: Monte Carlo options that determine what and how often data is written to the log file and how many structures are saved for the searches of the ligand-receptor complexes.

arg1 = 1: Print information for every search step to the log file.

arg3 = 10000: Use a large number here to avoid spurious lines in the eMBrAcE summary table at the end of the run.

arg6 = 2.0: Save only the two lowest energy conformations from the complexes for each ligand.

READ: Sequentially read the ligand structures.

MINI: Minimize the energy of each conformation generated in the search of each ligand and each ligand-receptor complex using the PRCG minimization technique, arg1=1, for up to arg3 (2000) iterations.

END: End the loop for the ligands.

MBAE: Turn off eMBrAcE and summarize the results. arg1 = -1

16.4.2.2 Low-Mode Conformational Search

This example performs a low-mode conformational search. The input file is similar to that of the MCMM search example above. A copy of this file is available at:

\$SCHRODINGER/macromodel-vversion/samples/Examples/MBAE_LMCS.com

The example is followed by descriptions of the opcodes. Opcode descriptions that are omitted may be read from the corresponding opcodes in the MCMM conformational search example given above, in [Section 16.4.2.1 on page 165](#).

SEARCH_LMCS.mae

SEARCH_LMCS-out.mae

FFLD	11	1	0	0	1.0000	0.0000	0.0000	0.0000
SOLV	3	1	0	0	0.0000	0.0000	0.0000	0.0000
EXNB	0	0	0	0	0.0000	0.0000	0.0000	0.0000
BDCO	0	0	0	0	89.4427	99999.0000	0.0000	0.0000
MSYM	0	0	0	0	0.0000	0.0000	0.0000	0.0000
MBAE	1	1	0	0	0.0000	0.0000	0.0000	0.0000
LMCS	50	0	0	0	0.0000	0.0000	3.0000	6.0000
MCSS	2	0	0	0	500.0000	0.0000	0.0000	0.0000
MCOP	1	0	1000	0	0.0000	1.0000	0.0000	0.0000
DEMX	0	0	0	0	500.0000	0.0000	0.0000	0.0000
SUBS	0	0	0	0	0.0000	0.0000	0.0000	0.0000
READ	0	0	0	0	0.0000	0.0000	0.0000	0.0000
AUTO	0	0	0	0	-1.0000	0.0000	0.0000	0.0000
CONV	2	0	0	0	0.0500	0.0000	0.0000	0.0000
MINI	1	0	2000	0	0.0000	0.0000	0.0000	0.0000
LMCS	50	0	0	0	0.0000	0.0000	0.0000	0.0000
BGIN	0	0	0	0	0.0000	0.0000	0.0000	0.0000
MCOP	1	0	10000	1	0.0000	2.0000	0.0000	0.0000
READ	0	0	0	0	0.0000	0.0000	0.0000	0.0000
MINI	1	0	2000	0	0.0000	0.0000	0.0000	0.0000

END	0	0	0	0	0.0000	0.0000	0.0000	0.0000
MBAE	-1	0	0	0	0.0000	0.0000	0.0000	0.0000

LMCS: Use the low-mode conformational search method. Arg1 defines the number of steps to use for the search of the receptor.

AUTO: Set up comparison atom lists, torsional constraints, and chirality checks for all the LMCS searches in this calculation. Arg5 should be set to -1 so that torsional moves (TORS) are not identified, and arg6 should be 0.0 or absent as the serial aspects of the searches are controlled by the next MCOP.

MINI: Minimize the energy of each conformation generated in the search of the receptor using the PRCG minimization technique, arg1=1, for up to arg3 (2000) iterations.

LMCS: The second LMCS causes low-mode conformational searching to be used in the search of each isolated ligand and each ligand-receptor complex. Arg1 defines the number of steps to use for each search. This need not match that used in the search of the receptor.

MCOP: Monte Carlo options that determine what and how often data is written to the log file and how many structures are saved for the searches of the ligand-receptor complexes.

arg1 = 1: Print information for every search step to the log file.

arg3 = 10000: Use a large number to avoid spurious lines in the eMBrAcE summary table at the end of the run.

arg4 = 1: Indicates that the searches are considered to be part of a serial calculation.

arg6 = 2.0: Saves only the two lowest energy conformations from the complexes for each ligand.

16.4.2.3 MCMM Conformation Search With COPY/ALGN

eMBrAcE conformational searches may also be used with COPY and ALGN to position ligands using a crystal structure of a complex of a closely related ligand. While this combination of commands may be useful, the positioning is crude and the searching of conformational space is slow and quite limited compared to that available in Schrödinger's docking program, Glide.

Below is an example command file for COPY and ALGN in combination with MBAE for an MCMM conformational search. Apart from the alignment step, the input file is similar to that of the MCMM search example above. A copy of the input file is available at:

`$SCHRODINGER/macromodel-vversion/samples/Examples/MBAE_ALGN.com`

The example is followed by descriptions of the opcodes. Opcode descriptions that are omitted may be read from the corresponding opcodes in the MCMM conformational search example given above, in [Section 16.4.2.1 on page 165](#).

```

ALGN_SEARCH_MCMC.mae
ALGN_SEARCH_MCMC-out.mae
FFLD      11      1      0      0      1.0000      0.0000      0.0000      0.0000
SOLV      3      1      0      0      0.0000      0.0000      0.0000      0.0000
EXNB      0      0      0      0      0.0000      0.0000      0.0000      0.0000
BDCC      0      0      0      0      89.4427 99999.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
WRIT      0      0      0      0      0.0000      0.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
COPY      0      0      0      0      0.0000      0.0000      0.0000      0.0000
BGIN      0      0      0      0      0.0000      0.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
ALGN      3      1      5      0      0.0000      0.0000      0.0000      0.0000
END        0      0      0      0      0.0000      0.0000      0.0000      0.0000
RWND      0      0      0      0      0.0000      0.0000      0.0000      0.0000
MSYM      0      0      0      0      0.0000      0.0000      0.0000      0.0000
MBAE      1      1      0      0      0.0000      0.0000      0.0000      0.0000
MCMC      50      0      0      0      0.0000      0.0000      0.0000      0.0000
MCNV      1      4      0      0      0.0000      0.0000      0.0000      0.0000
MCSS      2      0      0      0      500.0000      0.0000      0.0000      0.0000
MCOP      1      0     1000      0      0.0000      1.0000      0.0000      0.0000
DEMX      0      0      0      0     1000.0000      0.0000      0.0000      0.0000
SUBS      0      1      0      0      0.0000      0.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
AUTO      0      0      0      0      0.0000      0.0000      0.0000      0.0000
MINI      1      0     2000      0      0.0000      0.0000      0.0000      0.0000
MCMC      50      0      0      0      0.0000      0.0000      0.0000      0.0000
AUTO      0      0      0      0      0.0000      1.0000      0.0000      0.0000
BGIN      0      0      0      0      0.0000      0.0000      0.0000      0.0000
MCOP      1      0     1000      0      0.0000      2.0000      4.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
MINI      1      0     2000      0      0.0000      0.0000      0.0000      0.0000
END        0      0      0      0      0.0000      0.0000      0.0000      0.0000
MBAE     -1      0      0      0      0.0000      0.0000      0.0000      0.0000

```

READ: The first READ reads in the receptor.

WRIT: Record the receptor structure in the output structure file.

READ: The second READ reads in the first (reference) ligand.

COPY: Copy the current structure to the reference storage area for use in aligning the ligands.

BGIN: Loop over all remaining ligands in the input structure file.

READ: Read in the next ligand structure.

ALGN: Align the current structure with the reference structure by center of mass and moments of inertia (arg1 = 3), weighting atom positions by mass (arg2=1), and write all 4 such alignments to the output structure file (arg3 = 5).

END: End of loop for aligning the ligands.

RWND: Rewind the output structure file and use it as the input structure file for the remainder of the calculation.

DEMX: $\text{arg5} = 1000.0$ means retain only the conformations within 1000 kJ/mol of the minimum energy conformation found so far. We recommend using a large value to permit seed structures with different alignments to be restrained and sampled.

SUBS: Substructures may be used. All opcodes related to the substructure (SUBS, FXAT, FXDI, FXBA, and FXTA) should be defined in an .sbc file. The .sbc file should contain only lines referring to the receptor with the receptor atoms numbered starting from 1. In this case, we specify $\text{arg2} = 1$ to indicate that the name of the original *jobname*.sbc file is to be used.

READ: This READ command obtains the structure of the receptor.

MCMM: The second MCMM selects Monte Carlo Multiple Minimum searching for the search of each isolated ligand and each ligand-receptor complex. Arg1 defines the number of steps to use for each search. This need not match that used in the search of the receptor.

MCOP: Monte Carlo options that determine what and how often data is written to the log file and how many structures are saved for the searches of the ligand-receptor complexes.

$\text{arg1} = 1$: Print information for every search step to the log file.

$\text{arg3} = 10000$: Use a large number to avoid spurious lines in the eMBrAcE summary table at the end of the run.

$\text{arg6} = 2.0$: Save only the two lowest energy conformations from the complexes for each ligand.

$\text{arg7} = 4.0$: Seed the search for each ligand using all four alignments previously generated using the ALGN command.

16.4.3 Distributed eMBrAcE Calculations

Any eMBrAcE job may be distributed across a number of processors. To distribute an eMBrAcE job, add an NPRC line early in the .com file, such as immediately above the FFLD line. See the description of the NPRC opcode in the *MacroModel Reference Manual* for more information. For example,

```
NPRC      2      16      0      0      0.0000      0.0000      0.0000      0.0000
```

would run the job on two processors with 16 ligands processed in each subjob.

For energy difference eMBrAcE calculations, the receptor is processed first within the parent process prior to starting up any child processes.

There are two restrictions on distributed MBAE jobs. First, the output structure mode, which is controlled by `arg3` of the MBAE opcode, must be either 0 (complexes only) or 1 (ligands only) and not 2 (receptor, ligands and complexes). Second, distributed eMBrAcE runs may not contain COPY and ALGN commands to pre-position the ligands. Instead the pre-positioning of the ligands must be accomplished in a separate non-distributed calculation. The output from that calculation is then used as input for the distributed eMBrAcE conformational search calculation.

An example command file to align ligands prior to an eMBrAcE conformational search calculation is given below. This file is available at

`$SCHRODINGER/macromodel-vversion/samples/Examples/ALGN_pv.com`

COPY_ALGN.mae

COPY_ALGN-out.mae

READ	0	0	0	0	0.0000	0.0000	0.0000	0.0000
WRIT	0	0	0	0	0.0000	0.0000	0.0000	0.0000
READ	0	0	0	0	0.0000	0.0000	0.0000	0.0000
COPY	0	0	0	0	0.0000	0.0000	0.0000	0.0000
BGIN	0	0	0	0	0.0000	0.0000	0.0000	0.0000
READ	0	0	0	0	0.0000	0.0000	0.0000	0.0000
ALGN	3	1	5	0	0.0000	0.0000	0.0000	0.0000
END	0	0	0	0	0.0000	0.0000	0.0000	0.0000

READ: The first READ reads in the receptor.

WRIT: Record the receptor structure in the output structure file.

READ: The second READ reads in the first (reference) ligand.

COPY: Copy the current structure to the reference storage area for use in aligning the ligands.

BGIN: Loop over all remaining ligands in the input structure file.

READ: Read in the next ligand structure.

ALGN: Align the current structure with the reference structure by center of mass and moments of inertia (`arg1 = 3`), weighting atom positions by mass (`arg2=1`), and write all 4 such alignments to the output structure file (`arg3 = 5`).

END: End of loop for aligning the ligands.

An example command file is given below for a distributed eMBrAcE conformational search that uses the output of the separate COPY/ALGN job above. The file is available at:

`$SCHRODINGER/macromodel-vversion/samples/Examples/
dist_ALGNed_MBAE_MCOMM.com`

The descriptions of most of the opcodes may be read from the corresponding opcodes in the MCMM conformational search example given above, in [Section 16.4.2.1 on page 165](#).

COPY_ALGN-out.mae

MBAE_CSEARCH-out.mae

NPRC	2	16	0	0	0.0000	0.0000	0.0000	0.0000
FFLD	11	1	0	0	1.0000	0.0000	0.0000	0.0000
SOLV	3	1	0	0	0.0000	0.0000	0.0000	0.0000
EXNB	0	0	0	0	0.0000	0.0000	0.0000	0.0000
BDCO	0	0	0	0	89.4427	99999.0000	0.0000	0.0000
MSYM	0	0	0	0	0.0000	0.0000	0.0000	0.0000
MBAE	1	1	1	0	0.0000	0.0000	0.0000	0.0000
MCMM	50	0	0	0	0.0000	0.0000	0.0000	0.0000
MCNV	1	4	0	0	0.0000	0.0000	0.0000	0.0000
MCSS	2	0	0	0	500.0000	0.0000	0.0000	0.0000
MCOP	1	0	10000	0	0.0000	1.0000	0.0000	0.0000
DEMX	0	0	0	0	1000.0000	0.0000	0.0000	0.0000
SUBS	0	0	0	0	0.0000	0.0000	0.0000	0.0000
READ	0	0	0	0	0.0000	0.0000	0.0000	0.0000
AUTO	0	0	0	0	0.0000	0.0000	0.0000	0.0000
MINI	1	0	10	0	0.0000	0.0000	0.0000	0.0000
AUTO	0	0	0	0	0.0000	1.0000	0.0000	0.0000
MCMM	50	0	0	0	0.0000	0.0000	0.0000	0.0000
BGIN	0	0	0	0	0.0000	0.0000	0.0000	0.0000
MCOP	1	0	10000	0	0.0000	2.0000	4.0000	0.0000
READ	0	0	0	0	0.0000	0.0000	0.0000	0.0000
MINI	1	0	10	0	0.0000	0.0000	0.0000	0.0000
END	0	0	0	0	0.0000	0.0000	0.0000	0.0000
MBAE	-1	0	0	0	0.0000	0.0000	0.0000	0.0000

NPRC: Distribute the calculation. Arg1 is the number of processors to use and arg2 is the number of input structures to process in each subjob. Arg2 should be a multiple of 4 for ligands prepared with ALGN so that all four alignments of each ligand are sent to the same child process.

Free Energy Simulations

This chapter provides some of the theory behind the powerful simulation method of free energy perturbation calculations, and describes how to perform a free energy perturbation with MacroModel. If you intend to perform free energy perturbations, we strongly recommend that you become familiar with the current literature on this topic (see [Section 17.5 on page 186](#)).

Note: The MM2* and MM3* force fields should not be used for Free Energy Perturbation.

17.1 Free Energy Perturbation

The fundamental expression for free energy calculations as implemented in MacroModel is:

$$G_B - G_A = \Delta G = -RT \ln \left\langle \exp \left[\frac{-(H_B - H_A)}{RT} \right] \right\rangle_A \quad (1)$$

where H_A and H_B are the Hamiltonians for the two systems, A and B, and the $\langle \rangle$ notation represents an ensemble average over system A. This expression is valid only when there is a very small difference between the two systems, A and B (i.e., $H_B - H_A \approx RT$). To perform free energy calculations on meaningful systems, you generally perform a series of smaller simulations (windows), which can be summed to obtain a total free energy difference.

The coupling parameter λ is used to define each window in terms of the two endpoints, A and B, so that at any stage the Hamiltonian over which the ensemble average will be generated for the system is described in terms of the Hamiltonians for the endpoints:

$$H(\lambda) = \lambda H_B + (1 - \lambda) H_A \quad (2)$$

Then the overall expression for performing a simulation becomes:

$$\Delta G = \sum_{i=0}^n -RT \ln \left\langle \exp \left[\frac{-(H_{\lambda_{i+1}} - H_{\lambda_i})}{RT} \right] \right\rangle_{\lambda_i} \quad (3)$$

where n is the number of windows to be used in the simulation. At any window the “perturbation” is between a state at λ and that at $\lambda + d\lambda$. The method for performing this simulation as implemented in MacroModel is described as “single topology.” That is, at any value of λ , a set

of interactions is generated by mixing the interactions of the endpoints. For example, if an atom has a charge x in the starting point of the simulation and a charge y in the end point, at any given value of λ (i.e., at any window), the charge on the atom will be $z = y\lambda + (1-\lambda)x$. This process will be repeated for all the interactions in the system and the simulation will be performed with this new set of interactions. At points during the simulation, the energy, $H(\lambda)$, will be evaluated. Then the interactions corresponding to the “next” value of λ ($\lambda+d\lambda$) will be generated and $H(\lambda+d\lambda)$ will be calculated. These terms are used in the exponential of expression (3). In MacroModel we usually take advantage of the fact that for most values of λ , the energy at $\lambda-d\lambda$ can also be evaluated—this procedure is known as “double wide sampling,” and effectively allows a “forward” and “reverse” simulation to be performed in one simulation.

17.2 Setting Up FEP Calculations

Currently free energy calculations can be set up only by direct manipulation of the MacroModel command files. There is no Maestro interface for this procedure. The individual MacroModel opcodes are described in detail in [Chapter 3](#) of the *MacroModel Reference Manual*. This section illustrates how to combine these opcodes and produce a sensible command file. To do this we will use a simple example—a calculation of the free energy difference between D and L forms of the “alanine dipeptide”—this should give a free energy difference of zero in isolation. Though this is not a very useful simulation, it illustrates free energy perturbation calculations. Note that most of the perturbations we have tested involve small perturbations with changes of up to only three atoms.

17.2.1 The Structure File

When you perform a free energy perturbation calculation in MacroModel, the input structure file (.mae) must contain two separate MacroModel structures. These correspond to the start and end point for the calculation. The numbering of the structures is critical—all “equivalent” atoms must have the same number in each structure. Dummy atoms (MacroModel atom type 61) must be used as place holders for atoms created or destroyed during the perturbation. The numbering of the structures used in our example is shown in [Figure 17.1](#).

Here we are performing a FEP calculation between the L and D enantiomers of the alanine dipeptide using the AMBER* force field, so we have no explicit hydrogen on the alpha carbon. The perturbation involves changing dummy atom 13 in the starting (L) structure to a united atom methyl in the (D) structure. At the same time, atom number 9, which begins as a united atom methyl in the starting structure, is extinguished to a dummy atom in the structure at the end of the simulation. Numbering of the structures must be checked very carefully before a free energy simulation is attempted.




Figure 17.1. Numbering system for the alanine dipeptide.

17.2.2 The Command File

The command file used to perform the perturbation of the L to the D form is shown below.

```
fep.mae
fep-out.mae
EXNB      0      0      0      0      0.0000      0.0000      0.0000      0.0000
BDCO      0      0      0      0      0.0000      0.0000      0.0000      0.0000
FFLD      3      0      0      0      0.0000      0.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
FEIA      0      0      0      0      0.0000      0.0000      0.0000      0.0000
MCNV      1      2      0      0      0.0000      0.0000      0.0000      0.0000
MCSD      1      2      0      0      0.0000      0.0000      300.0000      0.0000
TORS      5      6      0      0      30.0000      180.0000      0.0000      0.0000
TORS      6      7      0      0      30.0000      180.0000      0.0000      0.0000
BGIN      0      0      0      0      0.0000      0.0000      0.0000      0.0000
FEAV      0      0      0      0      0.0000      0.0000      0.0000      0.0000
MINI      1      0      1000      0      0.0000      0.0000      0.0000      0.0000
MDIT      0      0      0      0      300.0000      0.0000      0.0000      0.0000
MDYN      0      0      1      0      1.5000      10.0000      300.0000      0.0000
FESA      0      0      0      0      0.0000      0.0500      0.0000      0.0000
MDYN      0      0      1      0      1.5000      100.0000      300.0000      0.0000
END      0      0      0      0      0.0000      0.0000      0.0000      0.0000
FESU      0      0      0      0      0.0000      0.0000      0.0000      0.0000
```

EXNB: Extended non-bonded cutoffs should be used in free energy calculations.

FFLD: This example used the AMBER* force field with the default distance-dependent dielectric treatment for electrostatics. We also recommend the use of the GB/SA solvation model in a real application.

READ: Required to read the first structure in the .mae file.

FEIA: Reads the second structure (the perturbation endpoint) from the .mae file and sets up the arrays for averaging of the interactions.

MCNV: Specifies that between 1 and 2 of the torsions specified by TORS commands will be varied at each MC step.

MCS D: In this example we use mixed-mode Monte Carlo/Stochastic Dynamics to perform the simulation. Because a successful FEP simulation requires a thorough sampling of the conformation space at each window, we strongly recommend using the mixed-mode simulation method during free energy calculations.

TORS: Defines the torsion angles used in the Monte Carlo part of the MCS D calculation. In this case it is the ϕ and ψ angles of the dipeptide.

BGIN: Free energy perturbation calculations are usually performed in a BGIN/END loop, which completes all the windows from $\lambda=0$ to $\lambda=1$.

FEAV: Generates the interactions corresponding to the value of λ given in arg5. The interpretation of this command can lead to some confusion. The first time it is encountered, it will generate an interaction array with the value of λ as given in arg5. When placed inside a BGIN/END pair, as in this example, the behavior is slightly different—when subsequently encountered, the FEAV command will increment the current value of λ (i.e., the middle of the double wide sampling window) by the difference between arg5 and arg6 in the following FESA command. In this case it will be incremented to 0.05 the first time through the BGIN/END loop.

MINI: Once the interactions are generated, it is necessary to minimize the structure with this set of interactions in order to remove any excess potential energy before beginning the dynamics part of the simulation.

MDIT: Sets initial random velocities corresponding to an initial temperature of the value of arg5 (in this case, 300 K).

MDYN: At any given value of λ (i.e., at the start of every window) it is good practice to first perform a short equilibration simulation before collecting data for the free energy calculation.

FESA: Begins the collection of free energy samples for the current value of (λ)—the middle of the window. The values of arg5 and arg6 indicate the values of the left side of the window and the right side of the window for double wide sampling. The difference between the values of arg5 and arg6 controls the number of windows which will be performed in the simulation—the relationship is $\text{\#windows} = 1/(\text{arg6} - \text{arg5}) + 1$. A maximum of 101 windows can be used in a MacroModel FEP calculation so the minimum permissible value of the difference between arg5 and arg6 is 0.01. When used inside a BGIN/END pair, as in this example, arg5 and arg6 will be incremented each pass through the loop, as illustrated in the table below.

Most simulations require at least 20 windows to ensure that the free energy difference per window (which must be <3 RT) remains sufficiently small.

Pass	$\lambda(\text{left})$	$\lambda(\text{middle})$	$\lambda(\text{right})$	
1	0.00	0.00	0.05	(no middle->left sampling)
2	0.0	0.05	0.10	(middle->left and left->right)
3	0.05	0.10	0.15	(middle->left and left->right)
4	0.10	0.15	0.20	(middle->left and left->right)
20	0.90	0.95	1.00	(middle->left and left->right)
21	0.95	1.0	1.0	(no middle->right sampling)

MDYN: Describes the dynamics simulation in which the free energy sampling is done. In this simple case we are only doing 100 ps of mixed-mode simulation at each window. In real cases > 500 ps simulation per window will probably be required.

END: The end of the BGIN/END loop. The loop terminates when the value of λ (middle) exceeds 1.0.

FESU: Prepare a summary of the free energy change during the simulation.

17.2.3 Other Possible Command Files

Performing specific windows:

Although it is useful to wrap free energy calculations in a BGIN/END loop and simulate with all windows, it is possible to run any specific windows. The windows are independent and free energy differences obtained from separate runs can be usefully combined. It is also sometimes useful to re-run specific windows from a simulation to obtain better statistics. The command file below performs double wide sampling for two windows with $\lambda(\text{middle})=0.5$ and $\lambda(\text{middle})=0.55$:

```
fep-diala.mae
fep-diala-out.mae
EXNB      0      0      0      0      0.0000      0.0000      0.0000      0.0000
BDCE      0      0      0      0      0.0000      0.0000      0.0000      0.0000
FFLD      3      0      0      0      0.0000      0.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
FEIA      0      0      0      0      0.0000      0.0000      0.0000      0.0000
MCNV      1      2      0      0      0.0000      0.0000      0.0000      0.0000
MCSD      1      2      0      0      0.0000      0.0000      300.0000      0.0000
TORS      5      6      0      0      30.0000      180.0000      0.0000      0.0000
TORS      6      7      0      0      30.0000      180.0000      0.0000      0.0000
FEAV      0      0      0      0      0.5000      0.0000      0.0000      0.0000
MINI      1      0      1000      0      0.0000      0.0000      0.0000      0.0000
```

MDIT	0	0	0	0	300.0000	0.0000	0.0000	0.0000
MDYN	0	0	1	0	1.5000	10.0000	300.0000	0.0000
FESA	0	0	0	0	0.4500	0.5500	0.0000	0.0000
MDYN	0	0	1	0	1.5000	100.0000	300.0000	0.0000
FEAV	0	0	0	0	0.5500	0.0000	0.0000	0.0000
MINI	1	0	1000	0	0.0000	0.0000	0.0000	0.0000
MDIT	0	0	0	0	300.0000	0.0000	0.0000	0.0000
MDYN	0	0	1	0	1.5000	10.0000	300.0000	0.0000
FESA	0	0	0	0	0.5000	0.6000	0.0000	0.0000
MDYN	0	0	1	0	1.5000	100.0000	300.0000	0.0000
FESU	0	0	0	0	0.0000	0.0000	0.0000	0.0000

Using the distributed MacroModel procedure to perform a free energy calculation:

For more information about distributed MacroModel, see [Section 4.3 on page 50](#). The NPRC command can be used to distribute the free energy calculation over a number of hosts in a manner similar to that described for Monte Carlo conformational searching. In this case the smallest job that can be distributed on any host is one complete FEP window, so the opportunities for load balancing are limited.

17.2.4 The Output File

Similar output for each window of the perturbation will be obtained in the .log file. The output obtained when the command file described above was run appears below. This output is for the value of $\lambda(\text{middle})=0.4$

```
FEP averaged force field, lambda = 0.400000
Loading FEP coordinates 1
Starting conjugate gradient minimization.
Minimization converged; gradient = 0.335E-01 .LT. 0.500E-01
Iterations = 120 out of 1000
Conf 1 E = -138.502 ( 0.034) kJ/mol

BatchMin V8.5      Stochastic Dynamics Simulation
                  MC/SD Mixed Mode

Time      Total E      0.5ps T      Ave H      Ave T      Ave H( 300)
(ps)      (kJ/mol)      (deg K)      (kJ/mol)      (deg K)      (kJ/mol)
1.000      -60.8      297.1      -97.2      295.4      -96.5
2.001      -72.2      172.3      -103.0      260.1      -96.5
3.000      -45.2      263.5      -102.1      265.5      -96.5
4.001      -62.5      308.8      -101.7      272.9      -97.3
5.001      -70.4      252.1      -101.4      269.1      -96.4
6.000      -73.9      272.6      -102.2      267.2      -96.9
7.001      -51.8      276.7      -101.6      265.1      -96.0
8.001      -60.0      271.4      -101.2      266.0      -95.7
9.000      -55.4      266.2      -101.5      268.3      -96.4
10.000     -54.3      261.5      -101.0      268.6      -95.9

Final potential and kinetic energy =      -93.46      39.77 kJ/mol
CPU Time =      1.2 sec

Average potential energy =      -101.02 kJ/mol
```

```

Acceptance Ratio for Mixed Mode Simulation = 0.02744863
Average kinetic energy = 43.55 kJ/mol (Av temperature = 268.6 deg K)
Average total energy = -57.47 kJ/mol (Std dev = 11.13 kJ/mol)

```

```

Average potential energy <H> scaled to 300.0 deg K = -95.92 kJ/mol
      Av stretch      12.72 kJ/mol
      Av bend         15.59 kJ/mol
      Av torsion       25.90 kJ/mol
      Av van der Waals 12.57 kJ/mol
      Av electrostatic -162.70 kJ/mol
      Av solvation 1   0.00 kJ/mol
      Av solvation 2   0.00 kJ/mol

```

```

BatchMin V8.5      Stochastic Dynamics Simulation
                   MC/SD Mixed Mode

```

```

Free energy perturbation calculation
Averaged interaction array, lambda = 0.400000

```

```
FEP Window 9:
```

Time (ps)	Total E (kJ/mol)	0.5ps T (deg K)	Average T (deg K)	<G>0.35-0.40 (kJ/mol)	<G>0.45-0.40 (kJ/mol)
10.000	-35.5	306.2	284.0	-1.11	0.89
20.001	-46.5	294.7	304.7	-1.16	0.98
30.000	-19.9	292.4	302.6	-1.04	0.82
40.000	-24.2	366.7	301.7	-0.98	0.77
50.001	-57.1	282.4	300.0	-0.95	0.72
60.000	-24.3	317.0	299.3	-0.94	0.70
70.001	-24.6	265.0	302.7	-0.96	0.72
80.001	-31.0	310.1	301.0	-0.94	0.71
90.002	-31.1	287.6	300.3	-0.95	0.71
100.001	-59.1	277.9	298.9	-0.97	0.74

```

Final potential and kinetic energy = -84.41      38.78 kJ/mol
CPU Time = 12.5 sec

```

```
Average potential energy = -96.48 kJ/mol
```

```

Acceptance Ratio for Mixed Mode Simulation = 0.03341983
Average kinetic energy = 48.46 kJ/mol (Av temperature = 298.9 deg K)
Average total energy = -48.02 kJ/mol (Std dev = 15.55 kJ/mol)

```

```

Average potential energy <H> scaled to 300.0 deg K = -96.29 kJ/mol
      Av stretch      15.40 kJ/mol
      Av bend         17.81 kJ/mol
      Av torsion       27.36 kJ/mol
      Av van der Waals 13.13 kJ/mol
      Av electrostatic -169.99 kJ/mol
      Av solvation 1   0.00 kJ/mol
      Av solvation 2   0.00 kJ/mol

```

```
Free energy perturbation window 9
```

```

Lambda (left,center,right) = 0.350 0.400 0.450
Number of samples in average = 3921

```

	Bonded	Nonbonded	Solvation	Total	+/-
G(left-center) =	-0.094	-0.867	0.000	-0.967	0.134 kJ/mol
G(right-center) =	0.059	0.694	0.000	0.737	0.166 kJ/mol

Similar output was obtained for the other 20 windows of this simulation. The first part of the output shows the result of the minimization at this value of λ and the equilibration simulation. Then the actual free energy sampling is performed and a free energy for the “forward” and

“reverse” simulation is reported. Notice that free energy reported should become stable to a few hundredths of a kJ/mol near the end of the sampling period. Note, however, that this is a necessary but insufficient indicator of convergence during the simulation, as this may only indicate that the simulation is stuck in one of a number of accessible potential energy minima.

After the actual sampling phase is over, a summary of the free energy change for this window will be printed as shown above. *Note that the “component” (e.g., Bonded, Nonbonded, and Solvation) free energies are only approximate and do not necessarily sum to the “Total” value.* After all windows are completed, a summary is printed of the perturbation over all windows (with the “forward” and “reverse” simulations shown side by side). At the end of the file is a summary of the complete perturbation:

```
Summary of total perturbation. Lambda: 0.000 -> 1.000
      Forward      Reverse      Average
G(bonded)      =      0.10200      0.07908      0.09054 kJ/mol
G(nonbonded)   =      1.08893     -1.12791     -0.01949 kJ/mol
G(solvation)   =      0.00000      0.00000      0.00000 kJ/mol
G(total)       =      1.03701     -0.90759      0.06471 kJ/mol
Standard dev   =      0.31498      0.29871           kJ/mol
```

This is probably the most important part of the output. The average of the “forward” and “reverse” total free energies gives an estimate of the free energy change for the entire simulation. In this case it is close to zero, as expected (sampling for 500 ps per window reduces the value to 0.02 kJ/mol). The standard deviation of the free energy is useful as an estimate of the error. In this case the free energy would be reported as 0.1 0.3kJ/mol. At 500 ps sampling it is 0.02±0.14 kJ/mol. Note that we use the larger of the “forward” or “reverse” standard deviation values and not the average. As mentioned above, only the “total” free energies are strictly meaningful—the *components are only estimates*. In addition to examining the final free energy summary, it is useful to plot the progress of the perturbation.

There should be a continuous, smooth change in free energy from window to window with no large “jumps” or sudden increases in the standard deviation. If discontinuities are observed, then the simulation should probably be run with more sampling per window and possibly more windows.

17.3 Criteria for a Successful Simulation

The answers to the following questions will help you perform successful simulations.

1. Do the input structures have the same numbering with dummy atoms used for atoms which are present in one structure and not in another?

One useful way to check the input file is to set up a simple command file as follows:

```
fep-diala.mae
fep-diala-out.mae
  EXNB      0      0      0      0      0.0000      0.0000      0.0000      0.0000
  BDCO      0      0      0      0      0.0000      0.0000      0.0000      0.0000
  FFLD      3      0      0      0      0.0000      0.0000      0.0000      0.0000
  READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
  FEIA      0      0      0      0      0.0000      0.0000      0.0000      0.0000
  FEAV      0      0      0      0      0.0000      0.0000      0.0000      0.0000
  MINI      1      0    1000      0      0.0000      0.0000      0.0000      0.0000
  FEAV      0      0      0      0      1.0000      0.0000      0.0000      0.0000
  MINI      1      0    1000      0      0.0000      0.0000      0.0000      0.0000
```

This command file performs a minimization with $\lambda = 0$ and $\lambda = 1$. Record the final minimized energies in each case. Now swap the structures in the input (.mae) file and re-run the calculation. For the second run, with the order of the structures reversed, the energy at the $\lambda = 1$ point should be the same as the energy from the $\lambda = 0$ point of the first run. Similarly, the energy of the $\lambda = 0$ point in the second run should be the same as the $\lambda = 1$ energy of the first. If these do not match exactly, then there is likely to be something wrong with the numbering of the input structure files.

2. Is sufficient sampling being performed?

While the expression (1) described above is exact, it is implicit that at each window complete sampling of all important conformations is achieved. In practice this is very difficult to achieve and we strongly recommend the use of the mixed-mode simulation methodology to help solve the sampling problem. Even so, you should also perform independent tests for convergence with simulations at least as long as the sampling period in the free energy calculations. Ideally you should be satisfied that you can achieve a reasonable degree of convergence for simulations of the starting structures in the length of time used for the sampling in the free energy calculations. Good tests for convergence are starting from different starting geometries and (for non-chiral molecules) monitoring the populations of “equivalent” conformations (e.g., +gauche/-gauche torsion angles). Another possibility is to compare structures sampled from a simulation with those obtained from an extensive conformational search. Most importantly note that stable values for the average free energy change over several tens of picoseconds of sampling, while required for a converged simulation, are not alone proof of convergence. If the simulation is “stuck” in one particular minimum, then the average quantities will look to be converged but the actual free energy may be incorrect.

3. Is the force-field adequate?

Free energy perturbation assumes that the potential energy of the system is well described by force field parameters. You cannot expect to get quantitative agreement with experimental free energies (or make useful predictions) unless you are using a high quality force field. You

should carefully check the quality of all the parameters in use in a free energy calculations. A summary of the parameter quality is printed at the top of the `.log` file.

4. Are the dummy atom parameters appropriate?

Natural bond lengths to dummy atoms may need to be set to approximate the length of the real counterpart of the dummy atom, e.g., the atom into which the dummy “grows.” In the past we have used very short (ca. 0.5 Å) natural bond lengths for dummy atoms. However, this can cause problems when the dummy atom has significant non-bonded interaction with other parts of the molecule early in the simulation. For example, in the case of perturbing axial methyl cyclohexane to equatorial methyl cyclohexane using the united atom model for AMBER*, if the dummy atoms start at 0.5 Å the free energy difference between the two conformations is greatly overestimated. Setting the dummy atom natural bond length to ca. 1.5 Å gives a much more reasonable free energy difference. Unfortunately there is not one “perfect” length for all systems, and you should experiment with the dummy atom natural bond length to ascertain that the results obtained are independent of the dummy atom parameters.

17.4 Other Types of Free Energy Calculations

Free energy calculations have been possible for some time—most of the simulation methods in MacroModel result in ensembles that actually sample the free energy surface of the molecule. In some cases actual values for free energy differences between well-defined conformational states can be obtained. This is best illustrated with two examples from our own laboratory.

17.4.1 Hydrogen Bonding Preference of a Glycyl Lactam in Organic Solution

Infrared studies of the glycyl lactam in CH_2Cl_2 [30] suggest that both hydrogen bonded and non-hydrogen bonded forms are significantly populated at room temperature. In this case conformational searching (which can be considered an “Enthalpy” at 0 K) was apparently at variance with experiment because the global minimum, a hydrogen bonded form, is 2.7 kcal/mol lower in energy than the most stable non-hydrogen bonded energy minimum. When molecular dynamics simulations were performed at 300 K [31] and the hydrogen bond population was monitored during the simulation, we found that only 56% of the conformations could be considered hydrogen bonding at this temperature. Thus the hydrogen bonded populations reflect the experimentally observed populations because the simulation is exploring the “free energy surface” of the molecule. The “minima” obtained from the conformational search are the lowest points on very broad energy wells, which include many states lacking well-defined hydrogen bonds. So even using a high quality force field and solvation model was not sufficient to reproduce experiment without the use of the appropriate simulation methodology.




Figure 17.2. The glycyl lactam structure.

17.4.2 Conformational Free Energies for a Diamide in Solution

This compound was one of a series studied by Gellman et al. [32] using variable temperature IR and NMR in organic solvent. If we consider that there are two states of the molecule (hydrogen bound and non-hydrogen bound, each representing many possible conformations), we can perform a series of simulations at different temperatures in the range 200-300 K and perform a van't Hoff analysis on the results, which is directly comparable to that reported from the experiments [33]. The results of this analysis are shown in the following figure and table:

ΔH (kcal/mol)				ΔS (e.u.)		
N	^1H NMR	IR	Simulation	^1H NMR	IR	Simulation
4	-1.5	-1.4	-1.65	-8.0	-6.8	-8.6

In this system, where there are two well-defined states, we have been able to obtain thermodynamic parameters (which agree well with experiment) from direct simulation without using any perturbation techniques.

The reason that both of these simulations were able to obtain quantitative agreement with experiment is that we closely followed this procedure:

- Careful parametrization of the gas-phase conformational energy differences against ab initio calculations at the 6-31G* level.
- Use of the GB/SA solvent model.




Figure 17.3. The diamide structure.




Figure 17.4. van't Hoff analysis of the diamide.

- Assuring converged simulations by long runs (ns time scale) using the mixed-mode MC/SD procedure. We were also careful to test for convergence by ensuring that we could obtained the same results from different starting conformations.

17.5 Literature

This section is not an exhaustive review of the free energy literature, but rather a set of representative reviews and papers covering both theory and applications of free energy calculations. Those marked with a “*” are particularly important for their critical analysis of the potential problems of free energy calculations, and we urge you to read and understand these reviews before attempting free energy calculations.

Kollman, P. A. *Chem Rev.* **1993**, 93, 2395.

Beveridge, D. L.; DiCapua, F. M. *Annu. Rev. Biophys. Biophys. Chem.* **1989**, 18, 431.

Jorgensen, W. L. *Acc. Chem. Res.* **1989**, 22, 184.

Jorgensen, W. L. *J. Am. Chem. Soc.* **1989**, 111, 755.

Cieplak, P.; Kollman, P. A. *J. Comput.-Aided Mol. Des.* **1993**, 7, 291.

* van Gunsteren, W. F., Weiner, P. K., Eds. *Computer Simulations of Biomolecular Systems*; ESCOM: Leiden, 1990.

* van Gunsteren, W. F., Mark, A. E. *Eur. J. Biochem.* **1992**, 204, 947.

* Mark, A. E.; van Helden, S. P.; Smith, P. E.; Janssen, L. H. M.; van Gunsteren, W. F. *J. Am. Chem. Soc.* **1994**, 116, 6293.

* Mitchell, M. J.; McCammon, J. A. *J. Comput. Chem.* **1991**, 12, 271.

Smith, P. E.; van Gunsteren, W. F. *J. Phys. Chem.* **1994**, 98, 13735.

Molecular Clustering with XCluster

XCluster is a powerful structural clustering tool that uses molecular similarity as the clustering criterion. XCluster calculations can be set up in Maestro on a set of conformers contained in the Project Facility. Details about XCluster can be found in the [MacroModel XCluster Manual](#).

The Maestro interface makes it easy to select atoms to be used in the XCluster analysis. Instead of typing in the atom numbers in pairs, you can pick the atoms in the Workspace or select the comparison atoms using the flexibility of selection available in the Atom Selection dialog box. This is an important feature because more meaningful clustering can be obtained by choosing the minimum number of comparison atoms that are relevant for the clustering to be performed.

XCluster calculations, including the selection of the comparison atoms or torsions, are prepared and started from the Maestro XCluster panel. When you start the calculation from Maestro, the job is transferred to the original XCluster interface, which controls the actual computation and the display of the results of the XCluster analysis. As in the past, XCluster jobs can still be prepared and submitted from the original interface or from the command line. This interface can be launched with the command `$SCHRODINGER/xcluster`.

18.1 The XCluster Panel

To open the XCluster panel, choose XCluster from the Applications menu.

This panel allows you to specify the source of structural input, select the distance criterion to be used, select the comparison atoms or torsions, choose the type of RMS calculation, and specify whether enantiomers are to be considered. The following sections describe the use of the XCluster panel to set up an XCluster calculation.

When you have set up the calculation, click Start. The original XCluster interface is launched, the computation transferred to its control, and the job is run. When the job is finished, you can investigate the clustering statistics and use the visualization tools from the Visualize menu of the XCluster interface.

18.1.1 Selecting the Structure Source

To perform an XCluster analysis, you must provide a collection of conformers as input. The input can be taken from Maestro's project facility (for instance, from the incorporated structural output of a conformational search or dynamics simulation) or from an external structure file. XCluster cannot perform clustering calculations on non-conformers.




Figure 18.1. The XCluster panel.

To use structures from the current project, select the entries in the Project Table, then select Selected entries in the Source of job input section. To use structures that are in a disk file, select From file and type a file name in the Input file text box, or click Open to open a file selector to choose a file. The first structure in this file is read into Maestro and replaces the contents of the Workspace, so you can use the picking tools to set up the comparison atoms or torsions.

18.1.2 Selecting the Clustering Criterion

In the Cluster by section, you can select Atomic RMS or Torsional RMS as the distance criterion for the clustering calculation. The atom selection tools displayed in the center of the panel depend on which of these options you select. The options in the lower section of the panel, Calculate RMS 'In-Place' (no superposition) and Compare enantiomers, affect how the clustering is performed.

Atomic RMS

This option selects the atomic RMS displacement of corresponding atoms as the distance criterion. By default, a rigid-body superposition is performed before evaluating RMS displacements. The default is overridden by the Calculate RMS 'In-Place' (no superposition) option.

Torsional RMS

This option selects the root-mean-square difference between corresponding torsion angles in pairs of structures as the distance criterion.

Calculate RMS 'in-place' (no superposition)

This option is not selected by default. If this option is selected, there is no rigid-body superposition, but the distance matrix is computed without translating or rotating the input conformers. This option is useful if the input conformers have a relationship to another structure such as a receptor, for which the relative coordinates between conformers should not be adjusted in the clustering analysis.

Compare enantiomers

This option is selected by default. When selected, the analysis generates the enantiomer of each input conformer, and uses the smallest of the distances computed for the input structure and its enantiomer as the distance entry in the distance matrix for the conformer. If the option is not selected, enantiomers of input conformers are not considered in the analysis.

18.1.3 Defining the Comparison Atoms or Torsions

The XCluster panel makes it easy to designate the comparison atoms quickly. If you are using project entries as the source of the set of conformers, ensure that one conformer (and only one) is included in the Workspace. If you are using an external file as input source, the first structure in the file is automatically imported and included in the Workspace when you open the file.

If the distance criterion is Atomic RMS, there are several ways to select the comparison atoms. The buttons Heavy Atoms and Heavy Atoms + O-H, S-H select all atoms of the type indicated. Atoms can also be defined by picking in the Workspace. The final way to select comparison atoms is with the Atom Selection dialog box, which you open with the Select button. With this tool, complex atom sets can be easily and quickly selected.

If you selected Torsional RMS, you can choose either Atom or Bond from the Pick menu in the Define comparison torsions section, and define the torsions by picking.

The Delete and Delete All buttons can be used to edit or clear the comparison list.

18.2 Command File Examples

A sample XCluster command file is given below for an Atomic RMS displacement calculation. A subset of heavy atoms has been chosen to serve as the comparison atoms in the RMS calculation. Five lead structures are written to the output structure file.

```
# start sample command file
Sfile: cluster.mae
Mmsym:
Enant:
Arms:
  1
  5
  8
 11
 14
 17
 33
 39
 42
 45
 48
 52
```

```
Cluster:
Writelead: -5 cluster-out.mae all
# end sample command file
```

Sfile: Specifies the input structure file.

Mmsym: The MMSYM symmetry library will be used in the RMS difference calculation. MMSYM automatically recognizes local and global molecular symmetry.

Enant: Enables consideration of enantiomers of each conformer.

Arms: Atom numbers used when generating the distance matrix.

Cluster: Builds clusters at all clustering levels and calculates simple statistics for the clusterings.

Writelead: Writes the lead structures for each of five clusters to the output structure file.

For more information on these commands, see [Chapter 4](#) and [Chapter 8](#) of the *MacroModel XCluster Manual*.

18.3 XCluster Output

XCluster analysis is mainly carried out using the visualization tools in the XCluster interface after the cluster calculation is finished. In addition, structural output can be generated from the Write panel, which you can open from the File menu or from the visualization panels. A log file is generated after a successful XCluster computation as *jobname.clg*.

Redundant Conformer Elimination

Collections of conformers can contain structures that are essentially the same, based upon energetic or geometric considerations. The redundant conformer elimination facility allows you to use MacroModel to remove the extra conformers rapidly without reminimizing or re-evaluating the energy.

One important use for this facility is in the complementary use of MacroModel and Jaguar to conduct a very high quality conformational search. In this approach, MacroModel is used to perform the initial conformational search, typically using GB/SA solvation. To increase the accuracy of the geometries and energies of the structures produced, Jaguar is used to reminimize the collection of conformers. Sometimes the molecular mechanics potential energy surface differs enough from the ab initio potential energy surface that structures considered distinct in the original conformational search minimize to the same structure in the Jaguar calculations. The redundant conformer elimination facility can then be used to remove these redundant conformers without adjusting the energy or geometry of the molecules.

Another use of the redundant conformer elimination facility is to reprocess the results of a MacroModel conformational search with less strict criteria for distinguishing conformers or a smaller energy range for retained conformers. As well, on occasions MacroModel conformational searches might not eliminate all of the redundant conformers. Running the structures through this facility provides a rapid mechanism for dealing with such cases.

19.1 Eliminating Redundant Conformers Using Maestro

You can set up and run jobs to eliminate redundant conformers from a set of conformers with the Redundant Conformer Elimination panel. To open this panel, choose Redundant Conformer Elimination from the MacroModel submenu of the Applications menu in the main menu bar. The upper and lower parts of the panel are common to all MacroModel panels. These components are described in detail in [Section 5.1 on page 53](#). Redundant conformer elimination performs no energy calculations, so the folders present in other MacroModel panels are not present. The remaining controls in this panel are for selecting atoms and settings for structural comparison.

Define comparison atoms

You can define comparison atoms by picking atoms in one of the conformers, which must be included in the Workspace, or by selecting one of the predefined sets of comparison atoms: Heavy Atoms + O-H, S-H or Heavy Atoms. Heavy atoms are defined as non-hydrogen atoms.




Figure 19.1. The Redundant Conformer Elimination *panel*.

The first option of the two predefined sets includes hydrogen atoms attached to O and S atoms. The comparison atoms are listed as atom numbers in the Comparison atoms list. To delete one or more atoms from the list, select them in the list and click Delete. You can delete all comparison atoms from the list by clicking Delete All.

Maximum distance between atoms in equal structures

This text box specifies the threshold for determining whether structures are equivalent. When the structures are compared, the maximum distance between pairs of corresponding atoms must be less than this threshold for the structures to be considered equivalent. The default is 0.25 Å.

Compare structures 'in place' (no superposition)

If you select this option, the conformers are not superimposed (translated and rotated to obtain the closest alignment of atoms) in the process of comparing the structures, but the comparison is made with the structures in their given locations. By default, the structures are superimposed. This option could be used, for example, for Glide poses, where the position of the ligands with respect to the receptor is important.

Distinguish enantiomers

By default, both enantiomers of a structure are considered when comparing with the other structures. If one enantiomer matches, the structure is eliminated. To ensure that both enantiomers of a structure are kept, select this option.

Source of energy

If you choose to include the energy as a criterion for comparison, structures whose energy differs by more than 1 kJ/mol are considered inequivalent. This test is applied before the geometries are compared and can considerably speed up the comparison. You can choose from the Jaguar energy or the energies for various force fields supplied with MacroModel, or you can disable energy comparison by selecting None. If the particular energy you selected is not available among the properties of the structure, the calculation stops. You should therefore choose an energy that has been calculated for all structures.

Energy window for saving structures

This is the threshold value for comparison of structures. Structures are kept only if their energy is less than this value above the current global minimum. Lowering this value results in fewer structures saved. The default value is 50 kJ/mol.

19.2 Command File Examples

To run MacroModel calculations, a molecular structure file and a command file are required. The molecular structure file contains the structures to be used as input in the calculation. The command file contains the name of the input structure file, the name of the output structure file, and an ordered list of operation codes (*opcodes*) for the calculations.

Once you set up a job in Maestro and click either Start or Write, Maestro writes out a molecular structure file and a command file. For many types of jobs, command files written this way are complete and adequate. For some types of jobs, however, you may need to adjust the Maestro-generated command file.

Redundant conformer elimination is performed in MacroModel by the [ADDC](#) opcode. Each time ADDC is executed, it tries to add the current structure to the collection of output conformers. The same checks for redundancy and transformations, such as net translation and rotation of the current structure, are performed as in a multiple minimization of conformers (see [COMP](#)), but no minimization or energy estimation is conducted. Energies from earlier MacroModel or Jaguar calculations may be used as part of the comparison process (see [DEMX](#) in the *MacroModel Reference Manual*). If MacroModel energies are to be used, then a [FFLD](#) line with the same force field used to calculate the energies must precede the ADDC line. See the *MacroModel Reference Manual* for information on [ADDC](#).

ADDC should be used much like MINI in a multiple minimization calculation. The command file below reads in conformers from the file `addc.mae`, orders them by increasing Jaguar energy, superimposes the structures on the lowest Jaguar energy structure, eliminates redundant conformers, and writes out the resulting structures to `addc-out.mae`.

```
addc.mae
addc-out.mae
DEMX      0      0      0      0      10.0000      10.0000      0.0000      0.0000
MSYM      1      0      0      0      0.0000      0.0000      0.0000      0.0000
BGIN      0      0      0      0      0.0000      0.0000      0.0000      0.0000
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
AUTO      1      0      0      0      0.0000      0.0000      0.0000      0.0000
ADDC      -1      0      0      0      0.0000      0.0000      0.0000      0.0000
END        0      0      0      0      0.0000      0.0000      0.0000      0.0000
```

DEMX: Arg5 sets the energy window to 10 kJ/mol.

BGIN: Begin the loop over conformers.

READ: Read the next conformer.

AUTO: Automatically set up comparisons. Arg1 = 1 sets up comparisons only once, and uses this set for the remaining conformers.

ADDC: Check the current conformer. Arg1 = -1 indicates that Jaguar energies should be present in the input structure file and that these energies are to be used in the process that eliminates redundant conformers.

Additional Features

MacroModel has some additional features that are not available from the Maestro interface. Command file examples of some of these types of calculations are provided in this chapter. Note that some of these command files do not have the first `MMOD` line as seen in command files generated from the interface. The `MMOD` monitoring command is not necessary because these files cannot currently be monitored from Maestro.

20.1 Geometry Calculations

Geometric properties such as atom positions, distances, angles, and dihedrals can be measured using the `GEOM` opcode. Below is an example of the command file for a geometry calculation and explanations of the opcodes contained in the file.

```
geom.mae
geom-out.mae
  BGIN
  READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
  GEOM      14      15      16      17      0.0000      0.0000      0.0000      0.0000
  END
```

BGIN/END: Geometries can be calculated for a number of input structures within this loop.

READ: Read the input file.

GEOM: Use the number of non-zero arguments in `arg1-arg4` to define whether atom coordinates, distances, bond angles, or dihedral angles should be calculated. For all of `arg1-arg4` non-zero, the value of the dihedral angle defined by the atom numbers listed in `arg1-arg4` is reported. `Arg5` is used for calculation of spin-spin coupling constants. For details of `arg1-arg5` of the `GEOM` opcode, see [Section 3.15](#) of the *MacroModel Reference Manual*.

20.2 Calculating Interaction Energies Using ASET

ASET can be used to examine the interaction energies between sets of atoms within a structure. In the sample below, the interaction energy between LYS 33 and the ligand in the CDK2 structure `1e1v` is computed.

```
aset.mae
aset-out.mae
  DEBG      1      0      0      0      0.0000      0.0000      0.0000      0.0000
```

FFLD	11	1	0	0	1.0000	0.0000	0.0000	0.0000
BDCO	0	0	0	0	89.4427	99999.0000	0.0000	0.0000
ASET	466	487	0	0	1.0000	2.0000	0.0000	0.0000
ASET	4599	4633	0	0	2.0000	2.0000	0.0000	0.0000
ASNT	1	2	1	1	0.0000	0.0000	0.0000	0.0000
READ	0	0	0	0	0.0000	0.0000	0.0000	0.0000
ELST	-1	0	0	0	0.0000	0.0000	0.0000	0.0000

DEBG: Different debugging flags can be set. In this example, `arg1=1` tells MacroModel to write out the set relationships to the log file.

ASET: Defines atom sets. Arg1-arg4 lists atom-numbers of atoms that participate in a set. The set number is defined in arg5. Arg6 allows for numerous alternatives of defining sets using arg1-arg4. For instance, a range of atoms can be added to a set by using arg1 and arg2 (first and last atom number of atoms to be included in the set) only, with `arg6=2`, as in this example. Arg7 and arg8 control ASET properties written to the output structure file. For details on the ASET opcode, see [Section 3.5](#) of the *MacroModel Reference Manual*.

ASNT: Turn on/off set interactions. Arg1 and arg2 give the set numbers (as defined in ASET). Arg3=1 turns force field interactions on and arg4=1 turns constraint interactions on.

READ: Read in the structure.

ELST: Calculate the interaction energies between the specified sets and record only limited information on the interactions in the system to the `.log` file.

20.3 Rewinding the Output File for Additional Minimization

The **RWND** opcode, when included in a MacroModel command file, uses the input or output file again for a second computation. Currently this is supported only for a conformation search followed by a multiple minimization of the resulting conformers. This enables a conformation search to be performed with a higher minimization gradient specified for generated conformers, followed by minimization to a lower gradient with duplicate elimination, all in a single MacroModel computation. This workflow can be more efficient than minimizing all generated structures during the conformation search itself. The same potential settings should be used for both segments of the calculation.

The command file below combines a low-mode conformation search with a subsequent multiple minimization of the generated conformers to efficiently produce a set of low-energy conformers with a small convergence gradient.

```
rwnd.mae
rwnd-out.mae
```

MMOD	0	1	0	0	0.0000	0.0000	0.0000	0.0000
FFLD	11	1	0	0	1.0000	0.0000	0.0000	0.0000
SOLV	3	1	0	0	0.0000	0.0000	0.0000	0.0000
EXNB	0	0	0	0	0.0000	0.0000	0.0000	0.0000
BDCO	0	0	0	0	89.4427	99999.0000	0.0000	0.0000
READ	0	0	0	0	0.0000	0.0000	0.0000	0.0000
LMCS	1000	0	0	0	0.0000	0.0000	3.0000	6.0000
MCSS	2	0	0	0	50.0000	0.0000	0.0000	0.0000
MCOP	1	0	0	0	0.0000	0.0000	0.0000	0.0000
DEMX	0	1666	0	0	50.0000	100.0000	0.0000	0.0000
COMP	0	0	0	0	0.0000	0.0000	0.0000	0.0000
MSYM	0	0	0	0	0.0000	0.0000	0.0000	0.0000
CONV	2	0	0	0	0.0500	0.0000	0.0000	0.0000
MINI	1	0	5000	0	0.0000	0.0000	0.0000	0.0000
RWND	0	0	0	0	0.0000	0.0000	0.0000	0.0000
BGIN	0	0	0	0	0.0000	0.0000	0.0000	0.0000
READ	0	0	0	0	0.0000	0.0000	0.0000	0.0000
CONV	2	0	0	0	0.0200	0.0000	0.0000	0.0000
MINI	1	0	1300	0	0.0000	0.0000	0.0000	0.0000
END	0	0	0	0	0.0000	0.0000	0.0000	0.0000

All opcodes above the **RWND** command are necessary for a standard low-mode search. The **COMP** command indicates that all heavy atoms are to be used for the redundant conformer comparison. The minimizations, after conformation generation, are continued only until a gradient of 0.05 is reached, as seen in **arg5** of the **CONV** command.

Arg1 of the **RWND** command indicates that the output file of the low-mode search is to be used as the input of the subsequent multiple minimization. The second argument indicates that the intermediate conformation search output structure file is to be discarded, and only the final results of the minimization are retained.

Below the **RWND** command are the opcodes necessary for the subsequent multiple minimization. Another **CONV** is included which indicates, with the fifth argument, a tighter convergence than was used in the conformation search.

20.4 Visualizing Vibrational Modes

Vibrational modes from MacroModel calculations can be visualized in Maestro using the ePlayer. Two opcodes, **VIBR** and **VBR2**, generate structures for a set of selected vibrational modes. **VIBR** is used for small molecules, whereas **VBR2** is used for larger molecules of the order of seven residues or larger. **VBR2** utilizes the same techniques as a large-scale low-mode conformation search, which include using ARPACK routines to solve the large eigenvalue problem. The arguments to **VIBR** and **VBR2** specify the first and last mode to visualize, as well as the number of frames to include for the period of each vibration. Rotational and translational modes are automatically disregarded. The structural output consists of a series of frames

(structures) that together animate the vibration when imported and played in Maestro's ePlayer.

VIBR and VBR2 computations are started from the command line with an appropriate command file and an input structure file. The input structure should be minimized to a low gradient. Below is a sample command file for a VIBR calculation.

```
vibr.mae
vibr-out.mae
FFLD      10      1      0      1      1.0000      0.0000      0.0000      0.0000
SOLV       3      1      0      0      0.0000      0.0000      0.0000      0.0000
EXNB       0      0      0      0      0.0000      0.0000      0.0000      0.0000
BDCC       0      0      0      0      89.4427 99999.0000      0.0000      0.0000
READ       0      0      0      0      0.0000      0.0000      0.0000      0.0000
CONV       2      0      0      0      0.0200      0.0000      0.0000      0.0000
MINI       1      0      500     0      0.0000      0.0000      0.0000      0.0000
VIBR       1      5      5      0      3.0000      0.0000      0.0000      0.0000
```

The command file specifies `vibr.mae` as the input structure file, and the output structural animations will reside in the file `vibr-out.mae`. The following four opcodes specify the potential settings, here MMFFs as the force field with solvation, extended non-bonded cutoffs, and BDCC usage for electrostatic interactions. After the structure is read, a minimization is requested (CONV, MINI). The `VIBR` opcode selects the first five non-trivial modes for animation, with five frames per quarter vibration cycle.

20.5 Alignment of Structures

The combination of the `COPY` and `ALGN` opcodes permits very rapid approximate alignment of ligands with a reference ligand using center of mass and moments of inertia. These commands may be used to preposition ligands crudely prior to an eMBRACE conformational search calculation based upon the position of a reference ligand already positioned in the active site. While this combination of commands may be useful the positioning is crude, and searching conformational space is slow and quite limited compared to that available in Schrödinger's docking program, Glide.

Note that `COPY/ALGN` may require the purchase of a separate license.

Below is an example command file for a `COPY/ALGN` run to align structures in the input structure file (`algn.mae`). Opcode descriptions follow.

```
algn.mae
algn_all-out.mae
READ      0      0      0      0      0.0000      0.0000      0.0000      0.0000
COPY      0      0      0      0      0.0000      0.0000      0.0000      0.0000
BGIN      0      0      0      0      0.0000      0.0000      0.0000      0.0000
```

READ	0	0	0	0	0.0000	0.0000	0.0000	0.0000
ALGN	3	1	5	0	0.0000	0.0000	0.0000	0.0000
END	0	0	0	0	0.0000	0.0000	0.0000	0.0000

READ: The first READ command reads in the first structure, which is the reference structure. The remaining structures are read in a BGIN/END loop with another READ command.

COPY: Saves a copy of the current structure in internal reference arrays for later use.

BGIN/END: Loop over the remaining structures in the file.

ALGN: Aligns each structure with the reference structure by the center of mass and moments of inertia (arg1=3, arg2=1). All four equivalent alignment are generated (arg3=5).

20.6 autoref: Restrained Minimizations

A script is available for performing restrained minimizations. The primary application is for the minimization of protein-ligand structures just after hydrogen atoms have been added. In this procedure, the heavy atoms in the system are restrained using a harmonic potential while the hydrogen atoms do not have restraining potentials applied.

Before you use `autoref`, you must convert the structure file from Maestro format (`.mae`) to MacroModel format (`.dat`) using the structure conversion utility `maemmod`, as shown below.

```
$SCHRODINGER/utilities/maemmod basename.mae basename.dat
```

The syntax of the `autoref` command is as follows:

```
$SCHRODINGER/utilities/autoref [options] basename.dat
```

The output file is named `basename_ref.dat`. The command options are as follows:

- l *rmsd* Specify maximum RMS deviation from the initial structure; default is 0.3.
- m Use the MMFFs force field instead of the default OPLS_2001.
- k Keep output structure files from intermediate stages.
- v Print version number and exit.

20.7 serial_split: Split Serial Search Output Structure Files

A serial search produces an output structure file that contains the collections of conformations produced from a number of different molecules. The `serial_split` utility can be used to produce separate files, each containing the conformers produced for a different molecule. The

`serial_split` utility can be used on the output files of eMBrAcE conformational searches if there are only ligands or only complexes in the file. The syntax of the command is:

```
SCHRODINGER/utilities/serial_split [options] input.mae basename
```

where *input.mae* contains the results of the previously performed serial conformational search. The options are given in [Table 20.1](#).

Table 20.1. Options for the serial_split utility.

Option	Description
-h	Print usage summary and exit.
-b	The serial number of the first molecule for which conformations should be extracted. If omitted, start with the first serial number present.
-e	The serial number of the last molecule for which conformations should be extracted. If omitted, continue to the last serial number present.
-v	Print the version number.

The files produced are named *basename*<serial number>.mae where <serial number> is the `i_mmod_Serial_Number` property listed for each structure in the input structure file. This number is usually a running count of the distinct structures present.

Getting Help

Schrödinger software is distributed with documentation in PDF format. If the documentation is not installed in `$SCHRODINGER/docs` on a computer that you have access to, you should install it or ask your system administrator to install it.

For help installing and setting up licenses for Schrödinger software and installing documentation, see the *Installation Guide*. For information on running jobs, see the *Job Control Guide*.

Maestro has automatic, context-sensitive help (Auto-Help and Balloon Help, or tooltips), and an online help system. To get help, follow the steps below.

- Check the Auto-Help text box, which is located at the foot of the main window. If help is available for the task you are performing, it is automatically displayed there. Auto-Help contains a single line of information. For more detailed information, use the online help.
- If you want information about a GUI element, such as a button or option, there may be Balloon Help for the item. Pause the cursor over the element. If the Balloon Help does not appear, check that Show Balloon Help is selected in the Help menu of the main window. If there is Balloon Help for the element, it appears within a few seconds.
- For information about a panel or the folder that is displayed in a panel, click the Help button in the panel. The Help panel is opened and a relevant help topic is displayed.
- For other information in the online help, open the Help panel and locate the topic by searching or by category. You can open the Help panel by choosing Help from the Help menu on the main menu bar or by pressing CTRL+H.

To view a list of all available MacroModel–related help topics, choose MacroModel from the Categories menu of the Categories tab. Double-click a topic title to view the topic.

If you do not find the information you need in the Maestro help system, check these sources:

- *Maestro User Manual*, for detailed information on using Maestro
- *Maestro Tutorial*, for a tutorial on the basic features of Maestro
- *Maestro Command Reference Manual*, for information on Maestro commands
- *MacroModel Quick Start Guide*, for a tutorial guide to using MacroModel
- *MacroModel Reference Manual*, for information on MacroModel commands
- Frequently Asked Questions pages, at https://www.schrodinger.com/MacroModel_FAQ.html

The manuals are also available in PDF format from the Schrödinger [Support Center](#). Information on additions and corrections to the manuals is available from this web page.

If you have questions that are not answered from any of the above sources, contact Schrödinger using the information below.

E-mail: help@schrodinger.com

USPS: Schrödinger, 101 SW Main St. Suite 1300, Portland, OR 97204

Phone: (503) 299-1150

Fax: (503) 299-4532

WWW: <http://www.schrodinger.com>

FTP: <ftp://ftp.schrodinger.com>

Generally, e-mail correspondence is best because you can send machine output, if necessary. When sending e-mail messages, please include the following information, most of which can be obtained by entering `$SCHRODINGER/machid` at a command prompt:

- All relevant user input and machine output
- MacroModel purchaser (company, research institution, or individual)
- Primary MacroModel user
- Computer platform type
- Operating system with version number
- MacroModel version number
- Maestro version number
- mmshare version number

References

1. Allinger, N. L. Conformational Analysis 130. MM2. A Hydrocarbon Force Field Utilizing V1 and V2 Torsional Terms. *J. Am. Chem. Soc.* **1977**, *99*, 8127.
2. Allinger, N. L.; Yuh, Y. H.; Lii, J-H. Molecular Mechanics. The MM3 Force Field for Hydrocarbons. *J. Am. Chem. Soc.* **1989**, *111*, 8551.
3. Ferguson, D. M.; Kollman, P. A. Can the Lennard-Jones 6-12 function replace the 10-12 form in molecular mechanics calculations? *J. Comput. Chem.* **1991**, *12*, 620.
4. McDonald, D. Q.; Still, W. C. AMBER torsional parameters for the peptide backbone. *Tetrahedron Lett.* **1992**, *33*, 7743.
5. Jorgensen, W. L.; Tirado-Rives, J. The OPLS Potential Functions for Proteins. Energy Minimization for Crystals of Cyclic Peptides and Crambin. *J. Am. Chem. Soc.* **1988**, *110*, 1657.
6. Jorgensen, W. L.; Maxwell, D. S.; Tirado-Rives, J. Development and Testing of the OPLS All-Atom Force Field on Conformational Energetics and Properties of Organic Liquids. *J. Am. Chem. Soc.* **1996**, *118*, 11225-11235.
7. Kaminski, G. A.; Friesner, R. A.; Tirado-Rives, J.; Jorgensen, W. L. Evaluation and Reparametrization of the OPLS-AA Force Field for Proteins via Comparison with Accurate Quantum Chemical Calculations on Peptides. *J. Phys. Chem. B* **2001**, *105*, 6474-6487.
8. Cornell, W. D.; Cieplak, P.; Bayly, C. I.; Gould, I. R.; Merz, K. M.; Ferguson, D. M.; Spellmeyer, D. C.; Fox, T.; Caldwell, J. W.; Kollman, P. A. A second generation force field for the simulation of proteins, nucleic acids and organic molecules. *J. Am. Chem. Soc.* **1995**, *117*, 5179.
9. Halgren, T. A. Merck Molecular Force Field. I. Basis, Form, Scope, Parameterization and Performance of MMFF94. *J. Comput. Chem.* **1996**, *17*, 490-519.
10. Halgren, T. A. Merck Molecular Force Field. II. MMFF94 van der Waals and Electrostatic Parameters for Intermolecular Interactions. *J. Comput. Chem.* **1996**, *17*, 520-552.
11. Halgren, T. A. Merck Molecular Force Field. III. Molecular Geometries and Vibrational Frequencies for MMFF94. *J. Comput. Chem.* **1996**, *17*, 553-586.
12. Halgren, T. A.; Nachbar, R. B. Merck Molecular Force Field. IV. Conformational Energies and Geometries. *J. Comput. Chem.* **1996**, *17*, 587-615.

13. Halgren, T. A. Merck Molecular Force Field. V. Extension of MMFF94 using Experimental Data, Additional Computational Data and Empirical Rules. *J. Comput. Chem.* **1996**, *17*, 616-641.
14. Halgren, T. A. MMFF VI. MMFF94s Option for Energy Minimization Studies. *J. Comput. Chem.* **1999**, *20*, 720-729.
15. Halgren, T. A. MMFF VII. Characterization of MMFF94, MMFF94s and Other Widely Available Force Fields for Conformational Energies and for Intermolecular Interaction Energies and Geometries. *J. Comput. Chem.* **1999**, *20*, 730-748.
16. Still, W. C.; Tempczyk, A.; Hawley, R. C.; Hendrickson, T. A General Treatment of Solvation for Molecular Mechanics. *J. Am. Chem. Soc.* **1990**, *112*, 6127.
17. Polak, E.; Ribiere, G. Note sur la Convergence de Méthodes de Directions Conjuguées. *Revue Francaise Informat. Recherche Operationelle, Serie Rouge*, **1969**, *16*, 35.
18. Ponder, J. W.; Richards, F. M. An Efficient Newton-like Method for Molecular Mechanics Energy Minimization of Large Molecules. *J. Comput. Chem.* **1987**, *8*, 1016.
19. Oren, S. S.; Spedicato, E. Optimal conditioning of self-scaling variable metric algorithms. *Mathematical Programming* **1976**, *10*, 70.
20. Chang, G.; Guida, W.; Still, W. C. An internal coordinate Monte-Carlo method for searching conformational space. *J. Am. Chem. Soc.* **1989**, *111*, 4379.
21. Saunders, M.; Houk, K. N.; Wu, Y.-D.; Still, W. C.; Lipton, M.; Chang, G.; Guida, W. C. Conformations of Cycloheptadecane: A Comparison of Methods for Conformational Searching. *J. Am. Chem. Soc.* **1990**, *112*, 1419.
22. Goodman, J. M.; Still, W. C. Searching Conformation Space. *J. Comput. Chem.* **1991**, *12*, 1110.
23. Kolossváry, I.; Guida, W. C. Low Mode Search. An efficient, automated computational method for conformational analysis: application to cyclic and acyclic alkanes and cyclic peptides. *J. Am. Chem. Soc.* **1996**, *118*, 5011.
24. Kolossváry, I.; Keseru, G. M. Hessian-Free Low-Mode Conformational Search for Large-Scale Protein Loop Optimization: Application to c-jun N-Terminal Kinase JNK3. *J. Comput. Chem.* **2001**, *22*, 21.
25. Keseru, G. M.; Kolossváry, I. Fully Flexible Low-Mode Docking: Application to Induced Fit in HIV Integrase. *J. Am. Chem. Soc.* **2001**, *123*, 12708.
26. Kolossváry, I.; Guida, W. C. Low-mode Conformational Search Elucidated. Application to C₃₉H₈₀ and Flexible Docking of 9-Deazaguanine Inhibitors to PNP. *J. Comput. Chem.* **1999**, *20*, 1671.

-
27. van Gunsteren, W. F.; Berendsen, H. J. C. A leap-frog algorithm for stochastic dynamics. *Mol. Simul.* **1988**, *1*, 173.
 28. Ryckaert, J.-P.; Ciccotti, G.; Berendsen, H. J. C. Numerical integration of the Cartesian equation of motion of a system with constraints: molecular dynamics of N-alkanes. *J. Comput. Phys.* **1977**, *23*, 327.
 29. Guarnieri, F.; Still, W. C. A Rapidly Convergent Simulation Method: Mixed Monte Carlo/Stochastic Dynamics. *J. Comput. Chem.* **1994**, *15*, 1302-1310.
 30. Dado, G. P.; Gellman, S. H. On the Use of AM1 Calculations for the Study of Intramolecular Hydrogen Bonding Phenomena in Simple Amides. *J. Am. Chem. Soc.* **1992**, *114*, 3138.
 31. McDonald, D. Q.; Still, W. C. An AMBER* study of Gellman's amides. *Tetrahedron Lett.* **1992**, *33*, 7747.
 32. Gellman, S. H.; Dado, G. P.; Liang, G-B.; Adams, B. R. Conformation-Directing Effects of a Single Intramolecular Amide-Amide Hydrogen Bond: Variable Temperature NMR and IR Studies on a Homologous Diamide Series. *J. Am. Chem. Soc.* **1991**, *113*, 1164.
 33. McDonald, D. Q.; Still, W. C. Conformational Free Energies from Simulation: Stochastic Dynamics/Monte Carlo Simulations of a Homologous Series of Gellman's Diamides. *J. Am. Chem. Soc.* **1994**, *116*, 11550.

Copyright Notices

C and C++ Libraries for Parsing PDB Records

The C and C++ libraries for parsing PDB records are Copyright (C) 1989 The Regents of the University of California. All rights reserved.

Redistribution and use in source and binary forms are permitted provided that the above copyright notice and this paragraph are duplicated in all such forms and that any documentation, advertising materials, and other materials related to such distribution and use acknowledge that the software was developed by the University of California, San Francisco. The name of the University may not be used to endorse or promote products derived from this software without specific prior written permission. THIS SOFTWARE IS PROVIDED “AS IS” AND WITHOUT ANY EXPRESS OR IMPLIED WARRANTIES, INCLUDING, WITHOUT LIMITATION, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

A

alignment of structures.....	200
AMBER*	32
compatible H-treatment	43
AMBER94	32
angles	
dihedral driving.....	89
monitored in dynamics calculations	132
arguments, opcode	47
atom flexibility	60
atom positions	
constraining	60
fixing and freezing	55
atom sets, interaction energies between	197
atoms, selecting.....	23–25
Auto-Help	30, 203
automatic setup	
in CSearch calculations	103
in MC/SD calculations	139

B

Balloon Help	30, 203
BDCO	36, 38
Build panel.....	21
building structures.....	20–23
button menu	9

C

charges, using from structure files	55
chiral atoms, selecting for CSearch	110
citing MacroModel in publications	4
command file.....	45
embedded in shell script	49
format	46
command file examples	
alignment of structures	200
ASET	197
ConfGen	124
CSearch.....	113
current energy	68
dihedral drive	90
dynamics.....	136
eMBrAcE.....	160
eMBrAcE conformational search	165
eMBrAcE, distributed.....	172
free-energy perturbation	177

geometry calculations.....	197
LOOP.....	149
MC/SD calculations.....	140
minimization	76
MINTA	145
multiple minimization	85
partition coefficient estimation	87
redundant conformer elimination	196
rewinding output file for minimization.....	198
vibrational animation.....	200
Command Script Editor panel.....	26
command scripts— <i>see</i> scripts	
comparison atoms	
in CSearch calculations	109
defining for XCluster	191
in multiple minimization calculations	81
ConfGen calculations.....	119
command file examples	124
number of steps.....	121
ConfGen panel— <i>see</i> Ligand Torsion Search panel	
conformational search— <i>see</i> CSearch calculations	
conformer generation.....	119
constraints	
on bond lengths using SHAKE.....	134
for fixing atoms.....	56
on positions, distances, angles, torsions	60
removing	60
Constraints folder.....	59
conventions, document.....	xiii
convergence guidelines	39
for CSearch jobs	39
for MCSD jobs	40
for minimization jobs	39
convergence threshold.....	76
CSearch calculations.....	97, 119
automatic setup of.....	103
chiral atoms	110
command file examples	113
comparison atoms	109
defining ring closures	105
defining torsion rotations	106
distance check.....	111
energy window.....	102
ligand bonds.....	113
markers	104
number of steps.....	101, 102
number of steps per rotatable bond	104

- resetting variables 103, 105
- torsion check 112
- translations and rotations 107
- variables 103
- CSearch methods
 - LMCS 97
 - MCMM 97
 - MCMM/LMCS 98
 - selecting 99
 - SPMC 97
- current energy calculations 67
- Current Energy panel 67
- current working directory 6
- D**
 - dielectric constant, modifying 55
 - dihedral drive calculations
 - command file examples 90
 - defining angles and increment 89
 - dihedral, defining for dihedral drive 89
 - directory
 - current working 6, 27
 - output 27
 - distances
 - atomic RMS for XCluster 190
 - checked in CSearch 111
 - constraining 60
 - monitored in dynamics calculations 131
 - distributed calculations 50
 - dynamics calculations 129
 - command file examples 136
 - equilibration time 135
 - method 134
 - monitored angles 132
 - monitored dihedrals 132
 - monitored distances 131
 - monitored H-bonds 133
 - monitored surface areas 131
 - SHAKE 134
 - simulation time 135
 - temperature 134
 - timestep 135
 - Dynamics folder 134
 - Dynamics panel 129
- E**
 - ECalc folder 67
 - electrostatic interaction truncation 36
 - electrostatic treatments 42
 - eMBrAcE calculations 153
 - association energy mode 155
 - command file examples 160
 - receptor entry 154
 - sample output 162
 - enantiomers, comparison of in XCluster 191
 - entries, Project Table 13
 - including, excluding, and fixing 18
 - selecting 17
 - sorting 15
 - specifying for receptor in eMBrAcE 154
 - environment variables
 - DISPLAY 6
 - SCHRODINGER 5–6, 45
 - SCHRODINGER_TMPDIR 45
 - ePlayer 15, 16
 - equilibration time 135
 - excluded entries 18
- F**
 - FEP calculations 175
 - hints for success 182
 - literature references for 186
 - preparing 176
 - sample output 180
 - structure file 176
 - file conversion 42
 - file I/O directory 27
 - files
 - atom type 43
 - force field 43
 - solvent 43
 - filters, project entry 17
 - fixed entries 18
 - Fixing and Freezing Atoms 55
 - FMNR minimization method 75
 - force fields 31
 - files 31
 - selecting 54
 - fragments, building structures from 20
 - free energy perturbation calculations—*see* FEP calculations

full screen mode 8, 13
 function key macros—*see* scripts

G

GB/SA 35
 geometry calculations 197
 grow bond 21

H

H-bonds monitored in dynamics calculations. 133
 Help panel 30, 203
 hydrogen treatment 43

I

included entries 18
 interaction energies using ASET 197
 iteration limit, minimization calculations 76

J

job files, writing 65
 job submission 65
 command line 47
 remote and distributed 50
 jobs
 general setup 53
 host 64
 incorporation preference 64
 input source 53
 monitoring 51, 65
 name 64
 remote login name 65
 jobs, running in Maestro 28–29

L

LBFGS minimization method 75
 ligand bonds checked in CSearch calculations 113
 Ligand Torsion Search panel 120
 LMCS 97
 log file, saving Maestro 30
 LOOP calculations 147
 command file examples 149
 example output 151
 input restrictions 147
 structure handling 148

low-mode conformational search methods 97
 low-mode settings 102
 .lsq file example 148
 LTSearch folder 122

M

macros—*see* scripts

Maestro

 main window 6, 7
 menus 8
 quitting 30
 running jobs from 28–29
 scratch projects 13
 starting 6
 undoing operations 28
 main window 7
 MBAE opcode 161, 164
 MC/SD calculations 139
 automatic setup 139
 command file examples 140
 MC to SD ratio 140
 MC/SD panel 139
 MCMM 97
 MCMM/LMCS 98
 menu button 9
 minimization calculations 73
 command file examples 76
 convergence parameters 75
 minimization methods 74
 FMNR 75
 LBFGS 75
 OSVM 75
 PRCG 74
 SD 75
 TNCG 74
 Minimization panel 73
 MINTA calculations 143
 command file examples 145
 energy evaluations per iteration 144
 hard limit 144
 input file 143
 number of iterations 143
 soft limit 145
 temperature 144
 MINTA panel 143
 mixed MCMM/Low-mode search method 98
 MM2* 31

MM3*	32	shortcut keys	19
MMFF	33	projects	13
Monitor folder	129	publications, citing MacroModel in	4
Monitor panel	29	Python scripts— <i>see</i> scripts	
Monte Carlo Multiple Minimum search method ...			
97			
mouse functions	5	Q	
Project Table panel	18–19	quitting Maestro	30
Workspace	12		
multi-ligand structure files, preminimization of	84	R	
multiple minimization calculations	79	receptor	
comparison atoms for	81	specifying entry for eMBrAcE	154
energy window	81	substructure in eMBrAcE	159
input file for	80	remote calculations	50
threshold for equivalent structures	82	results	
Multiple Minimization panel	79	dependence on task distribution	51
		results, checking and interpreting	
N		CSearch	118
nonbonded cutoffs	41	current energy	71
setting	55	dihedral drives	95
		dynamics	138
O		minimization	77
online help	30	MINTA	146
opcodes	47		
OPLS*	32	S	
OPLS_2001	32	Schrödinger contact information	204
compatible H-treatment	43	scratch entries	14
OPLS_2005	32	scratch projects	13
OPLSAA	32	scripts	
OSVM minimization method	75	function key macros	27
		macros	27
P		Maestro command	26
para_bmin	50	Python	25
parameter quality considerations	33	SD minimization method	75
partition coefficient estimation	82	selecting objects in the Workspace	9, 23
example	87	shells, defining	62
Potential folder	54	shortcut keys	
PRCG minimization method	74	main window	13
Preferences panel	27, 28	Project Table panel	19
premin script	84	simulation temperature	134
product installation	203	simulation time	135
project entries, <i>see</i> entries, Project Table		size limitations	41
Project Facility, introduction	13	solvation treatment	
Project Table panel	15	description	35
menus	16	selecting	54
mouse functions	18–19	SPMC	97
		Start dialog box	64

-
- structures
 - alignment of..... 200
 - building..... 20–23
 - displaying in sequence..... 15
 - substructure
 - defining..... 61
 - filling out residues in..... 62
 - for freezing or fixing atoms..... 56
 - shell radius..... 62, 63
 - shells..... 62
 - specifying for eMBrAcE..... 158
 - substructure files
 - format..... 58
 - reading and writing..... 63
 - Substructure folder..... 57, 61
 - superposition, rigid body..... 191
 - Systematic Pseudo-Monte Carlo search method.... 97
- T**
- technical support..... 30
 - temperature
 - dynamics calculations..... 134
 - MC/SD calculations..... 139
 - MINTA calculations..... 144
 - TNCG minimization method..... 74
 - toolbar
 - Build panel..... 22–23
 - main window..... 9–12
 - Project Table panel..... 15–16
 - torsional search..... 119
 - torsions
 - checked in CSearch..... 112
 - for comparison in XCluster..... 191
 - constraining..... 60
 - in MC/SD calculations..... 139
 - monitored in dynamics calculations..... 132
 - sampling amide, ester, azo, and imido..... 104
 - translations and rotations in MC/SD calculations.. 139
- U**
- undoing Maestro operations..... 28
- V**
- variable torsions in MC/SD calculations..... 140
 - vibrational modes, visualizing..... 199
- W**
- Workspace
 - description..... 6
 - full screen mode..... 8, 13
 - including, excluding, and fixing entries..... 18
 - mouse functions..... 12
 - scratch entries..... 14
- X**
- XCluster calculations..... 189

120 West 45th Street
32nd Floor
New York, NY 10036

101 SW Main Street
Suite 1300
Portland, OR 97204

3655 Nobel Drive
Suite 430
San Diego, CA 92122

Dynamostraße 13
68165 Mannheim
Germany

QuatroHouse, Frimley Road
Camberley GU16 7ER
United Kingdom

SCHRÖDINGER.